



One Touch Control



Operable with one hand, this ergonomic torque device is designed to enhance manipulation of the guide wire tip.



Sea Dragon™



TORQUE DEVICE

Saves time. No more tightening and untightening to reposition the device.

Simply squeeze the top cap with one hand to position, reposition or remove the device.

FEATURES & BENEFITS

PLASTIC COMPONENTS: Minimizes potential damage to the hydrophilic coating and the polymer jacket.

VERSATILITY: Device can be used for both hydrophilic and standard spring guide wires.

SMART DESIGN: Operable with one hand, squeeze the top cap to load, move and remove the device. Release the top cap to torque the wire.

EASY GRIP RIDGES: Enhances steering and torque performance.

QUICK RELEASE LOCKING MECHANISM: Device is quickly and easily removed from the guide wire facilitating wire/catheter placement and exchange.



SQUEEZE:

Squeeze the white top cap toward the green body to align lumens



LOAD:

Insert proximal guide wire tip through the device



RELEASE:

Release white top cap to activate locking mechanism



TORQUE:

Easy grip ridges enhances steering and torque performance

ORDERING INFORMATION

Cat. No.	Description
TD500	Guide wire torque device
	Accommodates guide wires from 0.018" (0.46mm) to 0.038" (0.097mm)
	Packaged individually, 10 units per box.



Merit Medical Systems, Inc. • 1600 West Merit Parkway, • South Jordan, Utah 84095 • 1-801-253-1600 • 1-800-35-MERIT

Merit Medical Europe • Horsterweg 24 • 6199 AC MAASTRICHT- AIRPORT • The Netherlands • +31 43 35 88 222

Merit Medical Ireland Ltd. • Parkmore Business Park West • Galway, Ireland • +353 (0) 91 703 733

Free phone from specific country:

Belgium 0800 72906 • France 0800 91 60 30 • Germany 0800 182 0871 • Luxembourg 0800 22522 • Netherlands 0800 022 8184 • UK 0800 973 115

Patents Pending

www.merit.com

400889001/A

ID 090707