

GRANULOMETRIC DISTRIBUTION OF PARTICLE SIZES:



A comparison between Merit Bearing nsPVA® Embolization Particles and Boston Scientific Contour™ Embolization Particles

Merit Medical Research & Development

Abstract

The standard manufacturing process of nonspherical polyvinyl alcohol particles starts with the fabrication of a polyvinyl alcohol sponge that is then ground into small nonuniform particles. The resulting shavings are passed through sieves with sequentially smaller holes to separate the particles into various size ranges. Due to the irregular shape of nsPVA particles, this common sieving method results in a percentage of particles that fall outside of specified size distribution ranges.

Objective

This granulometry analysis study was performed to compare the granulometric size distribution of Merit Bearing nsPVA Embolization Particles and Boston Scientific Contour Embolization Particles.

The following embolization particles and equipment were used:

- Merit Bearing nsPVA Embolization Particles (Fig. 1)
- Boston Scientific Contour Embolization Particles (Fig. 1)
- Camsizer® XT particle analyzer (Fig. 2)

Method

Seven size ranges (45-1180 μm) of Contour particles were purchased from Boston Scientific for this study. Merit tested the seven equivalent Bearing size ranges for this study.

The Camsizer XT particle analyzer performed an optical measurement analysis. Free-falling particles were illuminated by two pulsed LED light sources and their images were captured by two digital image processing cameras, each one specialized in a specific size range.

A method was specifically developed for each size range where image processing and calculations are relative to the dimensions, the velocity, and the sharpness of the particles. The Camsizer XT software stores and processes images of individual particles and can detect agglomerated particles and exclude them from the results.

The Camsizer XT performed an analysis on the number of particles that fell within and outside the stated size range.




Figure 1: Embolization Particles



Figure 2: Camsizer XT

RESULTS

The following charts indicate that the tested samples of Merit Bearing nsPVA[®] exhibited a greater percentage of particles consistently within the specified size ranges when compared to Boston Scientific Contour[™].




CONCLUSION


This study demonstrates that the size distribution of Bearing nsPVA particles is more tightly calibrated than Contour particles.

Disclaimer:


The tests described in this document were laboratory bench tests conducted by Merit Medical in 2012-2013. Granulometry analysis was performed with calibrated imaging equipment (Camsizer XT) on at least 6100 particles. These tests were conducted using the following lots of Contour[™] product: 14945627, 14959129, 14999768, 15228667, 15837033, 14792749, 15857804, 14860137, 15808410, 15008956, 15193505, 15855152, 14621269, 14996550, 15842238. Every effort was made to report accurate, verifiable results. No conclusions are drawn regarding the translation of bench data to the clinical model.



For all seven size ranges, on average, **Bearing nsPVA** had **33%** more particles in the specified size ranges than Contour.



In the 355-500 µm size range, **Bearing nsPVA** had **71%** and **Contour** had **51%** fall within range.



In the 500-710 µm size range, **Bearing nsPVA** had **72%** and **Contour** had **38%** fall within range.



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.208.4300
1.800.35.MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)
Amerikalaan 42, 6199 AE Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733

MeritEMEA.com

Austria
0800 295 374

Belgium
0800 72 906 (Dutch)
0800 73 172 (Français)

Denmark
80 88 00 24

Finland
0800 770 586

France
0800 91 60 30

Germany
0800 182 0871

Ireland (Republic)
1800 553 163

Italy
800 897 005

Luxembourg
8002 25 22

Netherlands
0800 022 81 84

Norway
800 11629

Portugal
308 801 034

Russia
+7 495 221 89 02

Spain
+34 911238406

Sweden
020 792 445

Switzerland
(Deutsch)
+41 225180252
(Français)
+41 225948000
(Italiano)
+41 225180035
UK
0800 973 115