BioSphere®
Medical

Biosphere Medical, Inc.
1050 Hingham Street
Rockland, MA 02370
USA

Tel: 800.394.0295
Fax: 781.792.2745

EC REP :

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<http://www.biospheremed.com>

Distributed in USA by:

Merit Medical Systems, Inc.
Customer Service 1-800-356-3748
781-681-7900

Instructions for use	English	2
Instructions d'utilisation.....	Français (French).....	4
Инструкции за употреба	Български (Bulgarian).....	6
Pokyny pro použití	Čeština (Czech).....	8
Brugsanvisning	Dansk (Danish).....	10
Gebruiksaanwijzing	Nederlands (Dutch).....	12
Kasutusjuhised.....	Eesti (Estonian).....	14
Käyttöohjeet.....	Suomi (Finnish).....	16
Gebrauchsanweisungen	Deutsch (German).....	18
Οδηγίες χρήσης.....	Ελληνικά (Greek).....	20
Használati útmutató	Magyar (Hungarian)	22
Istruzioni per l'uso.....	Italiano (Italian).....	24
Lietošanas norādījumi	Latviešu (Latvian)	26
Naudojimo instrukcijos.....	Lietuviškai (Lithuanian).....	28
Instruksjoner for bruk	Norsk (Norwegian)	30
Instrukcja stosowania	Polski (Polish)	32
Instruções de utilização.....	Português (Portuguese)	34
Instrucțiuni de utilizare	Română (Romanian)	36
Инструкции по применению	Русский (Russian)	38
Návod na použitie.....	Slovenčina (Slovak).....	40
Instrucciones de utilización	Español (Spanish).....	42
Brugsanvisning.....	Svenska (Swedish)	44
Kullanma talimatları	Türkçe (Turkish).....	46

DESCRIPTION

The EmboCath® Plus Infusion Microcatheter is a tapered 3.0-2.9F single lumen catheter designed to facilitate the access of distal vasculature over a guidewire. The catheter has a semi-rigid proximal shaft and becomes progressively more flexible toward the distal end. The shaft is reinforced, which provides improved torque transmission. The inner lumen is lined with a lubricious material to facilitate the movement of guidewires or other devices. The outer diameter of the catheter is coated with a hydrophilic surface to enhance the ability to navigate tortuous anatomy. The distal tip of the catheter has a single radiopaque marker to facilitate fluoroscopic visualization. The hub at the proximal end incorporates a standard luer adapter to facilitate the attachment of accessories. The catheter lumen is 0.028 inches and guidewires measuring up to 0.025 inches (0.635 mm) in diameter are recommended.

INDICATIONS FOR USE

The EmboCath Plus Infusion Microcatheter is intended for: infusion of various diagnostic, embolic and therapeutic agents into the body's vascular systems (peripheral); guidewire exchange/support; and superselective angiography of the peripheral vasculature. Diagnostic, embolic, or therapeutic agents are to be used according to specifications outlined by the manufacturer.

Note: Not intended for pediatric or neonatal use.

WARNING

- Contents supplied sterile.
- Do not use if pouch is opened or damaged. Verify that the package integrity has been maintained to ensure the sterility of the device.
- Do not use in case of any surface irregularities, bends, or kinks. Any damage of the micro-catheter may change its characteristics likely to affect its performance
- Use the device prior to the "Use Before" date noted on the package.
- For single patient use only. Do not reuse, reprocess, or resterilize. Reusing, reprocessing or resterilizing may compromise the structural integrity of the device and or lead to device failure, which in turn may result in patient injury, illness, or death. Reusing, reprocessing or resterilizing may also create a risk of contamination of the device and or cause patient infection or cross infection including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
- Discard the product after use.
- This device should be used only by physicians thoroughly trained in percutaneous, intravascular techniques and procedures in relevant areas of the anatomy.
- Do not attempt to move the catheter without observing the catheter tip. Always maintain visualization of the catheter under fluoroscopy.

- Do not advance or withdraw against resistance. If resistance is met, discontinue movement, determine the reason for resistance, and take appropriate action before continuing. Movement of the catheter or guidewire against resistance may result in separation of the catheter or guidewire tip, damage to the catheter, or vessel perforation.
- **Infusion dynamic pressure with this catheter should not exceed 800 psi (54.4 atm/bar). The static pressure with this catheter should not exceed 300 psi (20 atm/bar).** Pressure in excess of these maximum pressures may result in catheter rupture or device failure, possibly resulting in patient or user injury. If the flow of material through the catheter becomes restricted, do not attempt to clear the catheter lumen by infusion. Identify and resolve the cause of the blockage or replace the catheter with a new catheter before resuming infusion.
- Excessive tightening of a hemostatic valve onto the catheter shaft may result in damage to the catheter.
- The hydrophilic coating has an extremely lubricious surface only when properly hydrated.

PRECAUTIONS

- Do not expose Microcatheters to extreme temperatures or solvents.
- Prior to use, carefully examine the Microcatheter and packaging to verify proper function and integrity.

COMPLICATIONS

Procedures requiring percutaneous catheter introduction should not be attempted by physicians unfamiliar with the possible complications listed below. Possible complications may include, but are not limited to, the following:

- Puncture site hematoma
- Infection
- Ischemia
- Perforation of vessel or arterial wall
- Vasospasm
- Emboli

Flow Rate Table

EmboCath Plus Infusion Microcatheter Flow Rate (cc/second)

	Length (cm)	Dead Space Volume (cc)	Saline	100% Ionic Contrast	100% Non-Ionic Contrast
300 psi	100	0.5	3.8	1.1	1.6
	135	0.6	3.3	0.9	1.3
800 psi	100	0.5	9.3	2.9	3.7
	135	0.6	7.0	2.3	3.0

Flow rate values are approximate

PREPARATION FOR USE

- Open carefully the pouch
- Flush the hoop with sterile heparinized saline to activate the hydrophilic coating of the catheter. The luer fitting attached to the hoop will facilitate the flushing of the hoop.
- Flush the inner lumen of the catheter with heparinized saline by attaching a saline-filled syringe to the catheter hub.

- Remove the Microcatheter from the hoop. Repeat flushing of hoop if catheter is difficult to remove. Inspect it prior to use to verify that it is undamaged

The hydrophilic infusion catheter must remain hydrated in order to maintain lubricity.

INSTRUCTIONS FOR USE

- The EmboCath Plus Infusion Microcatheter should be utilized coaxially over a steerable guidewire. Follow the guidewire manufacturer's instructions for use. The catheter may be used with a guide catheter/sheath with a minimum internal diameter of 0.041 in. (1.04mm). A rotating hemostatic valve used in conjunction with the guiding catheter will provide a fluid-tight seal around the infusion catheter. Introduce the infusion catheter and wire assembly through the rotating hemostatic valve. Carefully tighten the valve around the hydrophilic infusion catheter to prevent backflow, allowing some movement through the valve by the hydrophilic infusion catheter.

Overtightening of a hemostatic valve may cause damage to the catheter.

- Advance the guidewire and hydrophilic infusion catheter to a selected vascular site by alternately advancing the guidewire and then tracking the hydrophilic infusion catheter over the guidewire.

Note: To facilitate catheter handling, the proximal portion of the hydrophilic infusion catheter does not contain the hydrophilic coating. Greater resistance may be encountered when this portion of the hydrophilic infusion catheter is advanced into the rotating hemostatic valve.

- To infuse, completely remove the guidewire from the hydrophilic infusion catheter. Connect a syringe with infusate to the hydrophilic infusion catheter hub and infuse as required. **Do not inject any fluid if air bubbles are visible.**

Note: Do not exceed Flow Rates/PSI as described in the Flow Rate Table.

- Aspirate infusion catheter upon removal and discard.

CAUTION:

Federal (USA) law restricts this device to use by or on the order of a licensed physician.

CONTRAINDICATIONS:

There are no known contraindications for the use of infusion microcatheters.

STORAGE:

Store the EmboCath Plus Infusion Microcatheter in a cool, dark, dry area.

COMPATIBILITY:

- The EmboCath Plus Infusion Microcatheter should be used with steerable guidewires that are up to 0.025 inches in diameter.
- The EmboCath Plus Infusion Microcatheter is compatible with microspheres of size range 700-900 µm and smaller.
- Catheters exhibited no detectable signs of material degradation during bench testing when exposed

between 60 and 75 minutes to a mixture of the following agents: 5mg mytomicin, 25mg doxorubicin, 50mg cisplatin, 0.25g fluorouracil, 5mL ethiodol, 5cc sterile contrast and 1cc sterile water.


WARRANTY

BioSphere Medical® warrants that reasonable care has been used in the design and manufacture of this instrument. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness. Handling, storage, cleaning and sterilization of this instrument, as well as factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond BioSphere Medical's control, directly affect the instrument and the results obtained from its use. BioSphere Medical's obligation under this warranty is limited to the repair or replacement of this instrument. BioSphere Medical is not liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this instrument. BioSphere Medical neither assumes nor authorizes any other person to assume for it any other or additional liability or responsibility in connection with this instrument.

BioSphere Medical assumes no liability with respect to instruments reused, reprocessed or re-sterilized, and makes no warranties, expressed or implied, including, but not limited to, merchantability or fitness for intended use with respect to such instrument.

All serious or life threatening adverse events or deaths associated with use of EmboCath Plus should be reported to the Competent Authority of the country where it occurred and to the device manufacturer.

Information on packaging:

Symbol	Designation
	Manufacturer: Name & Address
	Use by date: year-month
	Batch code
	Catalogue number
	Do not resterilize
	Do not use if package is damaged
	Keep away from sunlight
	Keep dry
	Do not re-use
	Caution - Refer to Instructions For Use
	Sterilized using ethylene oxide
	EC mark logo - Notified body identification : 0459
	Authorized Representative in European Community
	Maximum Guidewire
	Maximum pressure
	Radiopaque marker

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