SAFEGUARD RADIAL™

INSTRUCTIONS FOR USE
MODE D'EMPLOI
ISTRUZIONI PER L'USO
GEBRAUCHSANWEISUNG
INSTRUCCIONES DE USO
INSTRUÇÕES DE UTILIZAÇÃO
INSTRUÇÕES DE USO
GEBRUIKSAANWIJZING
BRUKSANVISNING
BRUGSANVISNING
ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
УКАЗАНИЯ ПО ПРИМЕНЕНИЮ
INSTRUKCJA UŻYTKOWANIA
NÁVOD K POUŽITÍ
ИНСТРУКЦИИ ЗА УПОТРЕБА
HASZNÁLATI UTASÍTÁS

 **MERTMEDICAL™**

SAFEGUARD RADIAL™

COMPRESSION DEVICE

INSTRUCTIONS FOR USE

PRODUCT DESCRIPTION

Merit Medical's SAFEGUARD Radial™ compression device is a sterile, single use disposable device. It has a clear medical grade polyurethane window and bulb (balloon) that facilitates visualization of the puncture site, a clear medical grade PVC flexible tube, and a pressure sensitive, self-adhesive peel backing. A valve on the end of the fill tube enables a standard luer or proprietary syringe to be connected to inflate and deflate the bulb with air to provide compression of the transradial puncture site. Product not made with natural rubber latex.

INDICATIONS FOR USE

SAFEGUARD Radial is a compression device to assist hemostasis of the radial artery after a transradial procedure.

CONTRAINDICATIONS

The adhesive portion of SAFEGUARD Radial should not be used over excoriated skin.

WARNINGS

- Prior to inflation of bulb, confirm that you are ONLY injecting air into SAFEGUARD Radial compression device and NOT to side port of sheath or other device.
- Over-inflation of bulb (above 7 mL of air) could cause bulb to excessively expand and cause pain, numbness, or radial artery occlusion.
- Under-inflation of bulb could compromise ability of device to assist hemostasis of radial artery.
- Maintain sterile field during application.

CAUTIONS

- Read instructions prior to use.
- RX Only: Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
- This device is intended for single use only. Do not reuse or resterilize.
- Do not use if package is damaged
- This device should be used by clinicians with adequate training in the use of the device.

REUSE PRECAUTION STATEMENT

For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

POTENTIAL COMPLICATIONS

Possible complications that may result from use of this device include, but are not limited to: hematoma; bleeding; pain or numbness; radial artery occlusion.

INSTRUCTIONS FOR USE

1. Ensure site is clean and dry.
2. After procedure, aspirate sheath, then withdraw sheath one inch (2-3 cm). (Figure 1)
3. Peel back adhesive backing of SAFEGUARD Radial compression device approximately halfway (both sides), then apply center of bulb over arteriotomy site. (Figure 2)

Note: The arrows on SAFEGUARD Radial should be pointing up the arm of the patient.

4. Remove remainder of adhesive backing (one side at a time) and completely secure around wrist. (Figure 3)
5. Attach and completely engage a standard luer or proprietary syringe to the tubing line labeled "AIR".

WARNING: Prior to inflation of bulb, confirm that you are ONLY injecting air into SAFEGUARD Radial compression device and NOT to side port of sheath or other device.

6. Inflate the bulb (balloon) with a maximum volume of 7 mL of air. (Figure 4) Remove syringe.








WARNING: Over-inflation of bulb above 7 mL of air could cause bulb to excessively expand and cause pain, numbness, or radial artery occlusion.

WARNING: Under-inflation of bulb could compromise ability of device to assist hemostasis of radial artery.

7. Remove sheath, then confirm there is no bleeding from puncture site by viewing site through bulb window. If bleeding is observed at any time, inject more air (not exceeding the maximum inflate volume) until bleeding stops. (Figure 5)
8. Per hospital protocol, ensure adequate distal perfusion is maintained. If necessary, adjust air volume in bulb. (Figure 6)

Note: Air volume and compression time may differ according to patient's condition, anticoagulant dosage, and size of puncture site.

9. Before removing device, confirm that bleeding has stopped. Attach syringe to tubing (holding plunger in place). Slowly deflate bulb. (Figure 7)
- Note:** Be sure not to create a vacuum by pulling back on the plunger too fast.
10. Carefully remove SAFEGUARD Radial compression device and apply sterile dressing per hospital protocol.

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Authorized Representative: Merit Medical Ireland Ltd, Parkmore Business Park West, Galway, Ireland