
Inflation Device and Fluid Dispensing Syringe

INSTRUCTIONS FOR USE

U.S. Patent Nos. 5,047,015 5,057,078 5,135,488
Other U.S. and Foreign Patents Pending.

DESCRIPTION:

The Monarch™ Inflation Device and Fluid Dispensing Syringe by Merit Medical is a 20ml disposable device with an integral pressure transducer, microcomputer, back-lit LCD, threaded plunger assembly with lock/release bar, a flexible high pressure extension tube, and a three-way medium pressure stopcock. The Monarch is designed to generate and monitor pressures over a range of -1 to +30.0 ATM/BAR (-7 to 441 PSI). The Monarch syringe dispenses .45ml of fluid ± .07ml for each 360° turn of the syringe plunger handle.

CAUTION: Federal (U.S.A.) law restricts this device to use by or on the order of a physician.

INSTRUCTIONS FOR USE:

Before use, inspect the device and packaging to verify that no damage has occurred as a result of shipping.

Press the green button next to the LCD display to activate the microprocessor. After three seconds, the LCD display and backlight will come on. The syringe will display 0.0 in the pressure area of the display and the time will begin incrementing.

The syringe will be set in the ATM/BAR mode when initially turned on. To change the pressure display to read in PSI, simply press the blue button. To change back to ATM/BAR, press the blue button once again.

After an inflation or pressure monitored injection has been made, pressing the green button will display last inflation information and an "←" indicator on the display.

CAUTION: If "ER" is displayed in the pressure area and a number appears in the time window, the syringe is defective. Please return the syringe to Merit Medical for credit.

NOTE: To conserve power the backlight will automatically turn off after ten minutes of inactivity. Pressing the green button or inflating the balloon will reactivate the backlight. The syringe's microprocessor will turn off after 90 minutes of inactivity. Pressing the green button will reactivate normal operation. After ten minutes at a constant pressure the backlight will turn off, however, the microprocessor will continue to monitor the pressure. Pressing the green button will reactivate the back light.

PRIMING WITH PRIMELOK™

The PrimeLok allows purging of air and excess fluid without squeezing the trigger. The Monarch is packaged with the PrimeLok in the disengaged position.

1. To engage PrimeLok, squeeze trigger and slide PrimeLok into slot.




2. To prep syringe, simply aspirate up to 20ml of contrast solution or fluid to be dispensed into the inflation syringe by pulling back on the plunger handle.

After an inflation or pressure monitored injection has been made, pressing the green button will display last inflation information and an "←" indicator on the display.

CAUTION: If "ER" is displayed in the pressure area and a number appears in the time window, the syringe is defective. Please return the syringe to Merit Medical for credit.

NOTE: To conserve power the backlight will automatically turn off after ten minutes of inactivity. Pressing the green button or inflating the balloon will reactivate the backlight. The syringe's microprocessor will turn off after 90 minutes of inactivity. Pressing the green button will reactivate normal operation. After ten minutes at a constant pressure the backlight will turn off, however, the microprocessor will continue to monitor the pressure. Pressing the green button will reactivate the back light.

CAUTION: Federal (U.S.A.) law restricts this device to use by or on the order of a physician.

INFLATION DEVICE INDICATIONS AND USAGE:

This inflation device is used to inflate and deflate an angioplasty balloon or other interventional devices, and to measure the pressure within the balloon.

ATTACHING THE INFLATION DEVICE TO THE BALLOON:

NOTE: Refer to the manufacturer's directions accompanying the balloon dilatation catheter or other interventional device for specific information on use, maximum inflation pressure, precautions, and warnings for that device.

1. Prepare and test the balloon catheter according to the catheter manufacturer's directions for use.

2. Create a fluid-fluid connection between the balloon and the syringe extension tube, connect the luer connectors securely.

3. Squeeze the trigger and pull back on the plunger handle to apply a vacuum to the balloon.

BALLOON INFLATION AND DEFLATION:

1. To inflate the balloon, squeeze the trigger allowing the plunger to return to a resting position (0 ATM/BAR or PSI). Release grip on the trigger, locking the plunger into position. To increase pressure, rotate handle clockwise until the desired pressure is reached. The lock mechanism maintains the pressure. Pressures above the maximum range will be indicated with flashing numbers and an "↑" indicator on the display.

NOTE: Loss of pressure may indicate a leak in the system.

FLUID DISPENSING INDICATIONS AND USAGE:

This device is intended for use by healthcare professionals to dispense fluids to the body from the Monarch syringe and monitor the pressure of that fluid.

DISPENSING FLUIDS AND MONITORING PRESSURES USING THE MONARCH SYRINGE

1. To slowly dispense fluids to the body rotate the plunger clockwise until the desired fluid is injected. The injection pressure will be displayed on the LCD and the timer automatically starts once the device generates a positive pressure. Pressures above the maximum range of the syringe will be indicated with flashing numbers and an "↑" indicator on the display.

2. To rapidly dispense fluids to the body squeeze the trigger while pushing the plunger forward. The injection pressure will be displayed on the LCD and the timer automatically starts once the device generates a positive pressure. Pressures above the maximum range of the syringe will be indicated with flashing numbers and an "↑" indicator on the display.

3. A negative pressure may be generated by squeezing the trigger and pulling back the plunger. Release grip on the trigger to lock the plunger in the negative pressure position.

4. The timer will stop once the pressure returns to zero or less. Data associated with the last injection will be displayed when the green button is depressed and held and the pressure is zero or less.

CAUTION: This syringe is capable of generating high fluid pressures in a closed fluid system. The volume change of fluid dispensed may not be accurate due to compliance of the plastic components as pressure changes.

RECOMMENDED SEPARATION DISTANCES BETWEEN PORTABLE AND MOBILE RADIO FREQUENCY (RF) COMMUNICATIONS EQUIPMENT AND THE MONARCH™ INFLATION DEVICE

The Monarch™ Inflation Device is intended for use in an electromagnetic environment in which RF radiated disturbances are controlled. The user of the Monarch™ Inflation Device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Monarch™ Inflation Device as recommended below, according to the maximum output power of the communications equipment.

DISPENSING FLUIDS AND MONITORING PRESSURES USING THE MONARCH SYRINGE

1. To slowly dispense fluids to the body rotate the plunger clockwise until the desired fluid is injected. The injection pressure will be displayed on the LCD and the timer automatically starts once the device generates a positive pressure. Pressures above the maximum range of the syringe will be indicated with flashing numbers and an "↑" indicator on the display.

2. To rapidly dispense fluids to the body squeeze the trigger while pushing the plunger forward. The injection pressure will be displayed on the LCD and the timer automatically starts once the device generates a positive pressure. Pressures above the maximum range of the syringe will be indicated with flashing numbers and an "↑" indicator on the display.

3. A negative pressure may be generated by squeezing the trigger and pulling back the plunger. Release grip on the trigger to lock the plunger in the negative pressure position.

4. The timer will stop once the pressure returns to zero or less. Data associated with the last injection will be displayed when the green button is depressed and held and the pressure is zero or less.

CAUTION: This syringe is capable of generating high fluid pressures in a closed fluid system. The volume change of fluid dispensed may not be accurate due to compliance of the plastic components as pressure changes.

RECOMMENDED SEPARATION DISTANCES BETWEEN PORTABLE AND MOBILE RADIO FREQUENCY (RF) COMMUNICATIONS EQUIPMENT AND THE MONARCH™ INFLATION DEVICE

The Monarch™ Inflation Device is intended for use in an electromagnetic environment in which RF radiated disturbances are controlled. The user of the Monarch™ Inflation Device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Monarch™ Inflation Device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (in watts) W	Separation distance according to frequency of transmitter (in meters) m		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
0.01	$d = [1.17] \sqrt{P}$	$d = [2.33] \sqrt{P}$	$d = [9.33] \sqrt{P}$
0.1	0.12	0.23	0.93
1	0.37	0.74	2.95
10	1.17	2.33	9.33
100	3.69	7.38	29.51
	11.67	23.33	93.33

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 – At 80 MHz and 800 MHz, the separation distance for the higher frequency applies.

NOTE 2 – These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Interference may occur.
Des interférences peuvent survenir.
Potrebbero presentarsi delle interferenze.
Es könnte zu einer Interferenz kommen.
Puede darse interferencia.
Podem ocorrer interferências.
Dit kan leiden tot hinderling.
Störning kan inträffa.
Interferens kan forekomme.
Μπορεί να παρουσιαστεί παρεμβολή.

GUIDANCE AND MANUFACTURER'S DECLARATION – ELECTROMAGNETIC IMMUNITY

The EQUIPMENT is intended for use in the electromagnetic environment specified below. The customer or the user of the EQUIPMENT should assure that it is used in such an environment.

Immunity Test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge (ESD) IEC 61000-4-20	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Increasing relative humidity will reduce the potential for ESD related difficulties

Users should follow local guidelines and practices regulating the disposal of infected waste products.

Rx Only - Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. / Attention: La loi fédérale américaine limite la vente de cet article à des médecins ou sur ordonnance. / Attenzione: La legge federale degli Stati Uniti limita la vendita di quest'articolo ai medici o dietro ricetta medica. / Achtung: Das Bundesgesetz der USA limitiert den Verkauf dieses Artikels an Ärzte oder gegen Rezept. / Precaución: La ley Federal de los Estados Unidos de América (EEUU) prohíbe la venta de este producto sin prescripción médica. / Atenção: A Lei Federal dos EUA restringe a venda deste dispositivo a, ou sob a prescrição de um médico. / Waarschuwing: de wetgeving van de Verenigde Staten bepaalt dat de verkoop van dit apparaat uitsluitend door of op recept van een arts kan geschieden. / Försiktighet: I USA förbjuder federala lagar att denna utrustning säljs av läkare eller på ordination av läkare. / Vigtigt: Foderale myndigheder i USA har bestemt, at denne artikel kun må sælges til læger eller på recept. / ΠΡΟΣΟΧΗ: Αμερικανικός ομοσπονδιακός νόμος επιτρέπει την πώληση αυτού του είδους μόνο σε ιατρούς ή έναντι ιατρικής συνταγής.

Manufacturer
Merit Medical Systems, Inc.
South Jordan, Utah 84095 U.S.A. 1-801-253-1600
U.S.A. Customer Service 1-800-356-3748

Authorized Representative
Merit Medical Ireland Ltd
Parkmore Business Park West, Galway, Ireland

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