OUTAKE™

CATHETER EXTRACTOR

INSTRUCTIONS FOR USE:

INTENDED USE: The Catheter Extractor is a catheter accessory device intended for minimally invasive removal of tunneled long term catheters. The Catheter Extractor is indicated for use with tunneled, double lumen hemodialysis catheters with cuff, sizes ranging from 10.8 French to 14.5 French.

WARNINGS AND PRECAUTIONS:

- It is recommended that this device be used by a physician who is familiar with the appropriate technique.
- Do not use any product if package is damaged or has been opened.
- Do not use if Extractor shows signs of damage.
- Sterile if package is unopened and undamaged.
- Single use only. Do not re-sterilize the Extractor or accessories by any method. The manufacturer will not be held liable for any damages caused by re-use or resterilization of the catheter or accessories.

RX Only: Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

INSTRUCTIONS FOR USE:

- Assemble all required supplies and bring to patient’s bedside. Ascertain if patient has Betadine allergy.
- Wash hands for at least 30 seconds with antimicrobial soap. Rinse hands and dry with sterile towel.
- Properly position the patient for the procedure. Have patient wear a cap if indicated (i.e. long hair).
- If sutures are present, they must be removed prior to procedure. Don sterile gloves. Using sterile suture scissors, clip the thread as close to the skin as possible. Lift the sutures out with sterile forceps.
- Don sterile gloves and begin preparation of the procedure site.
- Prep patient’s area around the catheter’s exit site and the exposed part of catheter with Providone-Iodine swab stick or equivalent antiseptic solution.
- The area to be prepped should extend from the intended site of catheter removal to a radius of 15cm. Cleanse with friction using a circular motion from center to periphery. Allow to dry (This provides a protective barrier against skin organisms)
- Don sterile mask, cap, and gown, and avoid contamination of site and equipment. Maintain aseptic technique throughout the procedure. Use a sterile drape as a barrier between patient’s mouth/nose/trach and insertion site.
- Inject lidocaine solution or equivalent anesthetic to the exit site and along the catheter up to the cuff.
- Snap the Extractor to the tube part of the exposed catheter.
- Insert the Extractor through the exit hole.
- Using a twisting, downward motion, break the cuff free from the fibrous tissue while pulling the catheter out.
- Pull catheter out of the exit hole slowly. Assure the whole catheter including its tip and cuff are removed without any parts remaining in the tissue.
- Hold pressure at the initial vein access site of either the Subclavian or Jugular veins and also at the exit site for approximately 6-10 minutes or until bleeding has subsided.
- Place sterile pressure dressing over the exit site.
- Examine the cuff to see if any parts of it are missing.
- Discard used supplies into a designated sharps or biohazardous waste container as appropriate.
- Wash hands.

Manufacturer: www.merit.com
Merit Medical Systems, Inc. South Jordan, Utah 84095
U.S.A. 1-801-253-1600 U.S.A. Customer Service 1-800-356-3748

EC REP
Authorized Representative:
Merit Medical Ireland Ltd, Parkmore Business Park West, Galway, Ireland

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