HYDROPHILIC

# Prelude SNAP™

Splittable Sheath Introducer

**INSTRUCTIONS FOR USE  
MODE D'EMPLOI  
ISTRUZIONI PER L'USO  
GEBRAUCHSANWEISUNG  
INSTRUCCIONES DE USO  
INSTRUÇÕES DE UTILIZAÇÃO  
INSTRUÇÕES DE UTILIZAÇÃO  
GEBRUIKSAANWIJZING  
BRUKSANVISNING  
BRUGSANVISNING  
ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ  
KULLANMA TALİMATLARI  
УКАЗАНИЯ ПО ПРИМЕНЕНИЮ  
NÁVOD K POUŽITÍ  
INSTRUKCJA OBSŁUGI  
ИНСТРУКЦИИ ЗА УПОТРЕБА  
HASZNÁLATI UTASÍTÁS**

## HYDROPHILIC

English

# Prelude SNAP™

## Splittable Sheath Introducer

### INSTRUCTIONS FOR USE

**This device is intended for one time use only. Read instructions prior to use.**

**Rx Only: Caution:** Federal (USA) law restricts this device to sale by or on the order of a physician.

#### For U.S.-California Only.

Proposition 65, a State of California voter initiative, requires the following notice:

**WARNING:** This product and its packaging have been sterilized with ethylene oxide. This packaging may expose you to ethylene oxide, a chemical known to the state of California to cause cancer or birth defects or other reproductive harm.

#### Indications

For the introduction of various types of pacing leads and catheters to the heart and coronary venous system.

#### Warnings

- This product is sensitive to light. Do not use if stored outside the protective outer carton.
- Store in a cool, dark, and dry place.
- Infusion through the sideport can be done only after all air is removed from the unit.

#### Precautions

- Do not alter this device in any way.
- **Single Use Devices:** This single-use product is not designed or validated to be reused. Reuse may cause a risk of cross-contamination, affect the measurement accuracy, system performance, or cause a malfunction as a result of the product being physically damaged due to cleaning, disinfection, re-sterilization, or reuse.
- Aspiration and saline flushing of the sheath and dilator should be performed to help minimize the potential for air embolism and clot formation.
- Indwelling introducer sheaths should be internally supported by a catheter, pacing lead, or dilator.
- Dilators, catheters, and pacing leads should be removed slowly from the sheath. Rapid removal may damage the valve resulting in flow of blood or air through the valve.
- Never advance or withdraw the guide wire or sheath when resistance is met. Determine the cause by fluoroscopy and take remedial action.
- When injecting or aspirating through the sheath, use the sideport only (if applicable).
- Ensure that the surface of the sheath is wet prior to insertion; the sheath should not be used in a dry state.

#### Adverse Events

Adverse events may include, but are not limited to the following:


- Air embolus
- Blood loss
- Vessel damage
- Infection
- Pneumothorax
- Hemothorax
- Pacing lead displacement

#### Use Sterile Technique *A suggested procedure:*

1. Peel open package and place contents in sterile field.
2. Prep skin and drape in area of anticipated veni-puncture.
3. Flush the needle, sheath, and dilator prior to use.
4. Insert the dilator into the sheath until the dilator cap attaches to the sheath hub.

5. Using sterile saline or water wet the outside of the sheath.
6. Insert the needle into vessel. Verify the needle position by observing venous blood return.
7. Aspirate the needle using the syringe.
8. Remove the syringe and insert soft tip of the guide wire through the needle into the vessel. Advance the guidewire to required depth. Leave an appropriate amount of the guidewire exposed. At no time should the guide wire be advanced or withdrawn when resistance is met. Determine the cause of resistance before proceeding. Fluoroscopic verification of the guide wire's entrance into the superior vena cava and right atrium is suggested.
9. Hold the guide wire in place and remove needle. Do not withdraw the guide wire back into the needle as this may result in separation of the guide wire.
10. Thread the dilator/sheath assembly over the guide wire.
11. Advance the dilator/sheath assembly with a twisting motion over the guide wire and into the vessel. Fluoroscopic observation may be advisable. Attaching a clamp or hemostat to the proximal end of the guide wire will prevent inadvertently advancing the guide wire entirely into the patient.
12. Once assembly is fully introduced into the venous system, remove the dilator cap from the sheath by rocking the dilator cap off the sheath hub.
13. Slowly retract the guide wire and dilator, leaving the sheath in position. The sheath will reduce the loss of blood and the inadvertent aspiration of air.
14. Aspirate all air from the sheath by using a syringe connected to the sideport (if applicable).
15. Flush the introducer with saline through the sideport (if applicable). If the introducer is to remain in place during lead positioning and testing, periodic flushing is advised.
16. Introduce pacemaker lead or catheter through the sheath and advance it into position.
17. Flush sheath with saline immediately before peeling sheath away in order to minimize backbleeding.
18. Split sheath by sharply snapping the sheath hub tabs and peeling sheath tube apart while withdrawing it from the vessel.

#### Graphic Symbols for Medical Device Labeling

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Manufacturer: Merit Medical Systems, Inc.
1600 W. Merit Pkwy, So. Jordan, Utah 84095 U.S.A.
U.S.A. Customer Service 1-800-356-3748
Authorized Representative: Merit Medical Ireland Ltd,