ANGIOGRAPHIC CATHETERS

INSTRUCTIONS FOR USE

Description: Angiographic catheters are single-use, sterile, non-pyrogenic, disposable intravascular catheters.

Indications: Angiographic catheters are designed to be used for delivering radiopaque media to selected sites in the vascular system in conjunction with routine diagnostic procedures. Angiographic catheters with marker bands may also be used for anatomical measurements.

Contraindications: None known.

Warnings:

- The Impress® angiographic catheter with hydrophilic coating should not be used in the coronary arteries
- Angiographic catheters are single use.
- Recommended guide wire size: See the individual catheter label.
- Guide wire use has been associated with greater incidence of thrombus formation.

Precautions:

- Do not use if package is opened or damaged.
- Do not use after expiration date.
- Do not autoclave.
- Carefully remove the catheter from the package by grasping the hub and slowly withdrawing it from the package.
- If resistance is felt when removing the guidewire from the catheter, remove the guide wire and catheter as a unit to prevent potential damage to the vessel wall.
- Flush all devices entering a blood vessel with sterile heparinized saline or similar isotonic solution before use.
- For catheters containing hydrophilic polymer coating on the distal surface of the catheter, the surface of the distal portion of the catheter must be completely wet with solution to remain lubricated. Proper caution must be taken to wet the catheter prior to and during use.

Adverse Events:

Merit Medical Angiographic Catheters should be used by physicians familiar with percutaneous catheter introduction. Complications which may be associated with the use of angiographic catheters include, but are not limited to, the following:

- Thrombus formation/emboli
- Hematoma at the puncture site
- Plaque dislodgement
- Cardiac arrhythmias
- Perforation of the vessel wall
- Vascular occlusion
- Stroke and death
- Air embolism
- Arterial wall damage
- Myocardial infarction
- Infection
- Allergic reaction to contrast media
- Vascular spasm

Directions for Use:

- A coated guide wire is recommended for use with Merit Angiographic Catheters.
- Introduce the catheter into the blood vessel using vascular entry technique of choice.
- PIGTAIL catheters. Pigtail catheters are fitted with a straightener. It is split for easy removal after the catheter tip is on the guide wire.
- Slide straightener forward until pigtail is straight.
- Place catheter tip on guide wire and remove straightener by pulling on the straightener tab.
- MARKER BANDS. Anatomical measurement accuracy is ± 4%.

Sterile and non-pyrogenic in unopened and undamaged inner package. Store in a cool dry place.

REUSE PRECAUTION STATEMENT

For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

Manufacturer

Merit Medical Systems, Inc. South Jordan, Utah 84095 U.S.A. 1-801-253-1600 U.S.A. Customer Service 1-800-356-3748

EC REP

Authorized Representative

Merit Medical Ireland Ltd, Parkmore Business Park West, Galway, Ireland European Customer Service by Country: Belgium 0800 72906; France 0800 916030; Germany 0800 1820871; Ireland 091 703700; Neth. 0800 0228184; U.K. 0800 973115

F FRENCH SIZE / Taille en French / Diametro esterno in French / Größe in French / Tamaño en French / Dimensão em “French” / Afmeting in French / Storlek i French / Storrelse i tommer / Μέγεθος σε French

PSI RECOM.PRESSURE LIMIT / Pression maximale recommandée / Limite pressione massima / Empt.Drucklimit / Presión máxima recomendada / Pressão maxima recomendada / Aanbevolen maximale druk / tryck-gräns / Anbef.tryklimit i psi / Προτεινόμενο όριο πίεσης

GW RECOM.GUIDEWIRE / Guide recommandé / Filo guida consigliato / Empf.Führungsdraht / Guia Recomendada / Fio Guia Recomendado / Aanbevolen voerdraad / Rekommederad ledäre / Anbef.guidewire / Προτεινόμενο οδηγό σύρμα

CT CONTENTS / Contenu / Contenido / Inhalt / Contenido / Conteúdo / Inhoud / Innehåll / Indhold / Περιεχόμενα

SH SIDE HOLES / Perforations laterals / Fori laterali / Seitenöffnungen / Orificios laterals / Orificios laterais / Zijgaten / Sidohål / Sidehuller / Πλευρικές οπές

FLO FLOW RATES and PRESSURES. FLOW RATES ACHIEVED AT10.4cP. INDIVIDUAL RESULTS MAY VARY. / Flux et pressions. Flux obtenus à 10,4cP. Les résultats individuels peuvent varier. / Velocità di flusso o pressioni. Velocità di flusso raggiunte a 10,4cP. I risultati individuali potrebbero variare. / Flussraten und Drücke. Flussraten erreicht bei 10,4cP.Einzelne Ergebnisse können variieren. / Caudal de flujo y presión. El caudal de flujo llegó e los 10,4cP. Los resultados individuales pueden variar. / Velocidades de flujo e pressões. Velocidades de fluxo obtidas a 10,4cP. Os resultados individuais podem variar. / StrömingsSNnelheid en druk. De stromingsnelheid komt tot stand bij 10,4cP. Afzonderlijke resultaten kunnen verschillend zijn. / Tillflödes-mätt och tryck. Uppnått tillflödesmätt vid 10,4cP. Individuella resultat kan variera. / Gennemstrømningshastigheder og tryk. Gennemstrømningshastigheder opnaet ved 10,4cP. Individuelle resultater kan variere. / Ταχύτητες ροής και Πίεσεις. Οι ταχύτητες ροής έφθασαν στα 10,4cP. Οι τιμές που επιτυγχάνει το κάθε άτομο μπορεί να ποικίλλουν.

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