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**CAUTION:** Federal (U.S.A.) law restricts this device to use by or on the order of a licensed physician.

#### INTENDED USE:

- Embolization of hypervascularized tumors.
- Embolization of arteriovenous malformations (AVM).

#### CLINICAL APPLICATIONS:

Scientific literature provides extensive documentation of embolization procedures using a wide variety of agents in both the neurological and peripheral circulatory systems, including the head, neck, spine, liver, genitourinary tract, uterus, gastrointestinal system, limbs, and lungs. A representative bibliography of this literature is provided following these instructions for use.

#### DESCRIPTION:

EmboGold Microspheres are part of a family of new embolic materials based on Merit Medical's proprietary microsphere technology. These colored spheres are designed to add easy visualization during handling, while still offering the controlled, targeted embolization characteristics of the first microsphere product, Embosphere® Microspheres.

EmboGold Microspheres are biocompatible, hydrophilic, non-resorbable, precisely calibrated microspheres produced from an acrylic polymer and impregnated with porcine gelatin. EmboGold Microspheres are available in a range of sizes.

#### DEVICE PACKAGING:

- Embolic contained in sterile, 20 mL syringe (pre-filled syringe), packaged in a peel-away tray.
- Each syringe contains approximately 1.0 mL or 2.0 mL of EmboGold Microspheres in pyrogen-free, sterile, physiological saline.

#### CONTRAINDICATIONS:

- Patients intolerant to occlusion procedures
- Vascular anatomy or blood flow that precludes catheter placement or embolic agent injection.
- Presence or likely onset of vasospasm.
- Presence or likely onset of hemorrhage.
- Presence of severe atherosomatous disease.
- Presence of feeding arteries smaller than distal branches from which they emerge.
- Presence of patent extra-to-intracranial anastomoses or shunts.
- Presence of collateral vessel pathways potentially endangering normal territories during embolization.
- Presence of end arteries leading directly to cranial nerves.

- Presence of arteries supplying the lesion not large enough to accept EmboGold Microspheres.
- Vascular resistance peripheral to the feeding arteries precluding passage of EmboGold Microspheres into the lesion.
- Do not use EmboGold Microspheres in the following applications:
  - a. Embolization of large diameter arteriovenous shunts (i.e. where the blood does not pass through an arterial/capillary/venous transition but directly from an artery to a vein).
  - b. In the pulmonary arterial vasculature.
  - c. Use in any vasculature where EmboGold Microspheres could pass directly into the internal carotid artery, vertebral artery, intracranial vasculature or the above-listed vessels.
- Patients who have an allergic response to gold.

#### WARNINGS:

- EmboGold Microspheres contain gelatin of porcine origin, and therefore, could cause an immune reaction in patients who are hypersensitive to collagen or gelatin. Careful consideration should be given prior to using this product in patients who are suspected to be allergic to injections containing gelatin stabilizers.
- Studies have shown that EmboGold Microspheres do not form aggregates, and as a result, penetrate deeper into the vasculature as compared to similarly sized PVA particles. Care must be taken to choose larger sized EmboGold Microspheres when embolizing arteriovenous malformations with large shunts to avoid passage of the spheres into the pulmonary or coronary circulation.
- Some of the EmboGold Microspheres may be slightly outside of the range, so the physician should be sure to carefully select the size of EmboGold Microspheres according to the size of the target vessels at the desired level of occlusion in the vasculature and after consideration of the arteriovenous angiographic appearance. EmboGold Microspheres size should be selected to prevent passage from artery to vein.
- Because of the significant complications of misembolization, extreme caution should be used for any procedures involving the extracranial circulation encompassing the head and neck, and the physician should carefully weigh the potential benefits of using embolization against the risks and potential complications of the procedure. These complications can include blindness, hearing loss, loss of smell, paralysis and death.
- The color of the EmboGold Microspheres could be visible through the skin if injected into arteries feeding superficial tissues.
- Serious radiation-induced skin injury may occur to the patient due to long periods of fluoroscopic exposure, large patient diameter, angled x-ray projections, and multiple image recording runs or radiographs. Refer to your facility's clinical protocol to ensure the proper radiation dose is applied for each specific type of procedure performed. Physicians should monitor patients that may be at risk.
- Onset of radiation-induced injury to the patient may be delayed. Patients should be counseled on potential radiation side effects and whom they should contact if they show symptoms.
- EmboGold Microspheres are not cleared for use in uterine fibroid embolization and are associated with delayed pain and/or rash when used in these procedures.

- Pay careful attention for signs of mistargeted embolization. During injection carefully monitor patient vital signs to include SAO2 (e.g. hypoxia, CNS changes). Consider terminating the procedure, investigating for possible shunting, or increasing microsphere size if any signs of mistargeting occur or patient symptoms develop.
- Consider upsizing the microspheres if angiographic evidence of embolization does not quickly appear evident during injection of the microspheres.

#### **Warnings about use of small microspheres**

- Careful consideration should be given whenever use is contemplated of embolic agents that are smaller in diameter than the resolution capability of your imaging equipment. The presence of arteriovenous anastomoses, branch vessels leading away from the target area or emergent vessels not evident prior to embolization can lead to mistargeted embolization and severe complications.
- Microspheres smaller than 100 microns will generally migrate distal to anastomotic feeders and therefore are more likely to terminate circulation to distal tissue. Greater potential ischemic injury results from use of smaller sized microspheres and consideration must be given to the consequence of this injury prior to embolization. The potential consequences include: swelling, necrosis, paralysis, abscess and/or stronger post embolization syndrome.
- Post-embolization swelling may result in ischemia to tissue adjacent to target area. Care must be given to avoid ischemia-intolerant, nontargeted tissue such as nervous tissue.

#### **PRECAUTION:**

- Do not use if the syringe, plunger seal or tray package appears damaged.
- For single patient use only - Contents supplied sterile - Never reuse, reprocess, or resterilize the contents of a syringe that has been opened. Reusing, reprocessing or resterilizing may compromise the structural integrity of the device and/or lead to device failure, which in turn may result in patient injury, illness or death. Reusing, reprocessing or resterilizing may also create a risk of contamination of the device and/or cause patient infection or cross infection including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. All procedures must be performed according to accepted aseptic technique.
- Do not connect the 20 mL syringe with EmboGold Microspheres directly to a microcatheter for embolic delivery, as a catheter occlusion may result.
- The syringe is intended for embolic use only. Do not use for another application.
- Select the size and quantity of EmboGold Microspheres appropriate for the pathology to be treated.
- Embolization with EmboGold Microspheres should only be performed by physicians who have received appropriate interventional occlusion training in the region intended to be embolized.
- Patients with known allergy to contrast medium may require corticosteroids prior to embolization.
- Additional evaluations or precautions may be necessary in managing perioperative care for patients with the following conditions:
  - Bleeding diathesis or hypercoagulative state

- Immunocompromise.

#### **POTENTIAL COMPLICATIONS:**

Vascular embolization is a high-risk procedure. Complications may occur at any time during or after the procedure, and may include, but are not limited to, the following:

- Paralysis resulting from untargeted embolization or ischemic injury from adjacent tissue edema.
- Undesirable reflux or passage of EmboGold Microspheres into normal arteries adjacent to the targeted lesion or through the lesion into other arteries or arterial beds, such as the internal carotid artery, pulmonary, or coronary circulations.
- Pulmonary Embolism due to arterial venous shunting.
- Ischemia at an undesirable location, including ischemic stroke, ischemic infarction (including myocardial infarction), and tissue necrosis.
- Capillary bed occlusion and tissue damage.
- Vessel or lesion rupture and hemorrhage.
- Neurological deficits, including cranial nerve palsies.
- Vasospasm.
- Death.
- Recanalization.
- Foreign body reactions necessitating medical intervention.
- Infection necessitating medical intervention.
- Clot formation at the tip of the catheter and subsequent dislodgment.
- Allergic response to gold.
- Complications related to catheterization (e.g., hematoma at the site of entry, clot formation at the tip of the catheter and subsequent dislodgment, and nerve and/or circulatory injuries, which may result in leg injury).
- Allergic reaction to medications (e.g. analgesics).
- Allergic reaction to contrast media or embolic material.
- Pain and/or rash, possibly delayed from the time of embolization.
- Blindness, hearing loss, loss of smell.
- Additional information is found in the Warnings section.

#### **STORAGE AND STERILITY:**

Inspect packaging prior to use to ensure seal integrity for maintenance of sterility.

- EmboGold Microspheres must be stored in a cool, dry and dark place in their original syringe and packaging.
- Use by the date indicated on the syringe label.
- Do not freeze.
- Do not resterilize.

#### **INSTRUCTIONS FOR USE:**

- Carefully evaluate the vascular network associated with the lesion using high resolution imaging prior to beginning the embolization procedure.
- EmboGold Microspheres are available in a range of sizes. Because of the potential for misembolization and the inherent variability in sphere sizes, the physician should be sure to carefully select the size of EmboGold Microspheres according to the size of the target vessels at the desired level of occlusion in the vasculature.
- When embolizing arteriovenous malformations (AVM), choose an

EmboGold Microspheres size that will occlude the nidus without passing through the AVM.

- Choose a delivery catheter based on the size of the target vessel and the microsphere size being used. EmboGold Microspheres can tolerate temporary compression of up to 33% in order to facilitate passage through the delivery catheter.
- EmboGold Microspheres are not radiopaque. It is recommended that the embolization be monitored using fluoroscopic visualization by adding the appropriate amount of contrast medium to the physiologic saline suspension fluid.
- To deliver EmboGold Microspheres:  
Match the total volume in the syringe with the same volume of undiluted contrast, which will result in a 50% microsphere/saline and 50% contrast solution. Remove all air from the syringe. To evenly suspend the EmboGold Microsphere/contrast solution, gently invert the 20mL syringe several times. Attach the 20mL syringe to one port of the luer-lock 3-way stopcock; attach a 1mL or 3mL injection syringe to another port on the stopcock, and, if desired, a delivery catheter may be attached to the remaining port on the stopcock. Wait several minutes to allow the EmboGold Microspheres to suspend properly. Draw EmboGold Microsphere/contrast solution into the injection syringe slowly and gently to minimize the potential of introducing air into the system. Purge all air from the system prior to injection. Inject the EmboGold Microsphere/ contrast solution under fluoroscopic visualization with the injection syringe, using a slow pulsatile injection, while observing the contrast flow rate. If there is no effect on the flow rate, repeat the delivery process with additional injections of EmboGold Microsphere/ contrast solution, or larger sized EmboGold Microspheres may be considered. If the EmboGold Microsphere/contrast solution requires resuspension, gently invert the 20mL syringe several times. Exercise conservative judgment in determining the embolization endpoint.
- Femoral puncture can result in arterial spasm. This may predispose to femoral thrombosis (e.g. leg injury). Femoral patency should be re-assessed prior to final catheter removal.
- Upon completion of the treatment, remove the catheter while maintaining gentle suction so as not to dislodge EmboGold Microspheres still within the catheter lumen.
- Apply pressure to the puncture site until hemostasis is complete.
- Discard any open, unused EmboGold Microspheres.

## REFERENCES:

- Bendszus M, Klein R, Burger R, et al: Efficacy of trisacryl gelatin microspheres versus polyvinyl alcohol particles in the preoperative embolization of meningiomas. AJNR, 21(2):255-61, Feb 2000.
- Deveikis JP: Endovascular therapy of intracranial arteriovenous malformations: materials and techniques. Neuroimaging Clin of N Am, 8(2):401-424, 1998.
- Frizzel RT, Fisher WS: Cure, morbidity, and mortality associated with embolization of brain arteriovenous malformations: A review of 1246 patients in 32 series over a 35-year period. Neurosurg, 37(6):1031-1040, Dec 1995.
- Gomes, A: Embolization therapy of congenital arteriovenous malformations; Use of alternative approaches. Radiology, 190:191-8, Jan 1994.
- Terada, T; Kinoshita, Y; Yokote, H; Tsuura, M; Itakura, T; Komai, N; Nakamura, Y; Tanaka, S; Kuriyama, T: Preoperative embolization of meningiomas fed by ophthalmic branch arteries. Surg Neurol, 45:161-6, 1996.

## Information on packaging:

Symbol	Designation
	Manufacturer: Name & Address
	Use by date: year-month
	Batch code
	Catalogue number
	Do not resterilise
	Do not use if package is damaged
	Keep away from sunlight
	Keep dry
	Do not re-use
	Caution - Refer to Instructions For Use
	Non-pyrogenic
	Sterilised using steam
	Lower limit of temperature

## Distributed in USA by:

Merit Medical Systems, Inc.

Customer Service                            1-800-356-3748

All serious or life-threatening adverse events or deaths associated with use of EmboGold Microspheres should be reported to the U.S. Food and Drug Administration under the MedWatch program and to the device manufacturer. Information about the MedWatch program and forms for reporting adverse events can be obtained at [www.fda.gov/safety/medwatch/howtoreport/ucm053074.htm](http://www.fda.gov/safety/medwatch/howtoreport/ucm053074.htm) or by calling toll free 888-463-6332. Reports to Merit Medical, Inc. can be made by calling toll free 800-394-0295.

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## RÉFÉRENCES :

1. Bendszus M, Klein R, Burger R, et al.: Efficacy of trisacryl gelatin microspheres versus polyvinyl alcohol particles in the preoperative embolization of meningiomas. AJNR, 21(2):255-61, Feb 2000.
2. Deveikis JP: Endovascular therapy of intracranial arteriovenous malformations: materials and techniques. Neuroimaging Clin N Am, 8(2):401-424, 1998.
3. Frizzel RT, Fisher WS: Cure, morbidity, and mortality associated with embolization of brain arteriovenous malformations: A review of 1246 patients in 32 series over a 35-year period. Neurosurg, 37(6):1031-1040, Dec 1995.
4. Gomes, A: Embolization therapy of congenital arteriovenous malformations; Use of alternative approaches. Radiology, 190:191-8, Jan 1994.
5. Terada, T; Kinoshita, Y; Yokote, H; Tsuura, M; Itakura, T; Komai, N; Nakamura, Y; Tanaka, S; Kuriyama, T: Preoperative embolization of meningiomas fed by ophthalmic branch arteries. Surg Neurol, 45:161-6, 1996.

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## BIBLIOGRAFIA:

- Bendszus M, Klein R, Burger R, et al.: Efficacy of trisacryl gelatin microspheres versus polyvinyl alcohol particles in the preoperative embolization of meningiomas. AJNR, 21(2):255-61, Feb 2000.
- Deveikis JP: Endovascular therapy of intracranial arteriovenous malformations: materials and techniques. Neuroimaging Clin N Am, 8(2):401-424, 1998.
- Frizzel RT, Fisher WS: Cure, morbidity, and mortality associated with embolization of brain arteriovenous malformations: A review of 1246 patients in 32 series over a 35-year period. Neurosurg, 37(6):1031-1040, Dec 1995.
- Gomes, A: Embolization therapy of congenital arteriovenous malformations; Use of alternative approaches. Radiology, 190:191-8, Jan 1994.
- Terada, T; Kinoshita, Y; Yokote, H; Tsuura, M; Itakura, T; Komai, N; Nakamura, Y; Tanaka, S; Kuriyama, T: Preoperative embolization of meningiomas fed by ophthalmic branch arteries. Surg Neurol, 45:161-6, 1996.

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**Biosphere Medical, S.A.**
Utah 84095 U.S.A. 1-801-253-1600
U.S.A. Customer Service 1-800-356-3758