


- To ensure fixation, suture at each eyelet on the wings of the Secondary Fixation Device.
- To open the device, hold the lock device down and pull the plate up, using the lever.

VENAGUIDE™

Note: This device is found in certain models only.



- Release guidewire by removing the guidewire cap.
- Straighten guidewire 'J' by retracting into introducer system with thumb.
- Insert into hub of introducer needle and advance guidewire into vein. Advance to required depth.

WARNING: The included guidewire does not have depth markings. Electrocardiogram (ECG) monitoring or Ultrasound guidance, or a combination of both is recommended to prevent under insertion/ over insertion of the guidewire.

GUIDING SYRINGE

Note: This accessory is found in certain models only.

This device is used for introducing guidewires for use with Merit central or peripheral venous catheters.

- Insert introducer needle attached to Guiding Syringe into vessel and aspirate.
- Feed the desired flexible end of the guidewire through the hole in rear of Guiding Syringe plunger into the vein. If a 'J' wire is to be used, straighten guidewire 'J' by retracting into introducer system with thumb.
- Hold the guidewire in place and remove introducer needle and Guiding Syringe.

WARNING: Do not aspirate with guidewire in place or air may enter the syringe

Caution: Do not reinfuse blood to minimize the risk of blood leakage from rear of syringe.

FLOSWITCH™ INTRODUCER

This device is used for introducing guidewires for use with Merit central or peripheral venous catheters.

Using a Floswitch™ introducer cannula, puncture the vein, aspirate and advance the cannula into the vein and remove the needle. At this stage the Floswitch™ may be switched off. The switch must be opened again to facilitate guidewire or catheter passage.

FLOSWITCH™ LUER LOCK ATTACHMENT

This device may be attached to a luer hub of a central or peripheral venous catheter.

- As standard practice the catheter should be secured to the skin.
- To attach the Floswitch™, insert the Floswitch™ luer into the catheter hub. Tighten the rotating collar ensuring that 'hand tight' connection has been made and that the catheter hub is fully advanced into the Floswitch™ collar (see figure 1.1).
- The Floswitch™ should be fixed in position by either suturing or taping the wings. The recommended procedure for taping is shown in figure 2.
- The black markings indicate flow status. When visible the catheter is open and there is free flow. When covered, the catheter is closed.

Figure 1. Floswitch™ Inspection Procedure





Figure 2. Floswitch Taping Procedure



WARNINGS - FLOSWITCH™ PRODUCTS

- The Floswitch™ should not be switched off (black marks covered) before the needle has been fully withdrawn. This applies to Floswitch™ introducers only.
- The Floswitch™ must not be switched off (black marks covered) until the guidewire has been fully withdrawn.
- When the catheter is not in use for infusion or aspiration, the Floswitch™ must be switched off and a suitable luer cap locked into the hub.
- When using an intermittent injection bung attached to a Floswitch™, only use short needles. Ensure Floswitch™ is in 'ON' position prior to injecting or aspirating. Do not insert the needle more than 8 mm into the Floswitch™.
- As standard practice, the security of the luer connection must be checked routinely. This is essential when lubricious substances such as Intralipids are being used. This applies to Floswitch™ luer lock attachment only.

STORAGE CONDITIONS

Store in cool dry place away from direct sunlight.

STERILE AND NON-PYROGENIC in unopened, undamaged package. For single use only. Check integrity of the individual package before use. Dispose of product after use. Do not resterilize.

Re-use may lead to infection or other illness/ injury.

For reordering information or assistance please contact local representative.

	Do Not Use If Package is Damaged
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Do not Resterilize
	Single Use
	Does not Contain DEHP, DIBP, DBP, BBP
	Non-Pyrogenic
	Sterilized using Ethylene Oxide
	Caution: Consult accompanying documents.



www.merit.com



Manufacturer:
Merit Medical Singapore Pte. Ltd.,
198 Yishun Avenue 7,
Singapore 768926



Authorized Representative:
Merit Medical Ireland Ltd
Parkmore Business Park West, Galway, Ireland
EC Customer Service +31 43 3588222