ONE Snare®

Endovascular Snare System Standard Snare Kit

EN	INSTRUCTIONS FOR USE
FR	MODE D'EMPLOI
IT	ISTRUZIONI PER L'USO
DE	GEBRAUCHSANWEISUNG
ES	INSTRUCCIONES DE USO
PT	INSTRUÇÕES DE UTILIZAÇÃO
NL	GEBRUIKSAANWIJZING
SV	BRUKSANVISNING
NO	BRUKSANVISNING
DA	BRUGSANVISNING
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
TR	KULLANMA TALIMATLARI
PL	INSTRUKCJA UŻYTKOWANIA
CS	NÁVOD K POUŽITÍ
BG	ИНСТРУКЦИИ ЗА УПОТРЕБА
HU	HASZNÁLATI UTASÍTÁS
RU	ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ
LV	LIETOŠANAS PAMĀCĪBA
LT	NAUDOJIMO INSTRUKCIJOS
RO	INSTRUCȚIUNI DE UTILIZARE
SK	NÁVOD NA POUŽITIE
HR	UPUTE ZA UPORABU
ET	KASUTUSJUHEND
SL	NAVODILA ZA UPORABO
FI	KÄYTTÖOHJEET

ONE Snare®

Endovascular Snare System

Standard Snare Kit

INSTRUCTIONS FOR USE

DESCRIPTION

The ONE Snare® Endovascular Snare System consists of the Snare, Snare Catheter, Insertion Tool and Torque Device. The Snare is constructed of nitinol cable and gold plated tungsten loop. The pre-formed Snare loop can be introduced through catheters without risk of snare deformation because of the snare's super-elastic construction. The snare catheter is constructed of polyether block amide and contains a platinum/iridium radiopaque marker band.

USER(S)/ PATIENT POPULATION

Placement of the devices in the ONE Snare Family should be performed by trained healthcare professionals. Clinician specialties typically include interventional radiologists and interventional cardiologists. For use on adult patients requiring retrieval or manipulation of foreign objects in the cardiovascular system, excluding the neurovasculature, where loop size is selected appropriately for the target vessel size.

CLINICAL BENEFITS

The intended clinical benefits of endovascular snares include retrieval or manipulation of foreign objects in the cardiovascular system, excluding the neurovasculature.

INTENDED USE

The ONE Snare® Endovascular Snare System is intended for use in the retrieval and manipulation of foreign objects in the cardiovascular system, excluding the neurovasculature.

INDICATIONS

The ONE Snare® Endovascular Snare System is indicated for use in patients who require retrieval or manipulation of foreign objects in the cardiovascular system.

CONTRAINDICATIONS

1. This device is not intended for the removal of foreign objects entrapped by tissue growth.
2. This device should not be used for fibrin sheath stripping in the presence of septal defect of Patent Foramen Ovale (PFO).
3. This device is not intended for removal of implanted pacing leads.
4. This device is not intended for use in the neurovasculature.

PERFORMANCE CHARACTERISTICS

The ONE Snare® Endovascular Snare System is designed with performance characteristics for use in the cardiovascular system.

- 90 degree pre-formed snare loop for a coaxial approach to foreign object capture
- Nitinol & gold-plated tungsten loop retains snare's circular shape to support capture of foreign objects, while visible under fluoroscopy
- Nitinol core wire provides shape-memory for support during vessel navigation.
- Radiopaque marker band on the snare delivery catheter to identify its precise location.
- These performance characteristics aid in the safe and effective completion of the procedure.

WARNINGS

1. Excessive force used to remove entrapped foreign objects may lead to device failure.
2. Do not use excessive force when manipulating the catheter through an introducer, or when manipulating the snare device. Excessive force may lead to device failure.
3. This device has been sterilized utilizing ethylene oxide and is considered sterile if the package is not opened or damaged. Do not use a device that has been damaged or if the package is open or damaged.
4. For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing, or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
5. After use this device may be a potential biohazard. Handle in a manner that will prevent accidental contamination. Dispose of device in a manner consistent with standard protocols for biohazard waste disposal.
6. Nitinol is a nickel titanium alloy. Possible reaction may occur for those patients who exhibit sensitivity to nickel.
7. There are insufficient safety and performance data to support use of the device in pediatric populations.
8. In the EU, any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the applicable Member State.
9. Use caution when retrieving foreign objects through the heart anatomy to avoid potential tissue / valve damage.

ONLY: CAUTION - Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

POTENTIAL COMPLICATIONS

1. Potential complications associated with foreign object retrieval devices in arterial vasculature include, but are not limited to:
 - Embolization
 - Stroke
 - Myocardial infarction (depending upon placement)
2. Potential complications associated with snare retrieval devices in venous vasculature include, but are not limited to:
 - Pulmonary embolism
3. Other potential complications associated with foreign object retrieval devices include, but are not limited to:
 - Vessel perforation
 - Device entrapment

DIRECTIONS FOR USE

PREPARE THE ONE SNARE® SYSTEM

Select the appropriate Snare diameter range for the site in which the foreign object is located. The Snare diameter range should approximate the size of the vessel in which it will be used.

1. Remove the Snare and Snare Catheter from their hoop holders and inspect for any damage.
2. Remove the Insertion Tool and Torque Device from the proximal end of the Snare shaft.
3. Flush catheter with normal saline.
4. Load the Snare into the Snare Catheter by inserting the proximal end of the snare into the distal (non-hubbed) end of the Snare Catheter, until the proximal end of the Snare shaft exits the hub and the loop can be retracted into the distal end of the Snare Catheter.
5. Test and inspect the device by extending and retracting the Snare loop through the distal end of the Snare Catheter 2-3 times, while carefully examining the Snare Catheter, radiopaque band and the device for any damage or defects.
6. When appropriate, the system (Snare and Snare Catheter) can be advanced to the desired site as a single unit assembled as described above.

ALTERNATIVE PREPARATION OF THE ONE SNARE® SYSTEM

If the Snare Catheter is already positioned within the vasculature, the provided Insertion Tool (located on the proximal end of the Snare and just distal to the Torque Device) may be used to position the snare in the indwelling Snare Catheter.

1. Remove the Snare from the protective holder and inspect for any damage.
2. Move the provided Insertion Tool (located on the proximal end of the Snare, just distal to the torque device) distally until the loop of the Snare is enclosed within the tubing portion of the Insertion Tool.
3. Insert the distal end of the Insertion Tool into the hub of the indwelling Snare Catheter until resistance is felt. This will indicate the tip of the Insertion Tool is properly aligned with the inner lumen.
4. Hold the Insertion Tool as straight as possible, grasp the shaft of the Snare just proximal to the hub of the Insertion Tool and advance the snare until it is secure within the lumen of the Snare Catheter. The Insertion Tool can be removed by grasping the blue tab and firmly peeling it away from the Snare shaft.

SNARE ASSISTED RETRIEVAL AND MANIPULATION SUGGESTIONS

If the Snare Catheter is already positioned within the vasculature, the provided Insertion Tool (located on the proximal end of the Snare and just distal to the Torque Device) may be used to position the snare in the indwelling Snare Catheter.

1. If present, remove the indwelling delivery catheter.
2. If a guidewire is in a patient at the location of a foreign object, advance a Snare Catheter over the guidewire to the desired location. Then remove the guidewire and advance the Snare through the Snare Catheter. Alternatively, cinch one loop of the Snare over the proximal end of the guidewire and advance the entire system (Snare and Snare Catheter assembly) into a guide catheter or introducer sheath until the distal end of the Snare Catheter is positioned proximal to the foreign object.
3. If a guidewire is not present, pull the Snare into the distal end of the Snare Catheter and advance through a guide catheter or introducer sheath until it is positioned proximal to the foreign object. Alternatively, collapse the Snare loop by pulling the device into the distal end of the Snare Insertion Tool. Place the tapered end of the Snare Insertion Tool into the proximal (hub) end of the Snare Catheter, guide catheter or sheath and advance the Snare forward maintaining constant contact between the Insertion Tool and Snare Catheter hub.

NOTE: When attempting to utilize guide catheters or sheaths not specifically manufactured for use with the ONE Snare® system, it is important to test product compatibility prior to use.
4. Gently push the Snare shaft forward to completely open the loop. The loop is then slowly advanced forward, and around the proximal end of the foreign object. Alternatively, the Snare may be advanced beyond the target location and the loops brought back around the distal end of the foreign object.
5. By advancing the Snare Catheter, the loop of the Snare is closed to capture the foreign object.

NOTE: Attempting to close the loop by pulling the Snare into the Snare Catheter will move the loop from its position around the foreign object.
6. To manipulate a foreign object, maintain tension on the Snare Catheter to retain the hold on the foreign object, and move the Snare and Snare Catheter together to manipulate a foreign object to the desired position.

7. To retrieve a foreign object, maintain tension on the Snare Catheter and move the Snare and Snare Catheter assembly together proximally to, or into a guide catheter or sheath. The foreign object is then withdrawn through or together with the guiding catheter or introducer sheath.
NOTE: Withdrawal of large foreign objects may require the insertion of larger sheaths, guiding catheters, or a cut-down at the peripheral site.

ONE SNARE® SYSTEM CONTAINS
(1) Snare, (1) Snare Catheter, (1) Insertion Tool and (1) Torque Device.

	Caution
	Do not use if package is damaged and consult instruction for use
	Catalog number
	Lot number
	Medical Device
	Unique Device Identifier
	Single use
	Do not resterilize
	Consult Instructions for Use For electronic copy scan QR Code, or go to www.merit.com/ifu and enter IFU ID Number. For printed copy available within seven calendar days, call U.S.A. or EU Customer Service.
	Sterilized using ethylene oxide
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Single sterile barrier system
	Use by date: YYYY-MM-DD
	Date of manufacture: YYYY-MM-DD
	Manufacturer
	Authorized Representative in European Community
	Non-pyrogenic
	Inner diameter
	Outer diameter

For a copy of this device’s current European Summary of Safety and Clinical Performance (SSCP), please go to the European database on medical devices (Eudamed), where it is linked to the basic UDI-DI. <https://ec.europa.eu/tools/eudamed>. Alternatively, download a copy of the SSCP from: <https://merit.com/sscp>.

Basic UDI-DI: 088445048822DX

ONE SNARE® ENDOVASCULAR SNARE SYSTEM

Catalogue	Snare			Catheter			
	Diameter Range	Length	Collapsed Snare Diameter (OD)	Size	Length	Catheter Inside Diameter (ID) (MIN)	Catheter Outside Diameter (OD) (MAX)
ONE500	5 mm	120 cm	0.040 in (1.02 mm)	4 Fr	100 cm	0.040 in (1.02mm)	0.057 in (1.45mm)
ONE1000	10 mm	120 cm	0.040 in (1.02 mm)	4 Fr	100 cm	0.040 in (1.02mm)	0.057 in (1.45mm)
ONE1500	15 mm	120 cm	0.050 in (1.27 mm)	6 Fr	100 cm	0.056 in (1.42mm)	0.082 in (2.08mm)
ONE1501	15 mm	80 cm	0.050 in (1.27 mm)	4 Fr	65 cm	0.040 in (1.02mm)	0.057 in (1.45mm)
ONE2000	20 mm	120 cm	0.050 in (1.27 mm)	6 Fr	100 cm	0.056 in (1.42mm)	0.082 in (2.08mm)
ONE2500	25 mm	120 cm	0.050 in (1.27 mm)	6 Fr	100 cm	0.056 in (1.42mm)	0.082 in (2.08mm)
ONE3000	30 mm	120 cm	0.050 in (1.27 mm)	6 Fr	100 cm	0.056 in (1.42mm)	0.082 in (2.08mm)
ONE3500	35 mm	120 cm	0.050 in (1.27 mm)	6 Fr	100 cm	0.056 in (1.42mm)	0.082 in (2.08mm)
ONE1001	10 mm	65 cm	0.040 in (1.02 mm)	4 Fr	48 cm	0.040 in (1.02mm)	0.057 in (1.45mm)
ONE2501	25 mm	65 cm	0.050 in (1.27 mm)	6 Fr	48 cm	0.056 in (1.42mm)	0.082 in (2.08mm)
ONE4000	-	-	-	4 Fr	100 cm	0.040 in (1.02mm)	0.057 in (1.45mm)
ONE6000	-	-	-	6 Fr	100 cm	0.056 in (1.42mm)	0.082 in (2.08mm)

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Merit Medical Systems, Inc.
South Jordan, Utah 84095 U.S.A.
U.S.A. Customer Service 1-800-356-3748
Merit Medical Ireland Ltd