ONE Snare[®]

Endovascular Snare System

Microsnare Kit

EN	INSTRUCTIONS FOR USE
FR	MODE D'EMPLOI
IT	ISTRUZIONI PER L'USO
DE	GEBRAUCHSANWEISUNG
ES	INSTRUCCIONES DE USO
PT	INSTRUÇÕES DE UTILIZAÇÃO
NL	GEBRUIKSAANWIJZING
SV	BRUKSANVISNING
NO	BRUKSANVISNING
DA	BRUGSANVISNING
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
TR	KULLANMA TALIMATLARI
PL	INSTRUKCJA UŻYTKOWANIA
CS	NÁVOD K POUŽITÍ
BG	ИНСТРУКЦИИ ЗА УПОТРЕБА
HU	HASZNÁLATI UTASÍTÁS
RU	ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ
LV	LIETOŠANAS PAMĀCĪBA
LT	NAUDOJIMO INSTRUKCIJOS
RO	INSTRUCȚIUNI DE UTILIZARE
SK	NÁVOD NA POUŽITIE
HR	UPUTE ZA UPORABU
ET	KASUTUSJUHEND
SL	NAVODILA ZA UPORABO
FI	KÄYTTÖOHJEET

ONE Snare®

Endovascular Snare System

Microsnare Kit

INSTRUCTIONS FOR USE

DESCRIPTION

The ONE Snare® Endovascular Microsnare System consists of the Snare, Snare Catheter, Insertion Tool and Torque Device. The Microsnare is constructed of nitinol cable and a gold plated tungsten loop. The pre-formed microsnare loop can be introduced through catheters without risk of snare deformation because of the snare's super-elastic construction. The microsnare catheter is constructed of polyether block amide (Pebax®) and contains a platinum/iridium radiopaque marker band.

USER(S)/ PATIENT POPULATION

Placement of the devices in the ONE Snare Family should be performed by trained healthcare professionals. Clinician specialties typically include interventional radiologists and interventional cardiologists. For use on adult patients requiring retrieval or manipulation of foreign objects in the cardiovascular system, excluding the neurovasculature, where loop size is selected appropriately for the target vessel size.

CLINICAL BENEFITS

The intended clinical benefits of endovascular snares include retrieval or manipulation of foreign objects in the cardiovascular system, excluding the neurovasculature.

INTENDED USE

The ONE Snare® Endovascular Microsnare System is intended for use in the retrieval and manipulation of foreign objects in the cardiovascular system, excluding the neurovasculature.

INDICATIONS

The ONE Snare® Endovascular Microsnare System is indicated for use in patients who require retrieval or manipulation of foreign objects in the cardiovascular system.

CONTRAINDICATIONS

1. This device is not intended for the removal of foreign objects entrapped by tissue growth.
2. This device should not be used for fibrin sheath stripping in the presence of septal defect of Patent Foramen Ovale (PFO).
3. This device is not intended for removal of implanted pacing leads.
4. This device is not intended for use in the neurovasculature.

PERFORMANCE CHARACTERISTICS

The ONE Snare® Endovascular Microsnare System is designed with performance characteristics for use in the cardiovascular system.

- 90 degree pre-formed snare loop for a coaxial approach to foreign object capture
- Nitinol & gold-plated tungsten loop retains snare's circular shape to support capture of foreign objects, while visible under fluoroscopy
- Nitinol core wire provides shape-memory for support during vessel navigation
- Radiopaque marker band on the snare delivery catheter to identify its precise location

These performance characteristics aid in the safe and effective completion of the procedure.

WARNINGS

1. Excessive force used to remove entrapped foreign objects may lead to device failure.
2. Do not use excessive force when manipulating the catheter through an introducer, or when manipulating the snare device. Excessive force may lead to device failure.
3. This device has been sterilized utilizing ethylene oxide and is considered sterile if the package is not opened or damaged. Do not use a device that has been damaged or if the package is open or damaged.
4. For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing, or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
5. After use this device may be a potential biohazard. Handle in a manner that will prevent accidental contamination. Dispose of device in a manner consistent with standard protocols for biohazard waste disposal.
6. Nitinol is a nickel titanium alloy. Possible reaction may occur for those patients who exhibit sensitivity to nickel.
7. There are insufficient safety and performance data to support use of the device in pediatric populations.
8. In the EU, any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the applicable Member State.

9. Use caution when retrieving foreign objects through the heart anatomy to avoid potential tissue / valve damage.

ONLY: CAUTION - Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

POTENTIAL COMPLICATIONS

1. Potential complications associated with snare retrieval devices in arterial vasculature include, but are not limited to:
 - Embolization
 - Stroke
 - Myocardial infarction (depending upon placement)
2. Potential complications associated with snare retrieval devices in venous vasculature include, but are not limited to:
 - Pulmonary embolism
3. Other potential complications associated with foreign object retrieval devices include, but are not limited to:
 - Vessel perforation
 - Device entrapment

DIRECTIONS FOR USE

PREPARE THE ONE SNARE® ENDOVASCULAR MICROSNARE SYSTEM

Select the appropriate Microsnare diameter range for the site in which the foreign object is located. The Microsnare diameter range should approximate the size of the vessel in which it will be used.

1. Remove the Microsnare and Microcatheter from their hoop holders and inspect for any damage.
2. Remove the Insertion Tool and Torque Device from the proximal end of the Microsnare shaft.
3. Flush catheter with normal saline.
4. Load the Microsnare into the Microcatheter by inserting the proximal end of the Microsnare into the distal (non-hubbed) end of the Microsnare Catheter, until the proximal end of the Microsnare shaft exits the hub and the loop can be retracted into the distal end of the Snare Microcatheter.
5. Test and Inspect the device by extending and retracting the Snare loop through the distal end of the Snare Catheter 2-3 times, while carefully examining the Snare Catheter, radiopaque band and the device for any damage or defects.
6. When appropriate, the system (Snare and Snare Catheter) can be advanced to the desired site as a single unit assembled as described above.

ALTERNATIVE PREPARATION OF THE ONE SNARE® ENDOVASCULAR MICROSNARE SYSTEM

If the Snare Microcatheter is already positioned within the vasculature, the provided Insertion Tool (located on the proximal end of the Microsnare and just distal to the Torque Device) may be used to position the Microsnare in the indwelling snare Microsnare Catheter.

1. Remove the Microsnare from the protective holder and inspect for any damage.
2. Move the provided Insertion Tool distally until the loop of the Microsnare is enclosed within the tubing portion of the Insertion Tool.
3. Insert the distal end of the Insertion Tool into the hub of the indwelling Snare Microcatheter until resistance is felt. This will indicate the tip of the Insertion Tool is properly aligned with the inner lumen.
4. Hold the Insertion Tool as straight as possible, grasp the shaft of the Microsnare just proximal to the hub of the Insertion Tool and advance the Microsnare until it is secure within the lumen of the Microcatheter. The Insertion Tool can be removed by grasping the blue tab and firmly peeling it away from the Microsnare Shaft.

POSITION MICROSNARE PROXIMAL TO FOREIGN OBJECT (MICROSNARE ASSEMBLED OUTSIDE THE BODY/ GUIDEWIRE PRESENT)

If the Snare Microcatheter is already positioned within the vasculature, the provided Insertion Tool (located on the proximal end of the Microsnare and just distal to the Torque Device) may be used to position the Microsnare in the indwelling snare Microsnare Catheter.

1. Remove the indwelling balloon or delivery catheter, if present. It may be necessary to change or extend the indwelling guidewire to facilitate balloon removal, and to upsize the indwelling guide catheter/sheath to accommodate the 2.3-3.0F (0.81-1.0 mm) Microsnare Catheter.
2. Place the Microsnare loop over the proximal end of the wire, and cinch until just snug by pulling the proximal end of the Microsnare shaft.
3. Advance the Microsnare, distal end first, into the guide catheter or sheath until the distal end of the Catheter is positioned just proximal to the foreign object.

MICROSNARE ASSEMBLED OUTSIDE THE BODY/GUIDEWIRE NOT PRESENT

1. Pull the Microsnare loop just into the distal end of the Microcatheter.
2. Advance the Microsnare in the guide catheter until it is positioned just proximal to the foreign object.

NOTE: When attempting to utilize guide catheters or sheaths not specifically manufactured for use with the ONE Snare® Microsnare system, it is important to test product compatibility prior to use.

MICROSNARE ASSEMBLED WITH AN INDWELLING CATHETER

1. Advance the Microsnare shaft in the Microcatheter until the loop is just proximal to the Microcatheter end hole.
2. Advance the Microsnare in the guide catheter until it is positioned just proximal to the foreign object.

MICROSNARE ASSISTED MANIPULATION/RETRIEVAL

1. Push the Microsnare shaft forward gently to completely open the loop. Slowly advance the loop forward, and around the proximal end of the foreign object.
2. Once the loop is around the foreign object, advance the Catheter to cinch the loop closed and grab the foreign object.
NOTE: Attempting to close the loop by pulling the Microsnare within the Catheter will move the loop from its position around the foreign object.
3. To manipulate a foreign object, maintain tension on the Catheter to retain the hold on the foreign object, and move the Microsnare and Microcatheter together to manipulate foreign object to the desired position.
4. To retrieve a foreign object, maintain tension on the Microcatheter and move the Microsnare and Catheter together (proximally) to, or into the guide catheter or sheath. The foreign object is then withdrawn through or together with the guiding catheter or vascular sheath.
NOTE: Withdrawal of large foreign objects may require the insertion of larger sheaths, or guiding catheters, or a cut down at the peripheral site.

ONE SNARE® ENDOVASCULAR MICROSNARE SYSTEM CONTAINS

(1) Microsnare, (1) Microsnare Catheter, (1) Insertion Tool and (1) Torque Device.

	Caution
	Do not use if package is damaged and consult instruction for use
	Catalog number
	Lot number
	Medical Device
	Unique Device Identifier
	Single use
	Do not resterilize
	Consult Instructions for Use For electronic copy scan QR Code, or go to www.merit.com/ifu and enter IFU ID Number. For printed copy available within seven calendar days, call U.S.A. or EU Customer Service.
	Sterilized using ethylene oxide
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Single sterile barrier system
	Use by date: YYYY-MM-DD
	Date of manufacture: YYYY-MM-DD
	Manufacturer
	Authorized Representative in European Community
	Non-pyrogenic
	Inner diameter
	Outer diameter

For a copy of this device’s current European Summary of Safety and Clinical Performance (SSCP), please go to the European database on medical devices (Eudamed), where it is linked to the basic UDI-DI. <https://ec.europa.eu/tools/eudamed>. Alternatively, download a copy of the SSCP from: <https://merit.com/sscp>.

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Merit Medical Systems, Inc.
South Jordan, Utah 84095 U.S.A.
U.S.A. Customer Service 1-800-356-3748
Merit Medical Ireland Ltd