INSTRUCTIONS FOR USE

MODE D'EMPLOI

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BRUKSANVISNING

BRUGSANVISNING

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NÁVOD NA POUŽITIE

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MERIT MEDICAL®

Guide Wire Extension**INSTRUCTIONS FOR USE****DESCRIPTION:**

The GO2WIRE™ Guide Wire extension device is a PTFE* coated stainless steel guide wire attachment that is 0.035" (0.089 cm) in diameter and 155 cm length. The extension device has a 152.5 cm long coating of white PTFE shrink tube. The non-coated stainless-steel attachment section is 2.5 cm in length. It is compatible with the shorter length (145cm & 175cm) GO2WIRES that have been modified for attachment to the extension device.
* Polytetrafluoroethylene (PTFE), a synthetic polymer which has lubricious properties.

INDICATIONS FOR USE:

The extension device is intended, when attached to the guide wire, to provide enough additional length to the original guide wire to allow catheter exchanges while maintaining position of the guidewire in the vasculature.

Carefully read all instructions prior to use. Observe all warnings and cautions noted throughout these instructions. Failure to do so may result in complications.

CONTRAINDICATIONS:

No known contraindications. Warnings and intended use of guidewires and interventional devices compatible with the GO2 Wire™ extension are described in the instructions supplied with the respective device.

POTENTIAL COMPLICATIONS:

Potential complications which may result from the use of the device include but are not limited to:

1. Refer to original guide wire instructions for use for potential complications of wire usage.
2. Separation of the extension device from the guide wire body
3. Potential loss of wire tip position in the vasculature
4. Other potential access site complications leading to bleeding, dissection, or perforation that may require intervention.

INSPECTION PRIOR TO USE:

Product is sterile if package is unopened and undamaged. Prior to use, carefully examine all devices to verify that the sterile package or product has not been damaged in shipment. Prior to and during use, inspect the guidewire and extension device carefully. Do not use a damaged device.

PRIOR TO USE:**CAUTION:**

1. The physician should be trained with the use of angiography and angioplasty products and the potential procedural complications.
 2. Confirm the compatibility of the extension wire with other interventional devices being used by testing the systems for any resistance prior to actual use. Free movement of the extension wire within the interventional device must be confirmed and maintained.
 3. The extension wire should be completely hydrated with saline or heparinized saline prior to removal from the flush dispenser hoop.
- NOTE:** Distal tip of extension wire may be positioned inside the dispenser hoop to protect the extension wire sleeve.
4. **CAUTION:** The safety and effectiveness of the Go2 steerable wire has not been established in the Coronary or Neurovasculature
 5. **WARNING:** Preclinical testing with this device showed the potential for clot formation in the absence of anticoagulation. Appropriate anticoagulation therapy should be considered to reduce the potential for thrombus formation on the device.

GUIDE WIRE EXTENSION DEVICE PREPARATION:

1. Employ aseptic technique during removal from the package and use.
 2. Remove the extension wire retaining J straightener from the dispenser hoop to expose the extension wire sleeve.
 3. Attach flush solution filled syringe to flush dispenser hoop luer.
 4. Inject saline until dripping out opposite end to completely fill the dispenser hoop.
 5. Detach syringe from flush dispenser hoop luer
 6. Inspect, prepare, and flush associated device(s) to be used according to the manufacturer's instructions (If applicable).
- NOTE:** To reduce the potential of clot formation, it is recommended that the guide wire extension device be flushed with saline or heparinized saline prior to use.

DIRECTIONS FOR USE:

1. Gently position the proximal end of the guide wire into the threaded extension wire sleeve.
2. Rotate the extension wire sleeve clockwise approximately 3 to 5 rotations until the extension device feels secure to the wire.
- NOTE:** Extension device is fully threaded when the extension portion contacts the white guide wire jacket.
3. Test connection security by grasping the wire and the extension device on each side of the joint and gently pulling to confirm secure connection.
4. Prepare the exchange catheter as directed in the device instructions for use. Withdraw the previously inserted catheter over the extended guide wire.
5. Examine connections to confirm integrity of the joint connection and confirm wire tip position in the vasculature under fluoroscopy before advancing another catheter.
6. Load and advance the catheter onto the proximal end of the guide wire extension device maintaining position of the guide wire tip and catheter under fluoroscopy.
7. When exchange is completed, remove the extension device. Gently grasp the extension wire sleeve and turn counter-clockwise 3 to 5 rotations.
8. Maintaining position of the guide wire in the vasculature, gently pull the extension portion away from the wire body to separate the extension device from the wire.

NOTE: Do not torque or manipulate the guide wire extension device, the extension device does not contain torqueing abilities like the guide wires. Attempted torqueing only results in limited maneuverability of the attached guide wire. The extension device is used for catheter exchanges only. Torqueing and other manipulations may be resumed after removal of the extension device from the guide wire.

DISPOSAL: After use, dispose of product and packaging in accordance with hospital protocol.

RX ONLY: Caution - Federal Law (USA) restricts this device to sale by or on the order of a physician.

REUSE PRECAUTION STATEMENT:

For Single Patient Use Only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

SYMBOL	DESIGNATION
	Use By: YYYY-MM-DD
	Lot Number
	Catalog Number
	Do Not Re-sterilize
	Do Not Use if Package is Damaged
	Single Use
	Caution: Consult accompanying documents. Read instructions prior to use.
	For electronic copy scan QR code or go to www.merit.com/ifu and enter IFU ID Number. For printed copy, call U.S.A or E.U Customer Service.
	Sterilized Using Ethylene Oxide
	Non-pyrogenic
	2D Barcode

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Merit Medical Systems, Inc.
South Jordan, Utah 84095 U.S.A.
U.S.A Customer Service 1-800-356-3748
Merit Medical Ireland Ltd,