MIGHTYTMwire

GUIDE WIRE

EN	INSTRUCTIONS FOR USE
FR	MODE D'EMPLOI
IT	ISTRUZIONI PER L'USO
DE	GEBRAUCHSANWEISUNG
ES	INSTRUCCIONES DE USO
PT	INSTRUÇÕES DE UTILIZAÇÃO
PTBR	INSTRUÇÕES DE USO
NL	GEBRUIKSAANWIJZING
SV	BRUKSANVISNING
NO	BRUKSANVISNING
DA	BRUGSANVISNING
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
TR	KULLANMA TALİMATI
PL	INSTRUKCJA UŻYTKOWANIA
CS	NÁVOD K POUŽITÍ
BG	ИНСТРУКЦИИ ЗА УПОТРЕБА
HU	HASZNÁLATI UTASÍTÁS
RU	ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ
LV	LIETOŠANAS PAMĀCĪBA
LT	NAUDOJIMO INSTRUKCIJOS
RO	INSTRUCȚIUNI DE UTILIZARE
SK	NÁVOD NA POUŽITIE
HR	UPUTE ZA UPORABU
ET	KASUTUSJUHEND
SL	NAVODILA ZA UPORABO

MIGHTY[™]wire

Guide wire

INSTRUCTIONS FOR USE

DEVICE DESCRIPTION

The Mighty Wire[™] Guide Wire will consist of 0.035" guide wire configurations available in 90, 145, 180, 230, 260 and 300cm. Wire lengths 230, 260 and 300cm have 4cm or 7cm of tip flexibility including a radiopaque coil at the distal tip for enhanced visibility. Wire lengths 90, 145, and 180cm have 4cm or 7cm of tip flexibility with no radiopaque coil. The shape of the distal section of the wire will be provided in straight, curved (7.5mm curve radius), double curved (75/15mm curve radius), and extended double curved configurations (55/15mm curve radius).

The 0.035" (0.89mm) wire is composed of a PTFE coated stainless-steel core wire with a flexible tip construction, PTFE coated stainless steel distal coil and radiopaque marker coil as outlined above. The PTFE coated coil and radiopaque coil, where present, are welded to the core on the very distal tip. The PTFE coating covers 100% of the guide wire surface. Across all guide wires within the range, the distal tip is radiopaque.

USER(S)

For use by physicians trained in diagnostic and interventional radiology, cardiology, nephrology, and vascular surgery procedures.

PATIENT POPULATION

The Mighty Wire Guide Wires are designed for use during diagnostic and interventional procedures by trained physicians. Using their education and experience, the physician determines based on the individual patient, the appropriate guide wire to support the associated devices to be used during the procedure. The guide wire navigates the anatomy and facilitates placement of the associated devices.

HOW SUPPLIED

The Mighty Wire Guide Wires are packaged in a plastic hoop, which is fitted with a flush port. This packaging is provided to facilitate compliance with the manufacturer recommended guidelines that the wire be flushed with saline or heparinized saline prior to use (See directions for use - Note). Individually packaged 5 units per box, refer to catalog for ordering information.

INTENDED USE/INDICATIONS FOR USE

The Mighty Wire Guide Wire is intended to be used to facilitate the placement of devices during diagnostic and interventional procedures in the central circulatory system, excluding the coronary arteries and neurovasculature.

CONTRAINDICATIONS

The Mighty Wire Guide Wire is not intended to be used in the coronary arteries and neurovasculature.

CLINICAL BENEFIT STATEMENT

The Mighty Wire Guide Wire has indirect clinical benefits for the patient since it assists other medical devices in achieving their intended purpose, without having a direct therapeutic or diagnostic function itself. It is used to deliver compatible diagnostic or therapeutic medical devices through the vasculature that have a direct therapeutic or diagnostic function.

PERFORMANCE CHARACTERISTICS

The Mighty Wire Guide Wires are designed with performance characteristics for use in a patient's vasculature system.

- Atraumatic distal flexible tips to facilitate introduction into the vasculature.
- Surface coating to allow smooth passage of the guide wire through the vasculature.
- Guide wire body stiffness that supports successful delivery of associated devices.
- These performance characteristics aid in a safe and effective completion of the procedure.

POTENTIAL COMPLICATIONS

Potential complications which may result from the use of the device include but are not limited to:

- Allergic Reaction
- Annular Rupture
- Aortic Regurgitation
- Arrhythmia
- Cardiac Tamponade
- Death
- Hemorrhage
- Hypotension
- Infection
- Thromboembolism
- Tricuspid Valve Injury
- Ventricle Perforation
- Ventricular Fibrillation
- Vessel Perforation

Other potential access site complications leading to bleeding, dissection, or perforation that may require intervention. Some of the stated potential adverse events may require additional surgical intervention.

INSPECTION PRIOR TO USE

Carefully read/observe all warnings and cautions noted in these instructions. Failure to do so may result in complications. Product is sterile if package is not opened or damaged, do not use if package integrity is suspect.

Carefully examine all guide wires to verify that the product has not been damaged in shipment. Inspect the guidewire for coil separation, bends or kinks which may have occurred. Do not use a damaged guide wire.

PRECAUTIONS

1. Physicians should be trained with the use of diagnostic and interventional products and any potential procedural complications.
2. The safety and effectiveness of the Mighty Wire has not been established in the Coronary arteries nor Neurovasculature.
3. Guide wire should be completely hydrated prior to removal from the dispenser hoop.
4. Confirm compatibility of guide wire outer diameter with associated device inner diameter by testing the systems together prior to actual use. Free wire movement within the device must be confirmed and maintained.
5. Physician should confirm the device in which they are intending to utilize in the procedure is the correct device and the device characteristics are as expected.
6. The guide wire is not for selective vessel cannulation or torquing and should not be torqued more than two full revolutions.

WARNINGS

1. Prediclinical testing with this device showed potential for clot formation in the absence of anticoagulation. Anticoagulation therapy should be considered to reduce potential for thrombus formation on the device.
2. Insert flexible end of guide wire into vasculature. Inserting from the rigid end may cause damage to vasculature, tissue, or associated device.
3. Advancement and withdrawal of the wire guide should be performed slowly and carefully. Use extreme caution when withdrawing PTFE coated guide wires back through a metal needle, the sharp edge of the needle may scrape the coating.
4. Withdrawal or manipulation of wire guide through a needle tip or against a sharp edge may result in damage, breakage, leading to embolization.
5. Guide Wire stiffness should be considered during manipulation through the vasculature.
6. If resistance is noted tactically or visually under fluoroscopy (buckling), determine cause, and take appropriate steps to relieve the resistance to avoid the risk of vessel perforation.
7. Always advance, manipulate, and remove the guide wire under fluoroscopic guidance to:
 - a. Ensure the distal tip is intraluminal and in the intended vessel.
 - b. Confirm the guide wire placement and location.
 - c. Prevent potential damage to the vasculature.
8. Hold guide wire in position while manipulating the catheter over the guide wire to prevent unintended distal wire tip movement.
9. During guide wire manipulation through valve leaflets, carefully monitor the leaflet position(s) throughout procedure.
10. In the event of a malfunction of the device and/or changes in the performance of the device, exercise caution as this may indicate a change that may affect the safety of the device.

ONLY: Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

INSTRUCTIONS FOR USE

The following instructions provide technical direction but do not obviate the necessity of formal training in the use of the Mighty Wire Guide Wires. The techniques and procedures described do not represent all medically acceptable protocols, nor are they intended as a substitute for the physician's experience and judgment in treating any specific patient.

PREPARATION AND DIRECTIONS FOR USE

1. Employ aseptic technique during removal from the package and use.
2. Upon removal from pouch, confirm no damage has occurred.
3. To avoid damaging the guide wire tip during removal from the flush hoop, remove the proximal guide wire body from the retention clip.
4. Slide the proximal end of the wire body forward towards the flush port in the flush hoop allowing the distal end of the wire to advance through the flush hoop.
5. Attach flush solution filled syringe to dispenser flush port. Inject solution until dripping out opposite end to completely fill the dispenser loop. Detach syringe.
6. Inspect, prepare, and flush the associated device(s) according to the manufacturer's instructions.
7. Gently grasp guide wire tip and J straightener together as a unit and gently pull forward to withdraw the fragile distal wire tip from the wire dispenser.
8. Carefully dispense guide wire into the flush port of the associated device. Remove J straightener by withdrawing it over the guide wire.
9. Under fluoroscopy and while maintaining position of the vascular access device, advance the guide wire to the targeted site.
10. With the wire guide secured in place, advance the associated device to the target site.
11. Once device is secured in place, remove guide wire, flush catheter following routine protocols.
12. Confirm wire position under fluoroscopy.
13. Inspect guide wire between utilizations to ensure no damage has occurred.

NOTE: To reduce the potential of clot formation, it is recommended that the guide wire be flushed with saline or heparinized saline prior to subsequent uses.

WARNING: RESISTANCE WHEN USING A RIGID GUIDE WIRE

1. Always advance or withdraw a wire slowly. Wire advancement with excessive force may cause coil penetration or other vessel damage.
2. Free movement of the guide wire within a vessel or catheter helps confirm suitable wire pathway.
3. Never push, twist, force or withdraw a guidewire that meets resistance. Resistance may be felt tactically or noted by tip buckling during fluoroscopy. Stop immediately to assess the tip under fluoroscopy to determine cause of resistance and/or the need for additional action to free the guide wire tip.
4. Continued manipulation of a guide wire when resistance is felt may cause guide wire damage, tip separation, and / or vascular injury. Extreme care should be taken when manipulating a catheter and wire combination to prevent possible intravascular tissue damage.
5. The guide wire and associated device should be moved as a unit when possible.
6. When reintroducing a guide wire into a catheter or device already within a vessel, confirm the catheter tip is free within the lumen (i.e., not against the vessel wall).

WARNING: After use, dispose of device in a manner consistent with standard protocols for biohazard waste disposal.

REUSE PRECAUTION STATEMENT

For Single Patient Use Only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

	Caution
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Medical Device
	Unique Device Identifier
	Do Not Use If Package is Damaged and Consult Instruction for Use
	Use by date: YYYY-MM-DD
	Catalog number
	Lot number
	Manufacturer
	Date of Manufacture: YYYY-MM-DD
	Consult Instructions for Use For electronic copy scan QR code, or go to www.merit.com/ifu and enter IFU ID. For printed copy, call U.S.A. or EU Customer Service
	Non-pyrogenic
	Single use
	Do not resterilize
	Keep Dry
	Single sterile barrier system
	Sterilized Using Ethylene Oxide
 e.g. 2D or 1D	2D Barcode

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Merit Medical Systems, Inc.
South Jordan, Utah 84095 U.S.A.
U.S.A. Customer Service 1-800-356-3748