# SAFEGUARD®

## 12 CM PRESSURE ASSISTED DEVICE

DISPOSITIF DE 12 CM ASSISTÉ PAR PRESSION  
DISPOSITIVO A PRESSIONE ASSISTITA SAFEGUARD 12 CM  
12 CM LANGE DRUCKUNTERSTÜTZTE VORRICHTUNG  
DISPOSITIVO ASISTIDO POR PRESIÓN DE 12 CM  
DISPOSITIVO 12 CM ASSISTIDO POR PRESSÃO  
DISPOSITIVO ASSISTIDO POR PRESSÃO DE 12 CM  
12 CM DRUKONDERSTEUNEND HULPMIDDEL  
12 CM TRYCKASSISTERAD ANORDNING  
12 CM TRYKBÆRENDE ANORDNING  
ΣΥΣΚΕΥΗ 12 CM ΜΕ ΥΠΟΒΟΗΘΗΣΗ ΠΙΕΣΗΣ  
12 CM BASINÇ DESTEKLİ CİHAZ  
ГЕМОСТАТИЧЕСКОЕ УСТРОЙСТВО 12 CM  
URZĄDZENIE UCISKOWE 12 CM  
PROSTŘEDEK S TLAKOVOU PODPOROU, 12 CM  
12 CM УСТРОЙСТВО ПОД НАЛЯГАНИЕ  
12 CM-ES ÖNTAPADÓ NYOMÓKÖTÉS

## INSTRUCTIONS FOR USE MODE D'EMPLOI

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## INSTRUCTIONS FOR USE

Read instructions prior to use.

Product not made with natural rubber latex.

### DEVICE DESCRIPTION

The SAFEGUARD 12 cm is a single use disposable device. SAFEGUARD has a clear medical grade polyurethane window and bladder, a clear medical grade PVC flexible fill tube, and a pressure sensitive, self-adhesive peel backing. A valve on the end of the fill tube enables a syringe to be connected to inflate the central bladder with air to provide pressure to the puncture site. The SAFEGUARD pressure assisted device has a sterile dressing with a clear window that facilitates visibility of the access site without removal or manipulation of the device.

### INDICATIONS

The indications for use for the SAFEGUARD 12 cm pressure assisted device are to assist in obtaining and maintaining hemostasis.

### CONTRAINDICATIONS

The adhesive portion of the SAFEGUARD device should not be used over excoriated skin.

### CAUTIONS

**R& Only. Caution:** Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

**CAUTION:** With over-inflation i.e., above 7 mL's of air, the bulb may begin to expand radially and could compromise the adhesive properties of the device.

**CAUTION:** Under-inflation of device could compromise the ability of the device to assist in obtaining and maintaining hemostasis.

### REUSE PRECAUTION STATEMENT

For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

### PRECAUTIONS

- Use proper aseptic techniques while handling product.
- Do not use if package is damaged
- Inspect device prior to use to verify that no damage has occurred during shipping.

### POTENTIAL ADVERSE EFFECTS

Possible adverse effects that may result from the use of this device:

- Hematoma
- Local bleeding
- Arterio-venous fistula or pseudoaneurysm

## PRE-HEMOSTASIS, or MANUAL ASSIST TECHNIQUE

### PLACEMENT OF SAFEGUARD







1. Before adhering SAFEGUARD to the patient, be sure skin is clean and dry. Determine the appropriate angle for SAFEGUARD placement to provide easy access to valve inflate/deflate port and to allow for easy sheath removal.

**Note:** Placement may require adjustment based on the patient's anatomy, angle of the puncture site, and the presence or absence of a procedural sheath.

2. Consider the point of maximum pulse, anatomy, angle of puncture and direction of flow to determine the appropriate SAFEGUARD position and verify.
3. Pull the procedural sheath back approximately 1" (2.5cm) so that when SAFEGUARD is adhered to the skin the sheath hub is outside the area of the SAFEGUARD adhesive.

**Note:** It is recommended that you aspirate the sheath prior to removal to prevent distal embolization from residual clot in sheath.

4. Remove the adhesive backing and place the bulb where you would position your fingers to hold manual compression (for example, in femoral artery procedures, typically the point of maximum femoral pulse). Make sure SAFEGUARD is completely adhered to the skin.
  5. Attach and completely engage the appropriately configured syringe to inflate the desired volume a maximum volume of 7 mL's of air into the bulb to apply pressure on the arteriotomy site. Syringe must be completely engaged in the valve to inflate/deflate the bulb. Remove syringe.
- Note:** Maintain pressure on the plunger while detaching syringe from the SAFEGUARD valve. Observe that the desired pressure is achieved and maintained.
6. Remove sheath, then immediately apply manual compression directly over inflated bulb.
  7. Hold manual compression until hemostasis has been achieved.
    - Slowly release manual compression.
    - Check distal/proximal pulses to assure flow is maintained.
    - Confirm hemostasis by viewing the site through the inflated bulb window.
  8. Per hospital protocol, periodically check the site through the clear window to confirm hemostasis and to manage the bulb volume and resultant pressure as needed. Continue to check distal/proximal blood flow to assure patency.
  9. Deflate bulb every two hours to allow for capillary refill and to assess the site. Re-inflate the bulb if necessary.


10. Deflate the bulb by attaching an appropriately sized and configured syringe to the valve, engage the valve and slowly depress the bulb allowing the syringe to fill with air. Alternatively, remove the plunger from the syringe, attach syringe and allow air to slowly release while gently depressing the bulb.

**Note:** Do not draw negative pressure in the syringe, as this will create a vacuum on the site.

11. Prior to discharge of the patient, remove SAFEGUARD and apply sterile dressing per hospital protocol.


#### POST-HEMOSTASIS TECHNIQUE

1. When hemostasis at the access site has been achieved, apply the SAFEGUARD device with the access site visible under the bulb window of the SAFEGUARD device. Consider the point of maximum pulse, anatomy, angle of puncture and direction of flow to determine the appropriate SAFEGUARD position and verify.


**Note:** Before adhering SAFEGUARD to the patient, be sure skin is clean and dry. Determine the appropriate angle for SAFEGUARD placement to provide easy access to valve inflate/deflate port.

2. Attach an appropriately sized and configured syringe to the valve of the SAFEGUARD device.

**Note:** Syringe must be completely engaged in the valve to inflate/deflate the bulb.

3. Inflate the bulb of the SAFEGUARD device with air to the desired volume of air (12 cm maximum of 7 mL's) to apply pressure on the arteriotomy site and remove the syringe. Check distal/proximal pulses to assure flow is maintained.

**Note:** Maintain pressure on the plunger while detaching syringe from the valve of the SAFEGUARD device. Observe that the desired volume is achieved and maintained.

4. Per hospital protocol, periodically check the site through the bulb window to assure hemostasis is maintained and the bulb maintains pressure.
5. Deflate bulb every two hours and assess the site. Re-inflate the bulb if necessary.
6. Deflate the bulb by attaching an appropriately sized and configured syringe to the valve, engage the valve and slowly depress the bulb allowing the syringe to fill with air. Alternatively, remove the plunger from the syringe, attach syringe and allow air to slowly release while gently depressing the bulb.

**Note:** Do not draw negative pressure in the syringe, as this will create a vacuum on the site.

7. Prior to discharge of the patient, remove the SAFEGUARD device and apply sterile dressing per hospital protocol.

Maintain sterile field during application.

Keep away from sunlight

Do not Resterilize

Do not Reuse

Do not use if package is damaged

Catalog Number

Lot Number

Use By

Sterilized Using Gamma

Caution: Consult accompanying document

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Manufacturer: Merit Medical Systems, Inc.
1600 West Merit Parkway, South Jordan, Utah 84095 U.S.A. 1-801-253-1600
U.S.A. Customer Service 1-800-356-3748
Authorized Representative: Merit Medical Ireland Ltd,