SAFEGUARD®

24 CM PRESSURE ASSISTED DEVICE

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INSTRUCTIONS FOR USE

Read instructions prior to use.

Product not made with natural rubber latex.

DEVICE DESCRIPTION

The SAFE GUARD 24 cm is a single use disposable device. SAFE GUARD has a clear medical grade polyurethane window and bladder, a clear medical grade PVC flexible fill tube, and a pressure sensitive, self-adhesive peel backing. A valve on the end of the fill tube enables a syringe to be connected to inflate the central bladder with air to provide pressure to the puncture site. The SAFE GUARD pressure assisted device has a sterile dressing with a clear window that facilitates visibility of the access site without removal or manipulation of the device.

INDICATIONS

The indications for use for the SAFE GUARD 24cm pressure assisted device are to assist in obtaining and maintaining hemostasis. The device is also indicated in the reduction of active compression time in femoral artery cannulation following diagnostic and interventional procedures with an ACT of 140 seconds or less, using a 6 Fr. and smaller sheath size.

CONTRAINDICATIONS

The adhesive portion of the SAFE GUARD device should not be used over excoriated skin.

CAUTIONS

Rx Only. Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

CAUTION: With over-inflation i.e., above 40 mL's of air, the bulb may begin to expand radially and could compromise the adhesive properties of the device.

CAUTION: Under-inflation of device could compromise the ability of the device to assist in obtaining and maintaining hemostasis.

REUSE PRECAUTION STATEMENT

For single patient use only. Do not reuse, reprocess or sterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

PRECAUTIONS

- Use proper aseptic techniques while handling product.
- Do not use if package is damaged
- Inspect device prior to use to verify that no damage has occurred during shipping.

POTENTIAL ADVERSE EFFECTS

Possible adverse effects that may result from the use of this device:

- Hematoma
- Local bleeding
- Arterio-venous fistula or pseudoaneurysm

	SAFE GUARD MAT N=100	Manual Compression Historical Control N=85
Avg. AC/TTH Manual Compression (minutes)	7.7 ± 3.3	28 Manual Compression Historical Control N=85

	Diagnostic	Interventional	Diagnostic	Interventional
Major Complications				
Access site bleeding requiring transfusion	0	0	0	0
Vascular repair or the need for vascular repair (via surgery, ultrasound-guided intervention, transcatheter embolization, or stent-graft)	0	1	0	2
Any new ipsilateral lower extremity ischemia	0	0	0	0
Surgery for access site related nerve injury	0	0	0	0
Access site-related infection requiring IV antibiotics or extending hospitalization	0	0	0	0
Minor Complications				
Access site hematoma >6 cm (after sheath pull)	1	1	-	-
Bleeding requiring > 30 minutes to re-establish hemostasis	0	0	-	-
Non-treated pseudoaneurysm	0	0	-	-
Non-treated arteriovenous (AV) fistula documented by ultrasound	0	0	-	-
Skin excoration at site of Safeguard after removal of dressing	0	0	-	-
Skin erythema	0	0	-	-
Allergic reaction to adhesive	0	0	-	-
Ipsilateral lower extremity arterial emboli, transient loss of lower extremity pulse, or deep vein thrombosis	0	0	-	-
Access site-related vessel laceration	0	0	-	-
Transient access-site related nerve injury	0	0	-	-
Access site wound dehiscence	0	0	-	-
Localized access site infection treated with intramuscular or oral antibiotics	0	0	-	-

*None Reported

There is no statistical difference in the rate of major complications between the control group and the treated group.

SAFETY AND EFFECTIVENESS RESULTS

A clinical trial was conducted to evaluate the safety and effectiveness of the SAFEGUARD Manual Assist Technique (MAT) in reducing active compression time compared to historical manual compression data. (Reduced Vascular Complications After Percutaneous Coronary Interventions with a Non-mechanical Suture Device: Results from the Randomized RACE Study, Sanborn, TA)

PRE-HEMOSTASIS, or MANUAL ASSIST TECHNIQUE (MAT)

PLACEMENT OF SAFEGUARD

1. Before adhering SAFEGUARD to the patient, be sure skin is clean and dry. Determine the appropriate angle for SAFEGUARD placement to provide easy access to valve inflate/deflate port and to allow for easy sheath removal.

Note: Placement may require adjustment based on the patient's anatomy, angle of the puncture site, and the presence or absence of a procedural sheath.







2. Consider the point of maximum pulse, anatomy, angle of puncture and direction of flow to determine the appropriate SAFEGUARD position and verify.
3. Pull the procedural sheath back approximately 1" (2.5 cm) so that when SAFEGUARD is adhered to the skin the sheath hub is outside the area of the SAFEGUARD adhesive.

Note: It is recommended that you aspirate the sheath prior to removal to prevent distal embolization from residual clot in sheath.

4. Remove the adhesive backing and place the bulb where you would position your fingers to hold manual compression (for example, in femoral artery procedures, typically the point of maximum femoral pulse). Make sure SAFEGUARD is completely adhered to the skin.
5. Attach and completely engage the appropriately configured syringe to inflate the desired volume a maximum volume of 40 mL's of air into the bulb to apply pressure on the arteriotomy site. Syringe must be completely engaged in the valve to inflate/deflate the bulb. Remove syringe.

Note: Maintain pressure on the plunger while detaching syringe from the SAFEGUARD valve. Observe that the desired pressure is achieved and maintained.

6. Remove sheath, then immediately apply manual compression directly over inflated bulb.
7. Hold manual compression until hemostasis has been achieved.*
- Slowly release manual compression.
- Check distal/proximal pulses to assure flow is maintained.
- Confirm hemostasis by viewing the site through the inflated bulb window.

* Recommendations (MAT only):

Diagnostic patients - minimum 5 minutes

Interventional patients - minimum 10 minutes

8. Per hospital protocol, periodically check the site through the bulb window to confirm hemostasis and to manage the bulb volume and resultant pressure as needed. Continue to check distal/proximal blood flow to assure patency.
9. Deflate bulb every two hours to allow for capillary refill and to assess the site. Re-inflate the bulb if necessary.
10. Deflate the bulb by attaching an appropriately sized and configured syringe to the valve, engage the valve and slowly depress the bulb allowing the syringe to fill with air. Alternatively, remove the plunger from the syringe, attach syringe and allow air to slowly release while gently depressing the bulb.

Note: Do not draw negative pressure in the syringe, as this will create a vacuum on the site.

11. Prior to discharge of the patient, remove SAFEGUARD and apply sterile dressing per hospital protocol.

POST-HEMOSTASIS TECHNIQUE

1. When hemostasis at the access site has been achieved, apply the SAFEGUARD device with the access site visible under the bulb window of the SAFEGUARD device. Consider the point of maximum pulse, anatomy, angle of puncture and direction of flow to determine the appropriate SAFEGUARD position and verify.

Note: Before adhering SAFEGUARD to the patient, be sure skin is clean and dry. Determine the appropriate angle for SAFEGUARD placement to provide easy access to valve inflate/deflate port.

2. Attach an appropriately sized and configured syringe to the valve of the SAFEGUARD device.

Note: Syringe must be completely engaged valve to inflate/deflate the bulb.

3. Inflate the bulb of the SAFEGUARD device with air to the desired volume of air (24cm maximum of 40 mL's) to apply pressure on the arteriotomy site and remove the syringe. Check distal/proximal pulses to assure flow is maintained.

Note: Maintain pressure on the plunger while detaching syringe from the valve of the SAFEGUARD device. Observe that the desired volume is achieved and maintained.




4. Per hospital protocol, periodically check the site through the bulb window to assure hemostasis is maintained and the bulb maintains pressure.

5. Deflate bulb every two hours and assess the site. Re-inflate the bulb if necessary.

6. Deflate the bulb by attaching an appropriately sized and configured syringe to the valve, engage the valve and slowly depress the bulb allowing the syringe to fill with air. Alternatively, remove the plunger from the syringe, attach syringe and allow air to slowly release while gently depressing the bulb.

Note: Do not draw negative pressure in the syringe, as this will create a vacuum on the site.

7. Prior to discharge of the patient, remove the SAFEGUARD device and apply sterile dressing per hospital protocol.

	Maintain sterile field during application.
	Keep away from sunlight
	Do not Resterilize
	Do not Reuse
	Do not use if package is damaged
	Catalog Number
	Lot Number
	Use By
	STERILE R Sterilized Using Gamma
	Caution: Consult accompanying document
	Manufacturer
	Authorized Representative

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Manufacturer: Merit Medical Systems, Inc.
1600 West Merit Parkway, South Jordan, Utah 84095 U.S.A.
U.S.A. Customer Service 1-800-356-3748
Authorized Representative: Merit Medical Ireland Ltd,