DISPOSABLE PRESSURE TRANSDUCER - INSTRUCTIONS FOR USE

DESCRIPTION

The Meritrans is a pre-calibrated, single use device for physiological pressure measurement. A separate reusable interface cable is used with this system to connect the transducer to a pressure monitor.

The Meritrans may be used with a manifold, pole mount organizer, or attached directly to the patient via a strap.

SET-UP

CONNECTIONS

1. Using aseptic technique, open the package containing the sterile transducer.
2. Set-up according to hospital protocol for pressure monitoring procedures. Connect to other monitoring equipment (e.g. manifolds, stopcocks, flush devices, tubing administration sets, etc.)

CAUTION: Make sure that all connections are tight. To prevent stripping, do not overtighten.

3. Ensure that the connectors are dry.
4. Connect the Meritrans disposable transducer cable to the reusable monitor cable.
5. Extract all air from the flush solution container.

CAUTION: If air is not extracted from the solution container, air may be forced into the system when the fluid is depleted.

FILL LINES

1. Remove caps (if any) from stopcock on transducer and open the system.
2. Fill all lines with sterile flush solution until free of air bubbles.

NOTE: Merit Medical recommends gravity filling the system rather than pressurizing the flush solution to avoid generating bubbles in the solution.

CAUTION: Verify that no air is trapped in any components of the fluid pathway. Air bubbles can cause serious patient harm and will negatively impact pressure wave forms.

3. Turn stopcock off.
4. Using sterile technique, replace vented caps with non-vented caps.
5. Repeat steps 1-4 for any additional stopcocks or ports.
6. Pressurize the I.V. solution source to 300 mmHg.

ZEROING THE SYSTEM

1. Zeroing should be performed after the system has been filled and mounted.
2. Turn the pressure monitoring system "off" to the patient.
3. Verify that the opening of the stopcock to be used for zeroing is positioned at the patient's mid-axillary level.
4. Being careful not to contaminate the zeroing port, turn the stopcock handle to open the port at the mid-axillary level to air.
5. Zero the transducer according to the monitor manufacturer's instructions.
6. Turn the stopcock handle "off" to the open port.
7. Turn the pressure monitoring system "on" to the patient.

CAUTIONS: Before injecting, turn the manifold handle to the transducer off in order to isolate the transducer from pressure.

CONNECTING TO THE CATHETER

1. Connect the monitoring kit pressure tubing carefully to the patient's catheter or sheath so that no air is introduced into the system.

CAUTIONS:

1. The Meritrans is a single use item. Contents sterile unless package is opened or damaged. Do not reuse or resterilize.
2. Federal (U.S.A.) Law restricts this device to sale by or on the order of a physician.
3. Do not autoclave the reusable cable.
4. If the transducer is used to measure left arterial pressure, an air eliminator filter must be installed between the solution source and the transducer.
5. Components of the pressure monitoring system in contact with the fluid path must be kept sterile.

MERITRANS SPECIFICATIONS

Excitation Voltage	1.0 to 10 Vdc-5kHz	Excitation Impedance	320-350 Ω
Signal Impedance	300 Ω ±15 Ω	Phase Shift	<5°
Sensitivity	5 µV/V/mmHg	Operating Temperature	15°C to 40°C
Storage Temperature	-25°C to 70°C	Maximum Half-Sine	
Zero Drift	1 mmHg/8 hours	Shock Acceleration	4500 G
Thermal Coefficient Offset	±0.3 mmHg/°C	Thermal Coefficient Span	±0.1%/°C
Light Sensitivity	<1 mmHg-		

Manufacturer:
Merit Medical Systems, Inc.
South Jordan, Utah 84095 U.S.A. 1-801-253-1600
U.S.A. Customer Service 1-800-356-3748

Authorized Representative:
Merit Medical Ireland Ltd
Parkmore Business Park West,
Galway, Ireland

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