DESCRIPTION

The Meritrans is a pre-calibrated, single use device for physiological pressure measurement. A separate reusable interface cable is used with this system to connect the transducer to a pressure monitor.

The Meritrans may be used with a manifold, pole mount organizer, or attached directly to the patient via a strap.

SET-UP

Connections

- Using aseptic technique, open the package containing the sterile transducer.
- Set-up according to hospital protocol for pressure monitoring procedures. Connect to other monitoring equipment (e.g. manifolds, stopcocks, flush devices, tubing administration sets, etc.)

CAUTION: Make sure that all connections are tight. To prevent stripping, do not overtighten.

- Ensure that the connectors are dry.
- Connect the Meritrans disposable transducer cable to the reusable monitor cable.
- Extract all air from the flush solution container.

CAUTION: If air is not extracted from the solution container, air may be forced into the system when the fluid is depleted.

FILL LINES

- Remove caps (if any) from stopcock on transducer and open the system.

- Fill all lines with sterile flush solution until free of air bubbles.

NOTE: Merit Medical recommends gravity filling the system rather than pressurizing the flush solution to avoid generating bubbles in the solution.

CAUTION: Verify that no air is trapped in any components of the fluid pathway. Air bubbles can cause serious patient harm and will negatively impact pressure wave forms.

- Turn stopcock off.
- Using sterile technique, replace vented caps with non-vented caps.
- Repeat steps 1-4 for any additional stopcocks or ports.
- Pressurize the I.V. solution source to 300 mmHg.

ZEROING THE SYSTEM

- Zeroing should be performed after the system has been filled and mounted.
- Turn the pressure monitoring system "off" to the patient.
- Verify that the opening of the stopcock to be used for zeroing is positioned at the patient's mid-axillary level.
- Being careful not to contaminate the zeroing port, turn the stopcock handle to open the port at the mid-axillary level to air.
- Zero the transducer according to the monitor manufacturer's instructions.
- Turn the stopcock handle "off" to the open port.
- Turn the pressure monitoring system "on" to the patient.

cautions: Before injecting, turn the manifold handle to the transducer off in order to isolate the transducer from pressure.

CONNECTING TO THE CATHETER

- Connect the monitoring kit pressure tubing carefully to the patient's catheter or sheath so that no air is introduced into the system.

CAUTIONS:

- The Meritrans is a single use item. Contents sterile unless package is opened or damaged. Do not reuse or resterilize.
- Federal (U.S.A.) Law restricts this device to sale by or on the order of a physician.
- Do not autoclave the reusable cable.
- If the transducer is used to measure left arterial pressure, an air eliminator filter must be installed between the solution source and the transducer.
- Components of the pressure monitoring system in contact with the fluid path must be kept sterile.

MERITRANS SPECIFICATIONS

Excitation Voltage	1.0 to 10 Vdc-5kHz	Excitation Impedance	240-350 Ω
Signal Impedance	300 Ω ±30 Ω	Phase Shift	<5°
Sensitivity	5 μV/V/mmHg	Operating Temperature	15°C to 40°C
Storage Temperature	-25°C to 70°C	Maximum Half-Sine Shock Acceleration	4500 G
Zero Drift	1 mmHg/8 hours	Thermal Coefficient Span	±0.1%/°C
Thermal Coefficient Offset	±0.3 mmHg/°C		
Light Sensitivity	<1 mmHg-		

REUSE PRECAUTION STATEMENT

For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

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