DeVos MANIFOLD™

MANIFOLD WITH INTEGRAL CHECK RELIEF VALVE

INSTRUCTIONS FOR USE

DESCRIPTION:

Manifold with integrated check relief valve. The check relief valve houses a flexible membrane which prevents the backflow of solutions and directs the flow of fluids in one direction only.

INDICATIONS AND USAGE:


The manifold with integral check relief valve is used in cardiology to control the flow of contrast without having to manipulate handles.

CAUTION: Federal (USA) law restricts this device to use by or on the order of a physician. This device is intended for single patient use. Read instructions prior to use.

INSTRUCTIONS FOR USE:

Inspect the device prior to use to verify that no damage has occurred during shipping and that the sterile barrier has not been invaded.

1. Remove device from package.
2. Ensure that all lines are firmly connected to the side ports of the manifold.
3. Firmly connect the contrast supply line to the CRV port.
4. Attach the male luer lock syringe to the proximal (female luer lock) end of the manifold.
5. Check connectors, manifold body and fluid paths for air bubbles.
6. If a pressurized contrast system is used (300 mmHg) the check relief valve will allow automatic refilling of the control syringe simply by aspirating to desired ml volume.
7. The system is now ready for use in angiographic procedures.

PRECAUTIONS:

Used correctly, this device will only allow aspiration and injection of contrast. It cannot be used to check patency of placement of the catheter by trying to aspirate blood.

Caution: This product can expose you to chemicals including Di(2-ethylhexyl)phthalate, (DEHP), which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

CE 0086

Manufacturer:
Merit Medical Systems, Inc. South Jordan, Utah 84095 www.merit.com
U.S.A. 1-801-253-1600 U.S.A. Customer Service 1-800-356-3748

EC REP Authorized Representative:
Merit Medical Ireland Ltd, Parkmore Business Park West, Galway, Ireland

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