DeVos MANIFOLD™

MANIFOLD WITH INTEGRAL CHECK RELIEF VALVE

INSTRUCTIONS FOR USE

DESCRIPTION:

Manifold with integrated check relief valve. The check relief valve houses a flexible membrane which prevents the backflow of solutions and directs the flow of fluids in one direction only.

INDICATIONS AND USAGE:


The manifold with integral check relief valve is used in cardiology to control the flow of contrast without having to manipulate handles.

CAUTION: Federal (USA) law restricts this device to use by or on the order of a physician. This device is intended for single patient use. Read instructions prior to use.

INSTRUCTIONS FOR USE:

Inspect the device prior to use to verify that no damage has occurred during shipping and that the sterile barrier has not been invaded.

1. Remove device from package.
2. Ensure that all lines are firmly connected to the side ports of the manifold.
3. Firmly connect the contrast supply line to the CRV port.
4. Attach the male luer lock syringe to the proximal (female luer lock) end of the manifold.
5. Check connectors, manifold body and fluid paths for air bubbles.
6. If a pressurized contrast system is used (300mmHg) the check relief valve will allow automatic refilling of the control syringe simply by aspirating to desired ml volume.
7. The system is now ready for use in angiographic procedures.

PRECAUTIONS:

Used correctly, this device will only allow aspiration and injection of contrast. It cannot be used to check patency of placement of the catheter by trying to aspirate blood.

This product is ethylene oxide sterilized product for single use only. Under the condition that the package is not opened or damaged, this product is an aseptic and non-pyrogenic product.

Shelf life: 3 years.

Storage Conditions: Room temperature, ventilation, dry, away from light.

Transportation Conditions: Do not expose to excessive heat or humidity during transportation.

www.merit.com

Manufacturing:

Merit Medical Systems, Inc. 1600 West Merit Parkway,
So. Jordan, Utah 84095 U.S.A. 1-801-253-1600
U.S.A. Customer Service 1-800-356-3748

Authorized Representative:

Merit Medical Ireland Ltd, Parkmore Business Park West, Galway, Ireland

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South Jordan, Utah 84095 U.S.A. 1-801-253-1600