# ACCESSPLUS™ / ACCESS-9 LARGE BORE HEMOSTASIS VALVE

## I N S T R U C T I O N S F O R U S E :

### DESCRIPTION:

Adjustable Merit AccessPLUS™/ Merit Access-9™ Hemostasis Valve.

### INDICATIONS AND USAGE:

The AccessPLUS/Access-9 Hemostasis Valve is recommended for maintaining a fluid-tight seal around percutaneous transluminal angioplasty catheters and guidewires.

**CAUTION:** Federal (U.S.A.) law restricts this device to use by or on the order of a physician. This device is intended for single patient use only. Read instructions prior to use. Read manufacturer's instructions for the use of catheters, guidewires and introducers.

### INSTRUCTIONS FOR USE:

Inspect the device prior to use to verify that no damage has occurred during shipping.

1. Connect the sideport of the hemostasis valve to the Merit manifold assembly. Flush and fill the assembly with saline. To fill the valve section, open the valve, place one finger over the luer fitting and continue to fill the assembly.
2. Connect the hemostasis valve to the guiding catheter. Aspirate the valve to remove any trapped air and flush thoroughly with saline. Purge blood by opening the valve while continuing to flush the assembly. Close the valve when the blood has been purged. Inspect carefully for air bubbles and reflux if necessary.

Manufacturer: Merit Medical Systems, Inc., So. Jordan, UT 84095 U.S.A. 1-801-253-1600 • U.S.A. Customer Service 1-800-356-3748  
Authorized Representative: Merit Medical Ireland Ltd, Parkmore Business Park West, Galway, Ireland

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3. Open the hemostasis valve and insert the guidewire or catheter. (If desired, insert the guidewire introducer and advance it through the valve.) Advance the guidewire or catheter an appropriate distance into the vasculature. Close the valve around the shaft of the guidewire, catheter or introducer. This forms a fluid-tight seal, yet does not inhibit movement of the guidewire, catheter or introducer.

**WARNING:** It is important that the valve be closed tight enough to prevent blood leakage, yet not so tight as to restrict function.

4. Refer to the dilatation catheter labeling for intended use, contraindications and potential complications associated with the use of dilatation catheters in PTCA and/or PTA procedure.
5. Withdraw the catheter to the point that the distal tip remains 10 to 20 cm inside the vasculature. Open the valve and withdraw the catheter completely. During insertion and removal through the valve, the balloon must be completely deflated to safely and easily traverse the valve.
6. Disconnect the hemostasis valve from the guiding catheter.

**PRECAUTIONS:** Do not inject any fluid if air bubbles are visible within the valve. First aspirate the valve to remove the air, then flush the valve as described above. The hemostasis valve must be completely closed during aspiration or injection. The narrowest portion of the hemostasis valve has an inner diameter of 0.120 inches (0.305 cm) (approximately 9 Fr).

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# ACCESSPLUS™ / ACCESS-9 VALVOLA EMOSTATICA LARGE LUMEN

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