
INTELLISYSTEM[®]25

Inflation Device and Fluid Dispensing Syringe

INSTRUCTIONS FOR USE

U.S. Patent Nos. 5,135,488 5,711,302.
Other U.S. and Foreign Patents Pending.

DESCRIPTION:

The Merit IntelliSystem25 Inflation Device and Fluid Dispensing Syringe is a 20ml disposable syringe with an integral pressure transducer, threaded plunger assembly with lock/release bar, a flexible high pressure extension tube, and a three-way medium pressure stopcock. The IntelliSystem25 is designed to generate and monitor pressures and dispense fluid with each turn of the handle.

INSTRUCTIONS FOR USE:

This device is intended to be used only with the IntelliSystem Monitor. Read IntelliSystem Monitor Instructions For Use before use.

Before use, inspect the device and packaging to verify that no damage has occurred as a result of shipping.

NOTE: To generate and monitor pressures over 20 ATM, 20 BAR or 300 PSI the IntelliSystem Monitor must have 6.0 version software or greater.

CAUTION: Federal (U.S.A.) law restricts this device to use by or on the order of a physician.

NOTE: The IntelliSystem25 syringe is packaged with a perforated bag over the cable connectors. Leave the bag in place until the syringe has been prepared with fluid. This will prevent fluid contamination on the cable connectors.

WARNING: This device should not be used in the presence of flammable anesthetics.

Manufacturer: Merit Medical Systems, Inc.
South Jordan, Utah 84095 U.S.A. 1-801-253-1600
U.S.A. Customer Service 1-800-356-3748

Authorized Representative:
Merit Medical Ireland Ltd
Parkmore Business Park West, Galway, Ireland

400646002_003 ID 2020-09-24

PRIMING WITH PRIMELOK™

The PrimeLok allows purging of air and excess fluid without squeezing the lock/release mechanism. The IntelliSystem25 is packaged with the PrimeLok and with the lock/release mechanism in the disengaged position.

1. To engage PrimeLok, squeeze handle and slide PrimeLok into plunger.




2. To prep syringe, simply aspirate up to 20ml of contrast solution or fluid to be dispensed into the syringe by pulling back on the plunger handle.

CAUTION: Inspect the syringe tubing and stopcock (if used) to insure that there is no air in the system.

3. Push handle against table to remove air in syringe.

4. To disengage squeeze handle and slide PrimeLok out of barrel. This will allow the plunger to lock in position and is now ready for use.

NOTE: The PrimeLok must be disengaged before pressure can be maintained by the lock/release mechanism.

ATTACHING THE INFLATION DEVICE TO THE BALLOON:

This inflation device is used to inflate and deflate an angioplasty balloon or other interventional device, in conjunction with the Merit IntelliSystem Monitor and to measure the pressure within the balloon. Refer to the IntelliSystem Monitor's User Manual for additional information on its function.

NOTE: Refer to the manufacturer's directions accompanying the balloon dilatation catheter or other interventional device for specific information on use, maximum inflation pressure, precautions, and warnings for that device.

1. Prepare and test the balloon catheter according to the catheter manufacturer's directions for use.

2. Create a fluid-fluid connection between the catheter and the inflation syringe extension tube and connect the luer connectors securely.

3. Squeeze the lock/release bar and pull back on the plunger handle to apply a vacuum to the balloon when indicated in the instructions for use of the interventional device.

BALLOON INFLATION AND DEFLATION:

1. To inflate the balloon, squeeze the lock/release bar allowing the plunger to return to zero pressure as displayed on the IntelliSystem Monitor. Release grip on the lock/release bar, locking the plunger into position. To increase pressure, rotate the plunger clockwise until the desired inflation pressure is reached. The lock mechanism maintains the pressure.

NOTE: Loss of pressure may indicate a leak in the system. Check connections.

2. To deflate the balloon, squeeze the lock/release bar and pull back to generate a negative pressure. Release grip on the lock/release bar to lock the plunger in the negative pressure position.

FLUID DISPENSING INDICATIONS AND USAGE:

This device is intended for use by healthcare professionals to dispense fluids to the body from the IntelliSystem25 Syringe and monitor the pressure of that fluid when used with the IntelliSystem Monitor.

1. To slowly dispense fluids to the body rotate the plunger clockwise until the desired fluid is injected. The injection pressure and time will be displayed on the Monitor .

2. To rapidly dispense fluids to the body squeeze the lock/release bar while pushing the plunger forward.

3. A negative pressure may be generated by squeezing the lock/release bar and pulling back the plunger. Release grip on the lock/release bar to lock the plunger in the negative pressure position.

CAUTION: This syringe is capable of generating high fluid pressures in a "closed" fluid system. The volume change of fluid dispensed may not be accurate due to compliance of the plastic components as pressure changes.

CAUTION: If applied pressure does not indicate on gauge/digital display, discontinue use immediately and replace it with a new unit.

SPECIFICATIONS :

Fluid Dispensed: 0.48 ±.07ml fluid dispensed for each 360° clockwise turn of the syringe plunger handle

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