PeriVac™ Kit

PERICARDIOCENTESIS KIT

INSTRUCTIONS FOR USE

English

INTENDED USE OF PRODUCT

INDICATIONS:
The Pericardiocentesis Kit is intended to withdraw fluid from the pericardial sac (membrane surrounding the heart) to evaluate the fluid for diagnostic purposes or to remove excess fluid causing compression of the heart.

CONTRAINDICATIONS:
There are no known contraindications for pericardial aspiration, although recurrent effusion or unresolved cardiac tamponade may warrant surgical intervention.


COMPLICATIONS/RISKS:
Complications and/or risks of pericardiocentesis procedures can include catheter occlusion, infection, cardiac tamponade from puncture, atrial and ventricular arrhythmias, laceration of coronary vessel or chamber wall, air injection into cardiac chambers, pneumothorax, hemotorax, and hypotension (perhaps reflexogenic). Perforation of the stomach, colon, right atrium and/or ventricle is also possible.

WARNINGS:
RX ONLY Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

- This device is intended for single use only.
- Store in a cool dry place.
- A guide wire should never be advanced or withdrawn against resistance. The cause of the resistance should be determined before continuing.

PRECAUTION:
Pericardiocentesis should only be performed by clinicians well-versed in the technique and aware of its risks and benefits.

DESCRIPTION OF KIT
The PeriVac Pericardiocentesis Kit is a sterile, one time use kit designed to assist the clinician in draining fluid from the pericardial sac. Pericardiocentesis procedures are ideally carried out in (but not limited to) the special procedures lab or cardiac cath lab. Pericardiocentesis is a sterile procedure which should be performed with fluoroscopic, ultrasonic, two-dimensional echocardiographic, or electrocardiographic monitoring and guidance.

INSTRUCTIONS FOR USE

1. Patient should be hemodynamically, electrocardiographically, and/or echocardiographically monitored throughout the procedure.
2. Pericardiocentesis should be performed as a sterile procedure. Using sterile technique, prep the selected puncture site.
3. Using sterile technique, open the PeriVac Pericardiocentesis Kit. Open the outer sterile wrap to prepare the sterile field.
4. Remove the backing from the plastic fenestrated drape and place it over the prepped puncture site.
5. Aspirate the preferred anesthetic, via the 18, 19 gauge needle into the 10ml syringe. Anesthetize the area with the 22, 25 gauge needle.
6. Make a small skin incision, approximately 3mm in length, with the #11 scalpel.
7. Attach one end of the alligator clip to the 18 gauge needle, and the other to the Y-lead on the EKG monitoring system.
8. Insert the 18 gauge needle through the incision. Local anesthesia can be instilled through the needle as it is advanced toward the pericardium. Aspiration during passage of the needle may block the needle with subcutaneous tissue. For this reason, tissue is removed by flushing the needle before passing through the pericardium.
9. After access into the pericardium is confirmed, pass the J guide wire through the needle and into the pericardial space. Do not force the guide wire; it should pass without resistance.
10. While maintaining wire position, remove the needle over the guide wire. Pass the dilator over the guide wire, dilating the subcutaneous tissue, diaphragm, and pericardium.
11. Remove the dilator and pass the drainage catheter over the guide wire, positioning the catheter into the pericardial space. Note that the guide wire should always extend beyond the catheter tip and a portion should always extend out the back of the catheter hub during the positioning of the catheter.
12. After the catheter is in position remove the guide wire leaving the catheter in place. Attach the 3-way stopcock to the female luer hub of the catheter.
13. Laboratory samples of pericardial fluid may be drawn at this time. After samples have been drawn, the drainage bag can be attached to the 3-way stopcock for fluid collection.
14. If the catheter is to be left in place for a period of time, secure the catheter to the skin and dress the site per hospital protocol.
15. If the catheter is to be removed, gently withdraw it and dress the site accordingly. Note: If the pigtail catheter is used, removal over the guide wire is recommended.
16. Catheter occlusion caused by coagulation of pericardial effusions can be minimized by interrupted drainage of the catheter every 4-6 hours.

PeriVac is a trademark of Boston Scientific Corporation or its affiliates.

SYMBOL	DESIGNATION
Rx ONLY	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician
LOT	Lot Number
(2)	Single use. Do not re-use.
!	Caution: Consult accompanying documents
(2)	Do not sterilize.
X	Non-pyrogenic
(2)	Do not use if package is damaged.
■	Manufacturer
STERILE EO	Sterile Using Ethylene Oxide
✉	Contents
DATE	Use By Date
REF	Catalog

This product can expose you to chemicals including Di(2-ethylhexyl)phthalate, (DEHP), which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

Manufactured For:
Boston Scientific Corporation
300 Boston Scientific Way
Marlborough, MA 01752 USA
U.S.A. Customer Service 1-800-356-3748

Manufactured For:
Boston Scientific Corporation
300 Boston Scientific Way
Marlborough, MA 01752 USA
U.S.A. Customer Service 1-800-272-1001

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