**CAUTION:** Rx Only: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

**DEVICE DESCRIPTION:** The SlipNot is a suture retention device with an integral threading mechanism.

**INTENDED USE:** The suture retention device is intended to temporarily secure sutures and then aid clinicians in locating and removing sutures efficiently.

#### CONTRAINDICATIONS

- There are no known contraindications for this device.

#### PRECAUTIONS

- The device is intended for single use only. Do not resterilize and/or reuse this device.
- Do not use the product if there is doubt as to whether the product is sterile.
- Store in a dark, dry, cool place. Avoid extended exposure to light.
- Upon removal from package, inspect the product to ensure no damage has occurred.
- Repositioning of the device is possible after deployment; use caution to insure device is in appropriate position to obtain best possible clinical results.

#### POTENTIAL ADVERSE EVENTS

- Undesired bleeding at the site could occur, if caution is not used when securing the device at deployment.
- Hematoma/Bleeding
- Patient Fever

#### PRODUCT RECOMMENDATIONS

For use with suture up to USP #0.

#### INSTRUCTIONS FOR USE

1. Purse-string suture is placed around the existing vascular sheath, puncture site, or skin area to be secured.
  2. The free ends of the purse-string suture are passed through the snare loop of the device.
  3. Depress the white button; pull the tab to thread the free ends of the purse-string through the device. Discard the snare with pull tab.
  4. Depress the white button and slide the device forward to produce tension around the site, enabling removal of the sheath if present. Release the white button to secure the suture within the device.
  5. To inspect the exit site to determine if bleeding is controlled, depress the white button and gently pull the device away from the exit site.
  6. Use caution not to completely pull the device off the purse-string suture.
  7. If bleeding is not controlled, apply tension by depressing the white button and sliding the device forward to produce tension at the site.
  8. To remove the suture, depress the white button and pull the device slightly away from the exit site. Cut one side of the purse-string suture near the skin. Grasp the device and pull back to remove suture.
- Alternate method: Remove the device completely from the purse-string suture. Cut one side of the purse-string suture and remove.

#### HOW SUPPLIED




Supplied sterile in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged.

#### REFERENCES

These instructions for use are based on physician experience and/or their published literature. Refer to your local Merit sales representative for information on available literature.

#### PATENTS

This product and its use are protected by one or more of the following patents: United States Patent No. 7,087,060. Other U.S. patents pending. Foreign patents issued and pending.

Manufacturer: Merit Medical Systems, Inc. South Jordan, Utah 84095  
U.S.A. 1-801-253-1600 U.S.A. Customer Service 1-800-356-3748

EC REP Authorized Representative:  
Merit Medical Ireland Ltd, Parkmore Business Park West, Galway, Ireland

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