MERIT[®] MANIFOLDS

EN	INSTRUCTIONS FOR USE
FR	MODE D'EMPLOI
IT	ISTRUZIONI PER L'USO
DE	GEBRAUCHSANWEISUNG
ES	INSTRUCCIONES DE USO
PT	INSTRUÇÕES PARA UTILIZAÇÃO
PTBR	INSTRUÇÕES DE USO
NL	GEBRUIKSINSTRUCTIES
SV	ANVÄNDNINGSSINSTRUKTIONER
NO	BRUKSANVISNING
DA	BRUGSANVISNING
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
TR	KULLANMA TALİMATLARI
PL	INSTRUKCJA UŻYTKOWANIA
CS	NÁVOD K POUŽITÍ
BG	ИНСТРУКЦИИ ЗА УПОТРЕБА
HU	HASZNÁLATI UTASÍTÁS
RU	ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ
LV	LIETOŠANAS INSTRUKCIJA
LT	NAUDOJIMO INSTRUKCIJOS
RO	INSTRUCȚIUNI DE UTILIZARE
SK	NÁVOD NA POUŽITIE
HR	UPUTE ZA UPORABU
ET	KASUTUSJUHEND
SL	NAVODILA ZA UPORABO

MERIT[®] MANIFOLDS

English


INSTRUCTIONS FOR USE

Product Description

The Merit manifold consists of a clear plastic body with colored handles that control the flow of fluids through multiple ports. The handles of the manifold are prominently marked with arrows indicating the direction of flow. The handles are easy to rotate but have sufficient torque to prevent inadvertent movement. Merit manifolds are available with various types of fitting connectors.

Catalog Numbers

Manifold catalog numbers are based on the following logic:

Merit's DeVos manifolds are a standard manifold with a check relief valve as the end port. The DeVos manifolds are available with:

- 2 or 3 port standard port spacing and 4 port wide
- Half body style
- Right port orientation
- 180°/360°
- 200 PSI (14 BAR) pressure rating
- OFF handles (ON handles)

Intended Use

Merit Manifolds are indicated for use in diagnostic and interventional applications to control or direct fluid flow between tubing, catheters, or other devices.

Contraindications

None known

Warning: Merit Manifolds are not recommended for use with lipids. Prolonged exposure to lipid solutions may result in stress cracking or leakage.

Cautions

- **R_X ONLY:** Federal (U.S.A.) law restricts this device to use by or on the order of a physician.
- After use, dispose of device in a manner consistent with standard protocols for biohazard waste disposal.

Precautions

- Carefully read instructions before using product. If product is being used in conjunction with other manufacturers' components, also read Instructions for Use.
- Use proper aseptic techniques while handling product.
- Inspect device prior to use to verify that no damage has occurred during shipping.

- Only use standard Luer connection devices. A standard Luer connection must conform to the harmonized standard ISO 594-2.
- **DO NOT OVER-TIGHTEN** connections.
- **DO NOT USE** any instrument to tighten connections.

REUSE PRECAUTION STATEMENT

For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or reesterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or reesterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

In the EU, any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the applicable Member State.

Instructions for Use

1. Using aseptic technique, open the package containing the sterile product.
2. Inspect for damage or improper assembly.
3. Check all connections before use and finger tighten.

Warning: To prevent stripping, do not over tighten.

4. Prime the manifold before use.

Note: Ensure all connections are securely tightened.

5. Inspect carefully for air bubbles and if necessary flush the lumen. Ensure that all air bubbles are removed.
6. Attach fluid devices/tubing. Ensure that all connections are secure.
7. Rotate the handles to the appropriate position to get the desired flow path.
 - a. The molded arrows on the manifold handles indicate the open port flow paths.
 - b. The molded "off" on the handle indicates a closed port preventing flow of fluid.

SYMBOL	DESIGNATION
	Catalog Number
	Lot Number
	Caution
	Consult Instructions for Use. For electronic copy scan QR Code, or go to www.merit.com/ifu and enter IFU ID Number. For printed copy, call U.S.A or E.U. Customer Service
	Sterilized using Ethylene Oxide
	Do not use if package is damaged.
	Single Use
	Use By
	Do Not Resterilize
	Non-pyrogenic
	Single sterile barrier system
	Manufacturer
	Authorized Representative in European Community
	Medical Device
	Date of Manufacture
	Unique Device Identifier

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Merit Medical Systems, Inc.
1600 West Merit Parkway, South Jordan, Utah 84095 U.S.A.
U.S.A Customer Service 1-800-356-3748
Merit Medical Ireland Ltd, Parkmore Business Park West, Galway, Ireland