

DiamondTOUCH™

English

INFLATION DEVICE AND FLUID DISPENSING SYRINGE

DiamondTOUCH™

DISPOSITIF DE GONFLAGE ET SERINGUE DE DOSAGE

INSTRUCTIONS FOR USE

DESCRIPTION: The DiamondTOUCH™ Inflation Device and Fluid Dispensing Syringe is a 30ml dispensing device with an inflation function, pressure transducer, microprocessor, barometric sensor, LCD touchscreen display, grayscale graphics, flexible tube or balloon, pressure extension bar, and a three way stopcock. The DiamondTOUCH™ is designed to generate positive and negative pressure, and monitor positive pressures over a range of zero to +35 ATM/BAR (zero to 514 PSI).

Rx Only: Caution: (U.S.A.) Law restricts this device to use by or on the order of a physician.

INDICATIONS FOR USE: The DiamondTOUCH™ Syringe is used to inflate and deflate balloon angioplasty catheters or other interventional devices and to measure the pressure and time of inflation within the balloon during the procedure. It is also used to dispense fluid and measure the pressure and time of inflation within the balloon during the procedure.

DEVICE INFORMATION: The DiamondTOUCH™ Syringe is used to inflate and deflate the balloon of the balloon catheter during the procedure. Press the blue button behind the device to turn the inflation function on. Turn the dial to increase or decrease the pressure until the desired pressure is reached. Then press the blue button again to turn off the inflation function.

PREPARATION OF THE DEVICE: The DiamondTOUCH™ Syringe is used to inflate and deflate the balloon of the balloon catheter during the procedure. Turn the dial to increase or decrease the pressure until the desired pressure is reached. Then press the blue button again to turn off the inflation function.

OPERATION OF THE DEVICE: The DiamondTOUCH™ Syringe is used to inflate and deflate the balloon of the balloon catheter during the procedure. Turn the dial to increase or decrease the pressure until the desired pressure is reached. Then press the blue button again to turn off the inflation function.

NOTICE: When in PSI mode, the tick marks on the left of the display that represent pressure will be limited to 350 PSI.

In the event the DiamondTOUCH™ is pressed in ATM/BAR mode, the grouping of tick marks on the left will flash. The maximum digital number at the center of the display will continue to show actual pressure throughout the device's pressure range (+4 to 514 PSI).

After a single or pressure measured injection has been made, a graph bar will remain to mark the highest point of pressure. Press the blue button once quickly to turn off the inflation function and, if you are representing the pressure, the pressure will be displayed in ATM/BAR mode. The DiamondTOUCH™ Syringe will accept a pressure greater than 350 PSI, the reprogramming of dispense fluid duration is up to the user.

After the inflation function has been started, the last inflation tick mark will disappear.

CAUTION: If the LCD displays anything besides the pressure and time windows shown above, the syringe is defective. Please return the syringe to Merit Medical for repair.

NOTE: To conserve power the backlight will automatically turn off after ten minutes of inactivity. However, the microprocessor will continue to monitor the pressure. Pressing the blue button or inflating the balloon will reactivate the backlight. The device will power itself off after 90 consecutive minutes at zero pressure.

ATTENTION: Inspect the syringe tubing and stopcock if used to ensure that there is no air in the system.

ATTACHING THE INFILATION DEVICE TO THE BALLOON: NOTE: To the manufacturer's directions concerning the balloon inflation catheter or other interventional devices for specific information on use, maximum inflation pressure, precautions, and warnings for that device.

1. Prepare and test the balloon catheter according to the catheter manufacturer's directions for use.

2. Create a fluid connection between the balloon and the syringe extension tube, connect the luer connector to secure.

3. Squash the blue handle or turn it clockwise to apply a vacuum to the balloon.

NOTES ON INFILATION AND DEFLECTION:

1. To inflate the balloon, squeeze the blue handle bar allowing the plunger to return to its resting position (0 ATM/BAR). Release grip on the luer handle bar, locking the plunger into position. To increase pressure, rotate handle clockwise until the desired pressure is achieved. Pressures above the maximum range will be indicated with flashing numbers.

The tick mark will return to the highest point of the last pressure reading. As the pressure decreases from the maximum pressure, the blue handle will begin to turn clockwise.

2. To deflate balloon, handle counter-clockwise to release pressure to 25 ATM or lower. Squeeze the blue handle bar to prevent a negative pressure. Release grip on the luer handle bar and pull back the plunger. Release grip on the luer handle to prevent a negative pressure. Release grip on the luer handle bar and pull back the plunger. Release grip on the luer handle to prevent a negative pressure.

3. The tick mark will reset to zero and the blue handle will turn clockwise to the original position.

4. To increase pressure, turn the blue handle clockwise to the original position.

5. To decrease pressure, turn the blue handle counter-clockwise to the original position.

6. To deflate balloon, handle counter-clockwise to release pressure to 25 ATM or lower. Squeeze the blue handle bar to prevent a negative pressure. Release grip on the luer handle bar and pull back the plunger. Release grip on the luer handle to prevent a negative pressure.

7. The tick mark will reset to zero and the blue handle will turn clockwise to the original position.

8. To increase pressure, turn the blue handle clockwise to the original position.

9. To decrease pressure, turn the blue handle counter-clockwise to the original position.

10. To deflate balloon, handle counter-clockwise to release pressure to 25 ATM or lower. Squeeze the blue handle bar to prevent a negative pressure. Release grip on the luer handle bar and pull back the plunger. Release grip on the luer handle to prevent a negative pressure.

11. The tick mark will reset to zero and the blue handle will turn clockwise to the original position.

12. To increase pressure, turn the blue handle clockwise to the original position.

13. To decrease pressure, turn the blue handle counter-clockwise to the original position.

14. To deflate balloon, handle counter-clockwise to release pressure to 25 ATM or lower. Squeeze the blue handle bar to prevent a negative pressure. Release grip on the luer handle bar and pull back the plunger. Release grip on the luer handle to prevent a negative pressure.

15. The tick mark will reset to zero and the blue handle will turn clockwise to the original position.

16. To increase pressure, turn the blue handle clockwise to the original position.

17. To decrease pressure, turn the blue handle counter-clockwise to the original position.

18. To deflate balloon, handle counter-clockwise to release pressure to 25 ATM or lower. Squeeze the blue handle bar to prevent a negative pressure. Release grip on the luer handle bar and pull back the plunger. Release grip on the luer handle to prevent a negative pressure.

19. The tick mark will reset to zero and the blue handle will turn clockwise to the original position.

20. To increase pressure, turn the blue handle clockwise to the original position.

21. To decrease pressure, turn the blue handle counter-clockwise to the original position.

22. To deflate balloon, handle counter-clockwise to release pressure to 25 ATM or lower. Squeeze the blue handle bar to prevent a negative pressure. Release grip on the luer handle bar and pull back the plunger. Release grip on the luer handle to prevent a negative pressure.

23. The tick mark will reset to zero and the blue handle will turn clockwise to the original position.

24. To increase pressure, turn the blue handle clockwise to the original position.

25. To decrease pressure, turn the blue handle counter-clockwise to the original position.

26. To deflate balloon, handle counter-clockwise to release pressure to 25 ATM or lower. Squeeze the blue handle bar to prevent a negative pressure. Release grip on the luer handle bar and pull back the plunger. Release grip on the luer handle to prevent a negative pressure.

27. The tick mark will reset to zero and the blue handle will turn clockwise to the original position.

28. To increase pressure, turn the blue handle clockwise to the original position.

29. To decrease pressure, turn the blue handle counter-clockwise to the original position.

30. To deflate balloon, handle counter-clockwise to release pressure to 25 ATM or lower. Squeeze the blue handle bar to prevent a negative pressure. Release grip on the luer handle bar and pull back the plunger. Release grip on the luer handle to prevent a negative pressure.

31. The tick mark will reset to zero and the blue handle will turn clockwise to the original position.

32. To increase pressure, turn the blue handle clockwise to the original position.

33. To decrease pressure, turn the blue handle counter-clockwise to the original position.

34. To deflate balloon, handle counter-clockwise to release pressure to 25 ATM or lower. Squeeze the blue handle bar to prevent a negative pressure. Release grip on the luer handle bar and pull back the plunger. Release grip on the luer handle to prevent a negative pressure.

35. The tick mark will reset to zero and the blue handle will turn clockwise to the original position.

36. To increase pressure, turn the blue handle clockwise to the original position.

37. To decrease pressure, turn the blue handle counter-clockwise to the original position.

38. To deflate balloon, handle counter-clockwise to release pressure to 25 ATM or lower. Squeeze the blue handle bar to prevent a negative pressure. Release grip on the luer handle bar and pull back the plunger. Release grip on the luer handle to prevent a negative pressure.

39. The tick mark will reset to zero and the blue handle will turn clockwise to the original position.

40. To increase pressure, turn the blue handle clockwise to the original position.

41. To decrease pressure, turn the blue handle counter-clockwise to the original position.

42. To deflate balloon, handle counter-clockwise to release pressure to 25 ATM or lower. Squeeze the blue handle bar to prevent a negative pressure. Release grip on the luer handle bar and pull back the plunger. Release grip on the luer handle to prevent a negative pressure.

43. The tick mark will reset to zero and the blue handle will turn clockwise to the original position.

44. To increase pressure, turn the blue handle clockwise to the original position.

45. To decrease pressure, turn the blue handle counter-clockwise to the original position.

46. To deflate balloon, handle counter-clockwise to release pressure to 25 ATM or lower. Squeeze the blue handle bar to prevent a negative pressure. Release grip on the luer handle bar and pull back the plunger. Release grip on the luer handle to prevent a negative pressure.

47. The tick mark will reset to zero and the blue handle will turn clockwise to the original position.

48. To increase pressure, turn the blue handle clockwise to the original position.

49. To decrease pressure, turn the blue handle counter-clockwise to the original position.

50. To deflate balloon, handle counter-clockwise to release pressure to 25 ATM or lower. Squeeze the blue handle bar to prevent a negative pressure. Release grip on the luer handle bar and pull back the plunger. Release grip on the luer handle to prevent a negative pressure.

51. The tick mark will reset to zero and the blue handle will turn clockwise to the original position.

52. To increase pressure, turn the blue handle clockwise to the original position.

53. To decrease pressure, turn the blue handle counter-clockwise to the original position.

54. To deflate balloon, handle counter-clockwise to release pressure to 25 ATM or lower. Squeeze the blue handle bar to prevent a negative pressure. Release grip on the luer handle bar and pull back the plunger. Release grip on the luer handle to prevent a negative pressure.

55. The tick mark will reset to zero and the blue handle will turn clockwise to the original position.

56. To increase pressure, turn the blue handle clockwise to the original position.

57. To decrease pressure, turn the blue handle counter-clockwise to the original position.

58. To deflate balloon, handle counter-clockwise to release pressure to 25 ATM or lower. Squeeze the blue handle bar to prevent a negative pressure. Release grip on the luer handle bar and pull back the plunger. Release grip on the luer handle to prevent a negative pressure.

59. The tick mark will reset to zero and the blue handle will turn clockwise to the original position.

60. To increase pressure, turn the blue handle clockwise to the original position.

61. To decrease pressure, turn the blue handle counter-clockwise to the original position.

62. To deflate balloon, handle counter-clockwise to release pressure to 25 ATM or lower. Squeeze the blue handle bar to prevent a negative pressure. Release grip on the luer handle bar and pull back the plunger. Release grip on the luer handle to prevent a negative pressure.

63. The tick mark will reset to zero and the blue handle will turn clockwise to the original position.

64. To increase pressure, turn the blue handle clockwise to the original position.

65. To decrease pressure, turn the blue handle counter-clockwise to the original position.

66. To deflate balloon, handle counter-clockwise to release pressure to 25 ATM or lower. Squeeze the blue handle bar to prevent a negative pressure. Release grip on the luer handle bar and pull back the plunger. Release grip on the luer handle to prevent a negative pressure.

67. The tick mark will reset to zero and the blue handle will turn clockwise to the original position.

68. To increase pressure, turn the blue handle clockwise to the original position.

69. To decrease pressure, turn the blue handle counter-clockwise to the original position.

70. To deflate balloon, handle counter-clockwise to release pressure to 25 ATM or lower. Squeeze the blue handle bar to prevent a negative pressure. Release grip on the luer handle bar and pull back the plunger. Release grip on the luer handle to prevent a negative pressure.

71. The tick mark will reset to zero and the blue handle will turn clockwise to the original position.

72. To increase pressure, turn the blue handle clockwise to the original position.

73. To decrease pressure, turn the blue handle counter-clockwise to the original position.

74. To deflate balloon, handle counter-clockwise to release pressure to 25 ATM or lower. Squeeze the blue handle bar to prevent a negative pressure. Release grip on the luer handle bar and pull back the plunger. Release grip on the luer handle to prevent a negative pressure.

75. The tick mark will reset to zero and the blue handle will turn clockwise to the original position.

76. To increase pressure, turn the blue handle clockwise to the original position.

77. To decrease pressure, turn the blue handle counter-clockwise to the original position.

78. To deflate balloon, handle counter-clockwise to release pressure to 25 ATM or lower. Squeeze the blue handle bar to prevent a negative pressure. Release grip on the luer handle bar and pull back the plunger. Release grip on the luer handle to prevent a negative pressure.

79. The tick mark will reset to zero and the blue handle will turn clockwise to the original position.

80. To increase pressure, turn the blue handle clockwise to the original position.

81. To decrease pressure, turn the blue handle counter-clockwise to the original position.

82. To deflate balloon, handle counter-clockwise to release pressure to 25 ATM or lower. Squeeze the blue handle bar to prevent a negative pressure. Release grip on the luer handle bar and pull back the plunger. Release grip on the luer handle to prevent a negative pressure.

83. The tick mark will reset to zero and the blue handle will turn clockwise to the original position.

84. To increase pressure, turn the blue handle clockwise to the original position.

85. To decrease pressure, turn the blue handle counter-clockwise to the original position.

86. To deflate balloon, handle counter-clockwise to release pressure to 25 ATM or lower. Squeeze the blue handle bar to prevent a negative pressure. Release grip on the luer handle bar and pull back the plunger. Release grip on the luer handle to prevent a negative pressure.

87. The tick mark will reset to zero and the blue handle will turn clockwise to the original position.

88. To increase pressure, turn the blue handle clockwise to the original position.

89. To decrease pressure, turn the blue handle counter-clockwise to the original position.

90. To deflate balloon, handle counter-clockwise to release pressure to 25 ATM or lower. Squeeze the blue handle bar to prevent a negative pressure. Release grip on the luer handle bar and pull back the plunger. Release grip on the luer handle to prevent a negative pressure.

91. The tick mark will reset to zero and the blue handle will turn clockwise to the original position.

92. To increase pressure, turn the blue handle clockwise to the original position.

93. To decrease pressure, turn the blue handle counter-clockwise to the original position.

94. To deflate balloon, handle counter-clockwise to release pressure to 25 ATM or lower. Squeeze the blue handle bar to prevent a negative pressure. Release grip on the luer handle bar and pull back the plunger. Release grip on the luer handle to prevent a negative pressure.

95. The tick mark will reset to zero and the blue handle will turn clockwise to the original position.