

REPLACEABLE PRESSURE INFUSOR BAG

INSTRUCTIONS FOR USE

DESCRIPTION OF THE PRODUCT

The pressure infusor bag (PIB) consists of a bag with a bladder, an envelope for holding the fluid bag, a pressure gauge, and tubing with an inflation bulb and a stopcock. Once the fluid bag is loaded in the envelope, the bag will then be inflated to the desired pressure.

INTENDED USE

A pressure infusor bag is used to apply pressure to a sealed bag of sterile fluid assisting in the infusion of the fluid.

CONTRAINDICATIONS

There are no Contraindications or Warnings identified for this device.

CAUTIONS

- Read instructions prior to use.
- Ensure adequate inflation pressure during use.
- **⚠ ONLY** Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

INSTRUCTIONS FOR USE

1. Insert intravenous bag between mesh and pressure bag.
2. Slide hanger loop through top of intravenous bag.
3. Inflate: Close valve on inflation bulb. Pump inflation bulb until gauge indicates the proper pressure.
4. Rotate stopcock to "Maintain Pressure" position.
5. Deflate: Open valve on inflation bulb.

Storage Conditions: Room temperature, ventilation, dry, away from light.

Transportation Conditions: Do not expose to excessive heat or humidity during transportation.

1. To deflate pressure bag turn the "off" handle on the stopcock/valve towards the bulb.

STORAGE CONDITIONS

Room temperature, ventilation, dry, away from light.

TRANSPORTATION CONDITIONS:

Do not expose to excessive heat or humidity during transportation.

	Caution
	Do not use if package is damaged and consult instruction for use
REF	Catalog number
LOT	Batch code
MD	Medical Device
UDI	Unique Device Identifier
	Consult Instructions for Use For electronic copy scan QR Code, or go to www.merit.com/ifu and enter IFU ID Number. For printed copy available within seven calendar days, call U.S.A. or EU Customer Service.
⚠ ONLY	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

	Use by date: YYYY-MM-DD
	Date of manufacture: YYYY-MM-DD
	Manufacturer
EC REP	Authorized Representative in European Community
	Contains DEHP
	Non-sterile



www.merit.com



Manufacturer:
Merit Medical Systems, Inc.
1600 West Merit Parkway,
South Jordan, Utah 84095 U.S.A.
1-801-253-1600
U.S.A Customer Service 1-800-356-3748



Authorized Representative:
Merit Medical Ireland Ltd,
Parkmore Business Park West
Galway, Ireland
European Customer Service +31 43 358 82 22