Description

A manifold with integral transducer is a pre-calibrated, single use device for physiological pressure measurement and fluid delivery. A separate reusable interface cable is used with this system to connect the transducer to a pressure monitor. A pressure monitor tubing line may be used as a compensation line with a pole mount organizer clip.

Indications

The Merit manifold with integrated transducer is used in cardiovascular diagnostic, surgical, and therapeutic applications to interconnect tubing, catheters, or other devices. Also used for measurement of blood pressure.

Rx Only. Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

Warnings

- Ensure that you are making secure connections when using this device to prevent the introduction of air into the system. All connections should be finger tightened. Over tightening can cause cracks and leaks to occur.

- Check for fluid leakage before and during the procedure. Leaks can result in the loss of sterility, fluid or blood loss, and/or air embolism. If a product leaks before or during use, retighten the leaking connection or replace the product.

- Prior to pressure injecting through the manifold, open the zeroing port. This product does not incorporate protection from accidental over pressurization. Over pressurizing may permanently impair the accuracy of the device.

Precautions

- Contents supplied sterile. Do not use if sterile barrier is damaged. If damage is found, contact a Merit Medical representative.

- Carefully read instructions before using product. If product is being used in conjunction with other manufacturers' components, also read instructions for use.

- Use proper aseptic techniques while handling product.

- The presence of air in the system will damp the transmission of the patient's pressure to the transducer. Be sure to eliminate all air bubbles within the system.

- The compensation line must be completely filled with fluid to produce accurate pressure readings.

- The female luer end of the compensation line (placed in the Merit organizer bracket) must be positioned at the mid axillary level in order to produce accurate readings.

- During fluid injection through the main lumen of the manifold, ensure proper orientation of the port handles so that fluid does not enter the side ports.

- Do not use yellow side port as an injection site for fluids.

- Inspect device prior to use to verify that no damage has occurred during shipping.

Instructions for Use

- Using aseptic technique, open the package containing the sterile product.

- Check all connections before use and finger tighten, if necessary. To prevent stripping, do not over tighten. Inspect for damage or improper assembly.

- Begin set-up according to hospital protocol for catheterization/pressure monitoring procedures. Purge the system of air bubbles.

- Ensure that all electrical connectors are dry. Connect the Meritrans disposable transducer cable to the reusable monitor cable. Align the connector pins, firmly join the connectors together.

- Attach the male end of the sterile compensation line to the transducer side port. Verify connection is tight but do not over tighten.

- Turn the manifold handles so the transducer is open to the fluid source. Open the yellow stopcock and flush the transducer free of air. Continue fluid filling through the yellow stopcock and out the compensation line until completely full. Debubble the compensation line. Turn the manifold handle and yellow stopcock off to the transducer.

Note: When filling the compensation line there will be an air space in the transducer housing between the yellow handle and the transducer element. This is normal and has no affect on pressure readings or frequency response.

- Mount a Merit organizer bracket on the pole. Place the compensation line mounting clip in the organizer bracket so that the slots of the plate are facing up and it snaps into the organizer. Place the female luer of the compensation line in the mounting clip. Adjust the bracket on the pole so that the top of the fluid column is at the mid axillary level.

- Inspect the entire length of the compensation line to ensure it is completely fluid filled and bubble free.

- Once the system has been fluid filled and the air is removed, the system is ready to be zero balanced. Ensure the yellow stopcock is open to atmosphere and the manifold transducer handle is closed to the patient.

- Balance and calibrate the system according to the monitor manufacturer's instructions.
- Once the system is balanced and has begun monitoring pressures, changes in elevation of the manifold will not affect the continued accuracy of pressure readings.

- To monitor patient pressure, ensure the yellow stopcock is off to atmosphere. Turn the manifold handles so that the transducer lumen is open to the catheter.

Note: Ensure all connections are securely tightened. Inspect carefully for air bubbles and flush the manifold lumen, if necessary.

- If the patient position is changed, adjust the height of the compensation mounting clip so that the top of the fluid column is maintained at mid axillary level.

REUSE PRECAUTION STATEMENT:

For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

MERITRANS SPECIFICATIONS			
Excitation Voltage	1.0 to 10 Vdc-5kHz	Excitation Impedance	240-350 Ω
Signal Impedance	300 Ω ±30Ω	Phase Shift	<5°
Sensitivity	5 μV/V/mmHg	Operating Temperature	15°C to 40°C
Storage Temperature	-25°C to 70°C	Maximum Half-Sine	
Zero Drift	1 mmHg/8 hours	Shock Acceleration	4500 G
Thermal Coefficient Offset	±0.3 mmHg/°C	Thermal Coefficient Span	±0.1%/°C
Light Sensitivity	<1 mmHg-		

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	
<p>Manufacturer: Merit Medical Systems, Inc. South Jordan, Utah 84095 U.S.A. 1-801-253-1600 www.merit.com U.S.A. Customer Service 1-800-356-3748</p> <p>Authorized Representative: Merit Medical Ireland Ltd, Parkmore Business Park West, Galway, Ireland European Customer Service by Country: Belgium 0800 72906; France 0800 916030; Germany 0800 1820871; Ireland 091 703733; Neth. 0800 0228184; U.K. 0800 973115</p>	

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