TRAM ®

TRAM  **P**™

EN	INSTRUCTIONS FOR USE
FR	MODE D'EMPLOI
IT	ISTRUZIONI PER L'USO
DE	GEBRAUCHSANWEISUNG
ES	INSTRUCCIONES DE USO
PT	INSTRUÇÕES DE UTILIZAÇÃO
PTBR	INSTRUÇÕES DE USO
NL	GEbruIKSAANWIJZING
SV	BRUKSANVISNING
NO	BRUKSANVISNING
DA	BRUGSANVISNING
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
TR	KULLANMA TALİMATLARI
RU	ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ
PL	INSTRUKCJA UŻYWANIA
CS	NÁVOD K POUŽITÍ
BG	ИНСТРУКЦИИ ЗА УПОТРЕБА
HU	HASZNÁLATI UTASÍTÁS
LV	LIETOŠANAS INSTRUKCIJA
LT	NAUDOJIMO INSTRUKCIJOS
RO	INSTRUCȚIUNI DE UTILIZARE
SK	NÁVOD NA POUŽITIE
HR	UPUTE ZA UPORABU
ET	KASUTUSJUHEND
SL	NAVODILA ZA UPORABO
FI	KÄYTTÖOHJEET

MERTMEDICAL®

INSTRUCTIONS FOR USE

DESCRIPTION

A manifold with integral transducer is a pre-calibrated, single use device for physiological pressure measurement and fluid delivery. A separate reusable interface cable is used with this system to connect the transducer to a pressure monitor.

A pressure monitor tubing line may be used as a compensation line with a pole mount organizer clip when using the TRAM-P.

INTENDED USE

The TRAM (Merit manifold with integrated transducer) and TRAM-P (Merit manifold with integrated transducer and compensation line) are used for cardiovascular diagnostic, surgical, and therapeutic applications to interconnect tubing, catheters, or other devices. They are also used for the measurement of blood pressure.

INDICATIONS

The TRAM and TRAM-P disposable pressure transducer is indicated for use in patients requiring invasive physiologic pressure monitoring.

CLINICAL BENEFITS

The indirect clinical benefit of the Merit TRAM and TRAM-P is physiologic pressure measurement to facilitate patient monitoring and/or medical management, as well as the control of fluid flow through an accessory device during procedures.

ONLY Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

USER/ PATIENT/ CLINICAL

User: Qualified nurses, clinicians, and physicians

Patient: Pediatric and adult applications

Clinical: Hospitals or appropriate clinical environments

PRECAUTIONS

- The manifold with integral transducer is an EO-sterilized, single use item. Do not use if sterile barrier is damaged. If damage is found, contact a Merit Medical representative.
- Carefully read instructions before using product.
- This product is to be used in conjunction with a monitor that is compliant to applicable IEC 60601 electrical safety standards. If this product is being used in conjunction with other manufacturers' components, their instructions for use must also be reviewed.
- Use proper aseptic techniques while handling product.
- The presence of air in the system will damp the transmission of the patient's pressure to the transducer. Be sure to eliminate all air bubbles within the system.
- During fluid injection through the main lumen of the manifold, ensure proper orientation of the port handles so that fluid does not enter the side ports.
- Do not use yellow side port as an injection site for fluids.
- Inspect device prior to use to verify that no damage has occurred during shipping.
- The device is to be kept away from sunlight.
- The device is to be kept dry from moisture.
- When using TRAM-P
 - The compensation line must be completely filled with fluid to produce accurate pressure readings.
 - The female luer end of the compensation line (placed in the Merit organizer bracket) must be positioned at the mid axillary level in order to produce accurate readings.

CONTRAINDICATIONS

None known.

REUSE PRECAUTION STATEMENT

For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

POTENTIAL COMPLICATIONS

- Abnormal Pressure Readings
- Air Emboli
- Clotted Catheter and Bleed-Back
- Inflammation
- Overinfusion
- Sepsis/Infection
- Soft Tissue Injury

INSTRUCTIONS FOR USE

- Using aseptic technique, open the package containing the sterile product.
- Check all connections before use and finger tighten, if necessary. To prevent stripping, do not over tighten. Inspect for damage or improper assembly.
- Begin set-up according to hospital protocol for catheterization/pressure monitoring procedures. Purge the system of air bubbles.
- Ensure that all electrical connectors are dry. Connect the transducer cable to the reusable monitor cable. Align the connector pins, firmly join the connectors together.
- TRAM**
 - To complete the set-up, open the stopcock to atmosphere and flush the transducer port free of air.
 - Once the entire system has been fluid filled and the air is removed, the system is ready to be zero balanced.
 - Ensure the yellow stopcock is open to air and the manifold transducer handle is closed to the patient.
 - Balance and calibrate the system according to the monitor manufacturer's instructions.
 - To monitor pressure, orient the port handles so that the transducer lumen is open towards the catheter or patient. (Inspect carefully for air bubbles and reflush, if necessary.)
- Note:** Ensure all connections are securely tightened. Inspect carefully for air bubbles and flush the manifold lumen, if necessary.
- The manifold should be returned to the original position used to zero balance the system. The reference point assures consistency as pressures can be affected by elevating the transducer.

6. TRAM-P

- Attach the male end of the sterile compensation line to the transducer side port. Verify connection is tight but do not over tighten.
- Turn the manifold handles so the transducer is open to the fluid source. Open the yellow stopcock and flush the transducer free of air.
- Continue fluid filling through the yellow stopcock and out the compensation line until completely full.
- Debubble the compensation line. Turn the manifold handle and yellow stopcock off to the transducer.
- Note:** When filling the compensation line there will be an air space in the transducer housing between the yellow handle and the transducer element. This is normal and has no effect on pressure readings or frequency response.
- Mount a Merit organizer bracket on the pole. Place the compensation line mounting clip in the organizer bracket so that the slots of the plate are facing up and it snaps into the organizer. Place the female luer of the compensation line in the mounting clip. Adjust the bracket on the pole so that the top of the fluid column is at the mid axillary level.
- Inspect the entire length of the compensation line to ensure it is completely fluid filled and bubble free.
- Once the system has been fluid filled and the air is removed, the system is ready to be zero balanced. Ensure the yellow stopcock is open to atmosphere and the manifold transducer handle is closed to the patient.
- Balance and calibrate the system according to the monitor manufacturer's instructions.
- Once the system is balanced and has begun monitoring pressures, changes in elevation of the manifold will not affect the continued accuracy of pressure readings.
- To monitor patient pressure, ensure the yellow stopcock is off to atmosphere. Turn the manifold handles so that the transducer lumen is open to the catheter.
- Note:** Ensure all connections are securely tightened. Inspect carefully for air bubbles and flush the manifold lumen, if necessary.
- If the patient position is changed, adjust the height of the compensation mounting clip so that the top of the fluid column is maintained at mid axillary level.

WARNINGS

- Ensure that you are making secure connections when using this device to prevent the introduction of air into the system. All connections should be finger tightened. Over tightening can cause cracks and leaks to occur.
- Check for fluid leakage before and during the procedure. Leaks can result in the loss of sterility, fluid or blood loss, and/or air embolism. If a product leaks before or during use, retighten the leaking connection or replace the product.
- Prior to pressure injecting through the manifold, open the zeroing port. This product does not incorporate protection from accidental over pressurization. Over pressurizing may permanently impair the accuracy of the device.

MR Conditional

Non-clinical testing has demonstrated the Merittrans portion of the device is MR Conditional. A patient with this device can be safely scanned in an MR system under the following conditions:

- Static magnetic field of 3-Tesla or less
- Maximum spatial gradient magnetic field of 4,000 Gauss/cm (40 T/m)
- The transducer and cable must not contact the patient during operation of the MR system
- The transducer and cable must not be placed inside the bore of the MR system during operation of the scanner
- The transducer is permitted in the MR system room, but must not be operational or connected to a monitoring system during the MRI procedure.

SAFE DISPOSAL

After use, dispose of device in a manner consistent with standard protocols for waste disposal.

STORAGE CONDITIONS

Store in a cool, dark and dry place.

TRANSPORTATION CONDITIONS

Do not expose to excessive heat or humidity during transportation.

In the EU, any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the applicable Member State.

MERITRANS SPECIFICATIONS	
Excitation Voltage	4-8 Vdc-5kHz
Excitation Impedance	240-350 Ω
Signal Impedance	300 Ω ±30 Ω
Phase Shift	<5°
Sensitivity	5µV/V/mmHg
Operation Temperature	15°C to 40°C
Storage Temperature	-25°C to 70°C
Shock Acceleration	4500 G
Zero Drift	1 mmHg//4 hours
Thermal Coefficient Offset	±0.3 mmHg/°C
Thermal Coefficient Span	±0.1%/°C
Light Sensitivity	<1 mmHg at 3000 ft-Candles

SYMBOL	DESIGNATION
	Catalog Number
	Lot Number
	Caution
	Consult Instructions for Use. For electronic copy scan QR Code, or go to www.merit.com/ifu and enter IFU ID Number. For printed copy available within seven days, call U.S.A or E.U. Customer Service
	Sterilized using Ethylene Oxide
	Do not use if package is damaged and Consult instructions For Use
	Single Use
	Use By
	Do Not Resterilize
	Non-pyrogenic
	Defibrillation-proof type CF applied part
	MR Conditional
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Single Sterile Barrier System
	Keep away from Sunlight
	Keep dry
	Manufacturer
	Authorized Representative in European Community
	Medical Device
	Date of Manufacture
	Unique Device Identifier

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Tram is a registered and unregistered trademark of Merit Medical Systems, Inc.
Merit Medical Systems, Inc.
U.S.A. 1-801-253-1600 U.S.A. Customer Service 1-800-356-3748
Merit Medical Ireland Ltd, Parkmore Business Park