FLOW CONTROL SWITCH

INSTRUCTIONS FOR USE

INTENDED USE

The flow control switch is an angiographic accessory intended for use as an ON/OFF device for angiography and other high pressure applications.

DEVICE DESCRIPTION

The flow control switch is a high pressure angiographic device with a male luer lock fitting at one end, a threaded female luer lock fitting at the other end, and a sliding, top-mounted ON/OFF switch.

The device is a one-handed actuation locking shutoff of fluid flow when installed in line. Luer lock fittings on both ends allow for connection to catheters and tubes. The front (male luer) portion of the device is connected to a catheter or other applicable device. The rear (female luer) portion of the device is connected to a compatible luer lock connecting device.

CONTRAINDICATIONS

None known.

PRECAUTIONS

- Before using this device, there should be a thorough understanding of the technical principles, clinical applications and risks associated with angiography and percutaneous interventional procedures.
- Prior to use carefully examine the unit to verify the device or sterile package has not been damaged in shipment.
- Do not exceed 1,200psi/82bar/8274kPa in the ON position and 150psi/10bar/1034kPa in the OFF position.

POTENTIAL COMPLICATIONS

| | |
|----------------------------|---|
| Arterial/Venous Thrombosis | Allergic Reaction (including anaphylaxis) |
| Cerebral Vascular Accident | Death |
| Embolism | Exposure to Biohazards |
| Hemorrhage | Infection |
| Myocardial Infarction | Transient Ischemic Attack (TIA) |

INSTRUCTIONS FOR USE

CONNECT DEVICE AS FOLLOWS

- Front portion (male luer) to catheter or other applicable device.
- Rear portion (female luer) to a compatible luer lock connecting device.
- Ensure all connections are tight.
- Place device in desired ON/OFF position.

REUSE PRECAUTION STATEMENT

For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

Single Use

Caution: Consult accompanying document

| | |
|----------------|-----------|
| STERILE | EO |
|----------------|-----------|

Non-pyrogenic
Sterile if package is unopened or undamaged

Rx Only: Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

www.merit.com

Manufacturer:
Merit Medical Systems, Inc. 1600 West Merit Pkwy, South Jordan, Utah 84095
U.S.A. 1-801-253-1600 U.S.A. Customer Service 1-800-356-3748

| | |
|-----------|------------|
| EC | REP |
|-----------|------------|

Authorized Representative:
Merit Medical Ireland Ltd, Parkmore Business Park West, Galway, Ireland

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