INSTRUCTIONS FOR USE


PRODUCT NAME

ReSolve® NL (Non-Locking) Drainage Catheter with Hydrophilic Coating

DESCRIPTION OF THE PRODUCT

The ReSolve NL Drainage Catheter with hydrophilic coating is a radiopaque catheter with multiple side holes used for percutaneous drainage. The components of the catheter allow for introduction and placement using a trocar stylette or over-the-wire method. The ReSolve NL Drainage Catheter may be packaged in a pouch with the following components:

- One (1) ReSolve NL Drainage Catheter with hydrophilic coating (A)
- One (1) Metal stiffening cannula (B)
- One (1) Trocar stylette (C)

INDICATIONS FOR USE

The ReSolve NL Drainage Catheter with hydrophilic coating is intended for percutaneous drainage of fluid from body cavities.

CONTRAINDICATIONS

The ReSolve NL Drainage Catheter is contraindicated for use where percutaneous drainage catheterization is unacceptable.

The ReSolve NL Drainage Catheter is contraindicated for intravascular use.

RESTRICTED DEVICE

Rx Only; Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician trained and/or experienced in the use of this device.

PRECAUTIONS

- Read manufacturer's instructions prior to use.
- Contents are sterile (via ethylene oxide) and non-pyrogenic.
- Do not use if packaging is opened, damaged or broken.
- For single patient use only. Do not reuse, reprocess, or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure, which in turn may result in patient injury, illness, or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.
- Follow universal precautions when inserting and maintaining this device. Due to the risk of bloodborne pathogens, healthcare professionals should always use standard blood and body fluid precautions in the care of all patients. Sterile technique should always be followed.
- Do not use after expiration date.
- Store in a cool, dry place.

WARNING: The ReSolve NL Drainage Catheter is not to be used to deliver nutritional supplements.

ADVERSE REACTIONS*

- Septic shock
- Bacteremia
- Hemorrhage
- Superinfection
- Bowel transgression
- Pleural transgression
- Vascular injury
- Pneumothorax
- Skin Infection
- Catheter Occlusion
- Catheter Dislodgement

*Brountzos EN. Quality improvement guidelines for percutaneous nephrostomies.[Internet]. 2006. CIRSE.org. http://cirse.org/files/File/05_qig.pdf. Accessed 04/11/13.

Wallace MJ, et al. Quality improvement guidelines for percutaneous drainage/aspiration of abscess and fluid collections. J Vasc Interv Radiol. 2010;21:431-435.

Ramchandani P, et al. Quality improvement guidelines for percutaneous nephrostomy. J Vasc Interv Radiol. 2003;14:S277-S281.

INSTRUCTIONS FOR USE

PLACEMENT TECHNIQUES OPTION 1: DIRECT PUNCTURE USING TROCAR STYLET

1. Remove the stiffening cannula and trocar stylette assembly from the ReSolve NL Drainage Catheter.

2. Ensure that the distal portion of the catheter is wet prior to placement.

CAUTION: To maximize the advantages of the hydrophilic coating on the surface of the distal portion of the catheter, wet the catheter prior to use with sterile water or saline. Keep catheter wet during placement.

WARNING: DO NOT wipe catheter with dry gauze or any solvents because it may damage the catheter coating.

3. Flush catheter prior to use.

4. Place the metal stiffening cannula into the catheter and tighten the Luer lock fittings. See Figure 1.

Figure 1

5. Remove the paper spacer from the trocar stylette. Advance the trocar stylette through the metal stiffening cannula and tighten the Luer lock fittings. See Figure 2.

Figure 2

6. Place the catheter/cannula/trocar assembly into the fluid collection site using standard insertion technique. Placement should be confirmed with diagnostic imaging.

7. After placement is confirmed, remove the trocar stylette and stiffening cannula.

8. The ReSolve NL Drainage Catheter is now ready to be connected to appropriate drainage bag or tubing.

9. A flush regimen should be designed for the circumstances of each patient and the protocol of the physician.

Note: Instruct patient or other healthcare personnel in appropriate device function and/ or maintenance.

INSTRUCTIONS FOR USE

PLACEMENT TECHNIQUES OPTION 2: SELDINGER ENTRY TECHNIQUE OR GUIDE WIRE EXCHANGE

1. Remove the stiffening cannula and trocar stylette assembly from catheter.

2. Ensure that the distal portion of the catheter is wet prior to placement.

CAUTION: To maximize the advantages of the hydrophilic coating on the surface of the distal portion of the catheter, wet the catheter prior to use with sterile water or saline. Keep catheter wet during placement.

WARNING: DO NOT wipe catheter with dry gauze or any solvents because it may damage the catheter coating.

3. Flush catheter prior to use.

4. Place the metal stiffening cannula into the catheter and tighten the Luer lock fittings. See Figure 3.


Figure 3

5. Place catheter over appropriate guide wire into the fluid collection site. The catheter accommodates an 0.038" (0.97 mm) wire. See Figure 4.

Figure 4

6. Advance the catheter into the site while holding the guide wire in place. Placement should be confirmed with diagnostic imaging. After placement is confirmed, remove cannula and guide wire.

7. The catheter is now ready to be connected to appropriate drainage bag or tubing.

8. A flush regimen should be designed for the circumstances of each patient and the protocol of the physician.

Note: Instruct patient or other healthcare personnel in appropriate device function and/ or maintenance.

CATHETER EXCHANGE OR REMOVAL

1. Disconnect catheter from drainage tubing or bag.

2. For catheter exchange or if access is to be maintained, advance appropriate guide wire through catheter; use diagnostic imaging to confirm wire placement. Guide wire will maintain access to drainage site, if necessary.

3. Carefully remove the catheter. Proceed with either catheter exchange or skin closure.

Storage Conditions: Do not expose to excessive heat or humidity.

Transportation Conditions: Do not expose to excessive heat or humidity during transportation.

Shelf Life: 3 years.

Body Contact Time: Less than 30 days.

Models:

RNL-6-038J, RNL-6-038S, RNL-8-038J, RNL-8-038S, RNL-10-038J, RNL-10-038S, RNL-12-038J, RNL-12-038S

Manufacturer: Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095 U.S.A. 1-801-253-1600
U.S.A. Customer Service 1-800-356-3748

Authorized Representative: Merit Medical Ireland Ltd,
Parkmore Business Park West, Galway, Ireland
www.merit.com

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