English

GUIDE WIRE INSERTION DEVICE

Step 1: Place the distal end of device (dilator, introducer, catheter, etc.) through the BowTie™. Device should not go past the mid part of the BowTie.

Step 2: Advance the BowTie, with the device, over the proximal end of the guide wire.

Step 3: When the guide wire is in the device, remove the BowTie by pulling down on the tab. The slit on the top of the BowTie facilitates easy removal.

INTENDED USE: The BowTie guide wire insertion device is used to facilitate placement of the proximal end of a micro guide wire $\leq 0.018''$ into the tip of another device such as a dilator, introducer or catheter.

Sterile if package is unopened and undamaged.

Non-pyrogenic / R_x only

REUSE PRECAUTION STATEMENT

For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

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Merit Medical Systems, Inc.
South Jordan, Utah 84095 U.S.A. 1-801-253-1600
U.S.A. Customer Service 1-800-356-3748
Merit Medical Ireland Ltd,