
INSTRUCTIONS FOR USE
MODE D'EMPLOI
ISTRUZIONI PER L'USO
GEBRAUCHSANWEISUNG
INSTRUCCIONES DE USO
INSTRUÇÕES DE UTILIZAÇÃO
GEBRUIKSAANWIJZING
BRUKSANVISNING
BRUGSANVISNING
ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
KULLANMA TALIMATLARI
УКАЗАНИЯ ПО
ПРИМЕНЕНИЮ

FLEX-NECK® CATHETER
EXTERNAL REPAIR KIT
CE-1400

For Flex-Neck® Peritoneal Dialysis Catheters ONLY


FLEX-NECK® CATHETER EXTERNAL REPAIR KIT CE-1400

English

For Flex-Neck Peritoneal Dialysis Catheters ONLY

PRODUCT DESCRIPTION

- Tube-to-Tube Connector for all Flex-Neck Adult, Adolescent, and Pediatric Catheters- Sterile Pouch (A)
- Tube-to-Tube Connector for all Flex-Neck Infant Catheters - Sterile Pouch (B)
- 20 cm of Repair Catheter Tubing (C)
- Catheter Connector (D)
- Cap (E)

INDICATIONS FOR USE

This device is intended for use with Flex-Neck catheters and Ash Advantage™ catheters in the event one of the two situations described below should occur:

1. Extender Application

In some patients, the external portion of the PD catheter may be too short for them to comfortably use. For example, obese patients may need an extender kit to be able to work with external part of the catheter. Other patients may require the external end of the catheter to be more centrally located than it normally might be.

2. Repair Application

If a peritoneal dialysis catheter develops a fracture in the external portion, the catheter must be trimmed to eliminate the fracture. Doing so, by definition, will shorten the catheter and may make it impractical or difficult or even unsuitable for use. Rather than replace the entire catheter, it is desirable to retain the catheter by extending the external part of the catheter with new tubing.

In either situation, this device enables qualified medical staff to extend and/or repair the external part of an existing Flex-Neck catheter or Ash Advantage catheter by adding a 20 cm long tube to the existing and/or remaining catheter tubing.

CONTRAINDICATIONS FOR USE

Do NOT use with other peritoneal dialysis catheter brands.

Rx Only; CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

PRECAUTIONS



- Read manufacturer's instructions prior to use.
- Contents are sterile (via ethylene oxide). Do not use if packaging is opened, damaged or broken.
- For single patient use only. Do not reuse, reprocess, or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure, which in turn may result in patient injury, illness, or death. Reuse, reprocessing, or

resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

- Do not use after expiration date.
- The medical techniques, procedures and potential complications stated herein do NOT give full and/or complete coverage or descriptions. They are not a substitute for adequate training and sound medical judgment by a physician.
- Use an aseptic procedure to open the package and to remove the contents.

POTENTIAL COMPLICATIONS

- Infections (exit-site or tunnel)
- Peritonitis
- Sepsis
- Leakage (initial or latent)
- Catheter separation from the connector

CAUTIONS

- Follow the specific order of connections as given in the Instructions for Use. Do not connect the various tubings and connectors in any other order.
- Use only sterile saline or sterile water as lubricants. Do not use any other lubricants.
- Catheter tubing can tear when subjected to repeated clamping, serrated-jaw forceps, excessive force, or rough tools.
- Do NOT use forceps with a serrated jaw.
- Do NOT use excessive force to lock the forceps closed.
- Use ONLY smooth-jawed forceps or equivalent.
- Do NOT clamp the catheter, or repair tubing repeatedly in the same area.
- Do NOT clamp near the connector.

INSTRUCTIONS FOR USE

1. Create a sterile field for the procedure:
 - Prepare the patient's skin and existing catheter with a cleaning or disinfecting solution as per hospital protocol.
 - Drape the catheter exit-site area in an appropriate manner.
2. If clamping the existing catheter or Repair Catheter Tubing (C) with sterile hemostats or forceps, use only smooth-jawed forceps or equivalent. Do not use forceps with a serrated jaw; see cautions.
3. Prepare the existing Flex-Neck Peritoneal Dialysis Catheter:
 - Cut catheter tubing just distal to the existing connector with sterile suture scissors, in a single, straight, perpendicular cut. Do not pull the connector off the catheter.
 - There must be at least 2.5 cm of catheter tubing remaining after trimming.
 - Verify that the cut is perpendicular to the tubing.
4. Select the appropriate Tube-to-Tube Connector, based on the size of the existing catheter:
 - Tube-to-Tube Connector (A) for all Flex-Neck Adult, Adolescent and Pediatric Catheters
 - Tube-to-Tube Connector (B) for all Flex-Neck Infant Catheter
5. Flush and rinse the Repair Kit components with sterile saline.
6. Wet the tapered tip of the appropriate Tube-to-Tube Connector (A or B) with sterile saline or sterile water and insert it into Repair Catheter Tubing (C).
 - Do not use a twisting motion to force the Tubing onto the Connector. Push the Connector into the Tubing with a single forward motion.
 - If Infant Tube-to-Tube Connector (B) is used, insert the wider end of the Connector into Repair Catheter Tubing (C). See Figure 2.
7. Advance the Repair Catheter Tubing (C) completely to the raised ring barrier of Tube-to-Tube Connector (A) or (B). See Figure 1 or 2, as appropriate. Pull carefully on Tube-to-Tube Connector (A or B) and Repair Catheter Tubing (C) to test the strength of the junction.
8. Wet the tapered tip of the Catheter Connector (D) with sterile saline or sterile water, and insert it into the other end of the Repair Catheter Tubing (C). Advance the tube completely to the Connector's raised shoulder ridge. The catheter tubing must completely pass over the Tapered Tip and to the Raised Shoulder Ridge, but not beyond that onto the Finger Grip. See Figures 3 and 4. Pull carefully on the Catheter Connector and Repair Tube to test the strength of the junction.
 - Do not use a twisting motion to force the Tubing onto the Connector. Push the Connector into the Tubing with a single forward motion.
9. Wet the remaining tapered tip of the Tube-to-Tube Connector (A or B) in sterile saline or sterile water and carefully insert it into the original, existing catheter. Advance the catheter end completely to the raised ring barrier of the Connector. See Figure 1 or 2, as appropriate. Pull carefully on the catheter and connector to test the strength of the junction.
 - **CAUTION:** Do NOT dislodge the cuffs.
 - Do not use a twisting motion to force the catheter onto the Tube-to-Tube connector. Push the Tube-to-Tube connector into the catheter with a single forward motion.

10. Attach either Cap (E), or a dialysis transfer set, to the threaded luer end of Catheter Connector (D).

CATHETER CLEANING AND CARE

All Flex-Neck Peritoneal Dialysis Catheters are made of silicone. Exit-site cleaning agents that are compatible with silicone catheters therefore may be acceptable for use on Flex-Neck Peritoneal Dialysis Catheters and repaired catheters. Such cleaning agents include:

- Electrolytically-produced sodium hypochlorite solutions (i.e., ExSept Plus[®])
- Normal (sterile) saline

Cleaning agents that are non-irritating, non-toxic, antibacterial, and in liquid form are generally recommended.

The following cleaning agents are not compatible with silicone catheters, and are not recommended for use with Flex-Neck Peritoneal Dialysis catheters:

- Acetone or acetone-based products
- Povidone-iodine or iodine-based products

Merit Medical Systems, does not provide specific recommendations or protocols for exit-site care and cleaning, whether by the healthcare professional or by the patient. Appropriate exit-site and catheter care treatment protocols should be individualized for each patient, and established by the patient's physician(s), nurse(s), dialysis center(s), and/or other relevant dialysis healthcare professionals.

ExSept Plus is a registered trademark of Alcavis HDC, LLC or one of its affiliates.

Copyright © Merit Medical Systems, Inc. All rights reserved.

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Copyright © Merit Medical Systems, Inc. Tous droits réservés.
Copyright © Merit Medical Systems, Inc. Tutti i diritti riservati.
Copyright © Merit Medical Systems, Inc. Alle Rechte vorbehalten.
Copyright © Merit Medical Systems, Inc. Reservados todos los derechos.
Copyright © Merit Medical Systems, Inc. Todos os direitos reservados.
Copyright © Merit Medical Systems, Inc. Alle rechten voorbehouden.