INSTRUCTIONS FOR USE
MODE D'EMPLOI
ISTRUZIONI PER L'USO
GEBRAUCHSANWEISUNG
INSTRUCCIONES DE USO
INSTRUÇÕES DE UTILIZAÇÃO
GEBRUIKSAANWIJZING
BRUKSANVISNING
BRUGSANVISNING
ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
KULLANMA TALİMATLARI
УКАЗАНИЯ ПО
ПРИМЕНЕНИЮ

Two-Piece
TITANIUM CATHETER
CONNECTOR
CC-2300

For use with Flex-Neck® Adult, Adolescent,
Pediatric and ExxTended™ Catheters ONLY

TWO-PIECE TITANIUM CATHETER CONNECTOR CC-2300

English

For use with Flex-Neck® Adult, Adolescent, Pediatric and ExxTended Catheters ONLY

PRODUCT DESCRIPTION

Each box contains one (1) pouch with two piece Titanium Connector. Includes main connector body and cap.

INDICATIONS FOR USE

The two-piece titanium peritoneal dialysis catheter connector is used to connect the peritoneal dialysis catheter to the dialysate line.

CONTRAINDICATIONS FOR USE

Do NOT use with Flex-Neck Infant Peritoneal Dialysis Catheters. Do NOT use with other peritoneal dialysis catheter brands.

MRI Safety Information:

MR Conditional

This product has been evaluated to determine MRI compatibility and have been found to be compatible. Non-clinical testing has demonstrated that the Peritoneal Dialysis Catheter is MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static Magnetic Field of 1.5T and 3.0T only
- Maximum spatial gradient magnetic field of 3,000 Gauss/cm (30T/m) or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode)

MRI-Related Heating

Under the scan conditions defined above, the Peritoneal Dialysis Catheter is expected to produce a maximum temperature rise of 1.5°C after 15 minutes of continuous scanning.

Artifact Information

In a non-clinical testing, the image artifact caused by the device extends approximately 10mm from the Peritoneal Dialysis Catheter when imaged with a gradient echo pulse sequence and a 3Tesla MRI system. The artifact does not obscure the device lumen.

Rx Only: CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

PRECAUTIONS

- Read manufacturer's instructions prior to use.
- Contents are sterile (via ethylene oxide). Do not use if packaging is opened, damaged, or broken.
- For single patient use only. Do not reuse, reprocess, or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure, which in turn may result in patient injury, illness, or death. Reuse, reprocessing or resterilization may also create a risk of contamination of

the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.


- Do not use after expiration date.
- The medical techniques, procedures and potential complications stated herein do NOT give full and/or complete coverage or descriptions. They are not a substitute for adequate training and sound medical judgment by a physician.
- Use an aseptic procedure to open the package and to remove the contents.

POTENTIAL COMPLICATIONS

- Infections (exit-site or tunnel)
- Peritonitis
- Sepsis
- Leakage (initial or latent)
- Catheter separation from the connector

CAUTIONS

- Catheter tubing can tear when subjected to repeated clamping, serrated-jaw forceps, excessive force, or rough tools.
- Do NOT use forceps with a serrated jaw.
- Do NOT use excessive force to lock the forceps closed.
- Use ONLY smooth-jawed forceps or equivalent.
- Do NOT clamp the catheter repeatedly in the same area.
- Do NOT clamp near the connector.

INSTRUCTIONS FOR USE

1. Create a sterile field for the procedure: Prepare the patient's skin and catheter with a cleaning or disinfecting solution as needed, as per hospital protocol. Drape the catheter exit-site area in an appropriate manner.
2. If clamping the catheter tubing with sterile hemostats or forceps, use only smooth-jawed forceps or equivalent. Do not use forceps with a serrated jaw. See cautions.
3. If replacing an existing connector, prepare the existing Flex-Neck Catheter:
 - Cut catheter tubing just distal to the existing connector with sterile suture scissors, in a single, straight, perpendicular cut. Do not pull the connector off the catheter.
 - There must be at least 2.5 cm of catheter tubing remaining after trimming.
 - Verify that the cut is perpendicular to the tubing.
4. Slide the Titanium Connector Cap over the catheter, with the tapered end first.
5. Insert the tip of the connector's Main Body into the catheter and advance until only the knob (A) of the Main Body is covered by the catheter as indicated by the arrow in diagram.
 - Do NOT advance past the knob (A) onto the neck (B).
 - If necessary, wet the tapered tip (A) of the Connector with sterile saline or sterile water before inserting it into the catheter.
 - Do not use any other lubricant.

- Do not use a twisting motion to force the catheter onto the Connector. Push the Connector into the catheter with a single forward motion.
6. Slide the Cap up and onto the main body and screw the two parts snugly together by rotating the cap. There must be no gap between the two parts of the connector, the Cap and the Main Body.
- Do NOT rotate the Main Body of the connector; the catheter may become twisted.
7. Pull carefully on the Replacement Connector and catheter to test the strength of the connection.
8. Attach either a connector cap, or a dialysis transfer set, to the threaded luer end of the Main Body.

CATHETER CLEANING AND CARE

All Flex-Neck Peritoneal Dialysis Catheters are made of silicone. Exit-site cleaning agents that are compatible with silicone catheters therefore may be acceptable for use on Flex-Neck Peritoneal Dialysis Catheters. Such cleaning agents include:

- Electrolytically-produced sodium hypochlorite solutions (i.e., ExSept Plus®)
- Normal (sterile) saline

Cleaning agents that are non-irritating, non-toxic, anti-bacterial, and in liquid form are generally recommended.

The following cleaning agents are not compatible with silicone catheters, and are not recommended for use with Flex-Neck Peritoneal Dialysis catheters:

- Acetone or acetone-based products
- Povidone-iodine or iodine-based products

Merit Medical Systems, does not provide specific recommendations or protocols for exit-site care and cleaning, whether by the healthcare professional or by the patient. Appropriate exit-site and catheter care treatment protocols should be individualized for each patient, and established by the patient's physician(s), nurse(s), dialysis center(s), and/or other relevant dialysis healthcare professionals.

ExSept Plus is a registered trademark of Alcavis HDC, LLC or one of its affiliates.

Copyright © Merit Medical Systems, Inc. All rights reserved.

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Copyright © Merit Medical Systems, Inc. Tous droits réservés.
Copyright © Merit Medical Systems, Inc. Tutti i diritti riservati.
Copyright © Merit Medical Systems, Inc. Alle Rechte vorbehalten.
Copyright © Merit Medical Systems, Inc. Reservados todos los derechos.
Copyright © Merit Medical Systems, Inc. Todos os direitos reservados.
Copyright © Merit Medical Systems, Inc. Alle rechten voorbehouden.