English

Czech

7.3F Hemostasis valve**7,3 F Hemostatický ventil****INSTRUCTIONS FOR USE****DESCRIPTION**

The PhD Hemostasis Valve™ has a 0.096" (2.44 mm) inside diameter. This device has two seals that operate independently: the compression seal and the secondary seal. The compression seal can be opened or closed by rotating the cap. Closing of the compression seal allows for pressure injections up to 400 psi and secures the diagnostic/interventional device in position within the vasculature.

The secondary seal is a diaphragm seal that forms around the diagnostic/interventional devices as they move into and out of the vasculature. This seal provides minimal fluid loss without restricting device movement. The secondary seal is open when the cap is pressed down, and closed when the cap is released. An open secondary seal allows air to be purged from the device and allows insertion/withdrawal of diagnostic/interventional devices.

INTENDED USE

The PhD Hemostasis Valve is intended to maintain hemostasis during the introduction/withdrawal and use of diagnostic and interventional devices.

INDICATIONS FOR USE

The PhD Hemostasis Valve is indicated for maintaining a seal around diagnostic/interventional devices with an outside diameter of ≤0.096" (2.44 mm) during interventional procedures.

CONTRAINDICATIONS

None known.

RX Only - Caution: Federal (USA) law restricts this device to sale by or on the order of a physician trained in angiography and percutaneous transluminal coronary angioplasty (PTCA), and/or percutaneous transluminal angioplasty (PTA).

PRECAUTIONS

- Read manufacturer's instructions for the use of catheters, guide wires, and introducers prior to use.
- Inspect the PhD Hemostasis Valve prior to use for any damage. Do not use if packaging is opened, damaged, or broken.
- Follow local guidelines and practices regulating the disposal of the PhD Hemostasis Valve after use.
- Contents are sterile (via ethylene oxide) and non-pyrogenic.

REUSE PRECAUTION STATEMENT

For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure, which in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

WARNING

- Do not inject any fluid if air bubbles are visible within the PhD Hemostasis Valve.
- To prevent damage to a balloon, it must be completely deflated during advancement/withdrawal through the PhD Hemostasis Valve.
- Failure to depress the cap prior to inserting a guide wire introducer could cause seal damage resulting in leakage and/or particulate embolization.
- Failure to open the secondary seal and compression seal prior to inserting or withdrawing a diagnostic/interventional device could cause damage to the device.
- Care should be taken to avoid over-tightening of the compression seal, as this may compromise the integrity of the diagnostic/interventional device lumen.
- Power injection at pressures greater than 400psi could result in leakage or detachment of components.

INSTRUCTIONS FOR USE:

1. Connect the sideport of the PhD Hemostasis Valve to the manifold assembly. With the compression valve open, flush and fill the assembly with saline. To fill the valve section, ensure the secondary seal is open by pushing down on the cap. Place one finger over the male Luer fitting and continue to fill the assembly until completely purged of air. Gently release the back end cap to engage the secondary seal.
2. Connect the rotating male Luer adapter of the PhD Hemostasis Valve to the guiding catheter. Ensure the connection is airtight.
3. To open the PhD Hemostasis Valve for catheter or guide wire insertion, push on the back end cap to open the secondary seal. The PhD Hemostasis Valve is now ready to accept any diagnostic/interventional device.
4. Insert the diagnostic/interventional device through an open secondary seal and an open compression seal.
CAUTION: PhD Hemostasis Valve will be fully open and allow for blood loss when secondary seal is open.
5. Advance the diagnostic/interventional device through an open compression seal. A fluid-tight seal, providing little to no blood loss, will form around the diagnostic or interventional device while not restricting movement of the device.
6. Once the diagnostic/interventional device is in place, if desired, close the compression seal to allow for pressure injections of up to 400 psi and/or secure the diagnostic/interventional device in position.
7. Remove the diagnostic/interventional device through an open compression seal. If resistance is felt, open the secondary seal as well.

Non-pyrogenic. Sterile if package is unopened and undamaged. Made in Ireland

Manufacturer:
Merit Medical Systems, Inc.
1600 West Merit Parkway, South Jordan, Utah 84095 U.S.A.
1-801-253-1600
U.S.A. Customer Service 1-800-356-3748

Authorized Representative:

Merit Medical Ireland Ltd, Parkmore Business Park West, Galway, Ireland

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