**StabiliT<sup>®</sup> VP**  
**Vertebroplasty System**  
**(For use with StabiliT<sup>®</sup> Bone Cement)**

# StabiliT® VP Vertebroplasty System (For use with StabiliT® Bone Cement)

## Instructions for Use (English) (EN)

### Important Information – Please Read Before Use


**CAUTION:** Federal law (USA) restricts this device to sale by or on the order of a physician.

### INDICATIONS

The StabiliT VP Vertebroplasty System is intended for percutaneous delivery of StabiliT Bone Cement in vertebroplasty procedures for the treatment of pathological fractures of the vertebrae. Painful vertebral compression fractures may result from osteoporosis, benign lesions (hemangioma), and malignant lesions (metastatic cancers, myeloma).

### DESCRIPTION

The StabiliT VP Vertebroplasty System is a system for the controlled delivery of StabiliT Bone Cement in the treatment of vertebral compression fractures. The StabiliT VP Vertebroplasty System consists of five (5) components (Figure 1).

**Figure 1. StabiliT VP Vertebroplasty System**

- StabiliT Delivery Syringe (#1) - Used for delivering bone cement into the vertebral body.
- StabiliT Delivery Instruments (#2, 3, 4) - Used to provide access to the vertebral body and facilitate delivery of PMMA bone cement in orthopedic procedures including vertebroplasty or kyphoplasty procedures.
  - Master Syringe (#2)
  - Locking Delivery Cannula (LDC) (#3)
  - StabiliT Introducer (#4)
    - Working Cannula
    - Stylet (diamond tip)
- StabiliT Bone Cement and Saturate Mixing System (#5)-The StabiliT Bone Cement and Saturate Mixing System is intended for:
  - The StabiliT Bone Cement is indicated for the treatment of pathological fractures of the vertebrae using a vertebroplasty or vertebral augmentation procedure. Painful vertebral compression fractures may result from osteoporosis, benign lesions (hemangioma), and malignant lesions (metastatic cancers, myeloma).
  - The Saturate Mixing System is intended for mixing of StabiliT Bone Cement.
    - Funnel (#5a)
    - Cement Syringe (#5b)
    - Filter (#5c)
    - Locking syringe (#5d)
    - Monomer (#5e)

### HOW SUPPLIED

- The StabiliT Delivery Syringe, StabiliT Delivery Instruments, and StabiliT Bone Cement and Saturate Mixing System are provided sterile. These devices are intended for single use only. DO NOT re-sterilize and/or reuse. DO NOT use if package is open or damaged and notify the manufacturer.

### CONTRAINDICATIONS

- The use of this product is contraindicated in patients with coagulation disorders, or with severe pulmonary insufficiency.
- The use of this product is contraindicated in patients with a compromise in the posterior column of the vertebral body or the walls of the pedicles.
- The use of PMMA bone cement is contraindicated in the presence of active or incompletely treated infection at the site where the bone cement is to be applied.
- This product should not be used in patients with sensitivity to any of the components of the PMMA bone cement.

### WARNINGS

- This device should only be used by qualified physicians with training in the clinical procedure in which it is being used.
- For safe and effective use of the StabiliT VP Vertebroplasty System, the physician should have specific training, experience, and thorough familiarity with the use and application of this product.
- Percutaneous vertebroplasty procedures should only be performed in medical settings in which emergency decompressive surgery is available.
- It is essential to maintain strict sterile technique during the vertebroplasty procedure and during all phases of handling this product.
- Precise Working Cannula placement is required for this procedure. Incorrect device placement could result in patient injury.
- The Working Cannula (part of the StabiliT Introducer) is not intended for delivering bone cement. Always use the LDC to deliver bone cement to the vertebral body.
- Insertion of the Working Cannula must be performed with the Stylet in place inside the Working Cannula.
- The Stylet or the Locking Delivery Cannula must be within the Working Cannula during manipulation or repositioning of the Working Cannula.
- Removal of the Working Cannula must be performed by rotation and axial motion. DO NOT bend the cannula sideways, patient injury may occur.
- Use appropriate imaging techniques to confirm correct Working Cannula placement, absence

of damage to surrounding structures, and appropriate location of delivered bone cement. Imaging, such as venography, can be used to assess the ability of the vertebra to contain the delivered bone cement.

- Thoroughly read the IFU that is packaged with the StabiliT Bone Cement and Saturate Mixing System before use.
- Dispose used product per Local, State and Federal Blood borne pathogen controls including Biohazard sharps container and disposal procedures.
- All devices are provided sterile. These devices are intended for single use only. DO NOT re-sterilize. DO NOT use if package is open or damaged. Reconditioning, refurbishing, repair, modification, or reesterilization of the device(s) to enable further use is expressly prohibited as may result in loss of function and/or infection.
- For devices penetrating bone, DO NOT use if dense bone is encountered. Device damage resulting in patient injury may occur. Breakage of the device may require intervention for retrieval.
- IFUs for each device (if packaged separately or along with this IFU) must be followed to perform a procedure using the StabiliT VP Vertebroplasty System.

### CAUTIONS

- Examine all packaging prior to opening. DO NOT use device if damaged or the sterile packaging is breached.
- Wear safety glasses or a face shield when delivering the bone cement.
- Ensure that all luer lock connectors are tightened securely. Improperly secured connections could result in disconnection during injection.
- Exercise caution in cases involving extensive vertebral destruction and significant vertebral collapse (i.e., the vertebral body is less than 1/3 of its original height). Such cases may lead to a technically difficult procedure.

### ADVERSE EVENTS

- Serious adverse events, some with fatal outcome, associated with the use of polymethylmethacrylate (PMMA) include:
  - Myocardial infarction
  - Cardiac arrest
  - Cerebrovascular accident
  - Pulmonary embolism
  - Anaphylaxis
  - Diffusion of the bone cement outside the vertebral body: in the peripheral veins (pulmonary embolism), in the epidural plexus (myelopathy, radiculopathy), in the intervertebral disc
- The most frequent adverse reactions reported with PMMA are:
  - Transitory fall in blood pressure
  - Thrombophlebitis
  - Hemorrhage and hematoma
  - Superficial or deep wound infection
  - Bursitis
  - Short-term cardiac irregularities
  - Heterotopic bone formation
- Other potential adverse events reported for PMMA include:
  - Pyrexia
  - Hematuria
  - Dysuria
  - Bladder fistula
  - Transitory worsening of pain due to heat released during polymerization
  - Nerve entrapment and dysphasia due to extrusion of the bone cement beyond its intended application
  - Adhesions and stricture of the ileum due to heat released during polymerization
- Potential adverse events associated with vertebroplasty include:
  - Pneumonia
  - Intercostal neuralgia
  - Collapse of a vertebra adjacent to the one injected, due to an osteoporotic disease
  - Pneumothorax
  - Extravasation of bone cement into soft tissue
  - Fracture of a pedicle
  - Rib fracture in patients with diffuse osteopenia, especially during thoracic vertebroplasty procedures, due to the significant downward force exerted during Working Cannula insertion
  - Compression of the spinal cord with paralysis or loss of feeling

### PREPARATION AND USE

- 1) Check packaging for damage prior to placing contents in sterile field.
- 2) Remove product from package using standard sterile technique.
- 3) Check all components for damage. If the pressure gauge needle on the StabiliT Delivery Syringe is resting outside the "0" box, do not use.
- 4) Mix bone cement per the StabiliT Bone Cement and Saturate Mixing System IFU.
- 5) Access to the vertebral body
  - a) Under fluoroscopic guidance place the Introducer into the vertebral body. Alternately, the LDC can be used to access the vertebral body directly instead of using the Introducer.
  - b) Once the Introducer is positioned in the vertebral body, remove the Stylet with a counterclockwise turn leaving the Working Cannula in place.
- 6) Prepare the StabiliT Delivery Syringe
  - a) Remove stopcock from StabiliT Delivery Syringe packaging
  - b) Connect the StabiliT Delivery Syringe extension tube to the female connector on the stopcock (opposite from the male connector).
  - c) Ensure the male port on the stopcock is set to the closed position.
  - d) Advance the plunger with enough force to completely remove any air present in the syringe.
  - e) Submerge the stopcock in sterile water (or saline).
  - f) Squeeze the trigger on the StabiliT Delivery Syringe and pull back the handle to fill the Syringe with fluid. Do so until the entire Syringe is filled.
  - g) Invert the StabiliT Delivery Syringe, squeeze the trigger, and advance the plunger to remove any remaining air in the syringe and extension tube.
  - h) Squeeze the trigger on the StabiliT Delivery Syringe and pull back fully to aspirate with sterile water (or saline).

**Caution:** Inspect StabiliT Delivery Syringe tubing and stopcock to ensure that there is no air in the system.
- 7) Assemble system components

- a) One minute prior to bone cement delivery (see Table 1) remove the Filter and Funnel Assembly. Clean the Cement Syringe of excess bone cement and completely thread the Cement Syringe into the Locking Delivery Cannula.
- b) Thread the Master Syringe onto the Cement Syringe.  
**Caution:** Ensure Master Syringe is FULLY threaded onto Cement Syringe before proceeding. Not doing so can cause user injury or device malfunction.
- c) Attach the StabilIT Delivery Syringe to the Master Syringe. Securely connect the luer of the Master Syringe to the stopcock at the end of the StabilIT Delivery Syringe extension tube.  
**Caution:** Do not begin bone cement delivery until saturation and preparation time is complete (see Table 1).
- 8) Cement delivery
- a) Confirm that the StabilIT Delivery Syringe trigger is released to ensure that the plunger is locked in position.
- b) Prime the LDC with bone cement by rotating the StabilIT Delivery Syringe handle in the CLOCKWISE direction. Once bone cement exits LDC tip, stop cement flow by squeezing the trigger on the StabilIT Delivery Syringe. Release trigger to lock the plunger in the withdrawn position. Wipe LDC tip clean.
- c) Under fluoroscopic guidance, stabilize the Working Cannula and insert the LDC until the rotating wheel contacts the Working Cannula luer. Rotate the LDC wheel to lock the LDC to the Working Cannula.
- d) When prepared to deliver bone cement squeeze the StabilIT Delivery Syringe trigger and push the handle forward until resistance is met and release the trigger. Under fluoroscopic guidance, deliver bone cement by rotating the handle in the CLOCKWISE direction.
- e) To stop bone cement delivery, squeeze the trigger on the StabilIT Delivery Syringe. Release trigger to lock the plunger in the withdrawn position. To re-engage, squeeze the trigger and push the handle forward until resistance is met, then release trigger. Continue delivering bone cement by rotating the handle in the CLOCKWISE direction.  
**Caution:** To protect the threads of the lock release handle, the gauge must indicate 25 ATM or lower before the quick release mechanism is used to pull back.  
**Caution:** The clutch will disengage (signaled by a clicking sound) if the operator exceeds the maximum pressure for the StabilIT Delivery Syringe. Once this has occurred, the clutch may disengage at lower pressures during subsequent attempts to increase pressure.  
**Caution:** Following completion of bone cement delivery, remove the LDC from the Working Cannula within 1 minute and immediately insert and lock the Stylet in the Working Cannula. If no additional bone cement delivery is required, remove Introducer (Cannula with Stylet).  
**Warning:** Removal of the Working Cannula should only be performed after insertion of the Stylet by rotation and axial motion. DO NOT bend Working Cannula sideways, patient injury may occur.

Table 1: Timing of Various Activities at Different Ambient Temperatures

Activity	Approximate Cumulative Time From Initiation of Saturation (minutes)	
	@ 18-19 °C (65-67 °F)	@ 20-23 °C (68-74 °F)
Bone Cement Saturation and Preparation (See Bone Cement IFU)	0-10 minutes	0-5 minutes
Bone Cement Delivery	10-45 minutes (35 minutes working time)	5-40 minutes (35 minutes working time)

Table 2: Effect of Ambient Temperature on Cement

Temperature °F (°C)	Minimum Setting Time (minutes)
81 (27)	35
73 (23)	51
66 (19)	82

#### STORAGE & HANDLING

Handle with care.  
 Store in original packaging in a clean, cool, and dry location.  
 Avoid exposure to temperature and humidity extremes.

#### SYMBOL GLOSSARY

	Caution		Use By
	Consult Instructions For Use		Authorized Representative in the European Community
	Sterilized using Irradiation		Keep away from sunlight
	Sterilized using Ethylene Oxide		Keep away from moisture
	Sterilized using Aseptic Processing Techniques		Flammable
	Lot Number		Single Use Device, DO NOT REUSE
	Catalog Number		Store below 25°C
	No Latex		Manufacturer
	Do not use if package is opened or damaged		Diamond tip
	Short		Bevel Tip
	Long		Cannula Gauge
	Device Length in Centimeters		Introducer Gauge/Locking Delivery Cannula Gauge
	Do not Resterilize		Quantity

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Merit Medical Systems, Inc.
South Jordan, Utah 84095 U.S.A.
U.S.A. Customer Service 1-800-356-3748
Merit Medical Ireland Ltd.