STAR™ Tumor Ablation System

STAR™ Tumor Ablation System

Instructions for Use (English) (EN)

Instructions for Use

Do not Reuse

Important Information – Please Read Before Use

CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician.

INDICATIONS

The STAR Tumor Ablation System is indicated for palliative treatment in spinal procedures by ablation of metastatic malignant lesions in a vertebral body.

CONTRAINDICATIONS


- SpineSTAR Ablation Instrument is contraindicated in patients with heart pacemakers or other electronic device implants.
- SpineSTAR Ablation Instrument is contraindicated in vertebral body levels C1-C7.

DESCRIPTION

The SpineSTAR Ablation Instrument is a sterile single use device for the ablation of tumor within the vertebral body. It is to be used with the StabiliT® Introducer (Stylet and Working Cannula), the AE Cable, Hand Switch Cable and the MetaSTAR® RF Generator.

The STAR Tumor Ablation System consists of seven (7) components.

StabiliT® Introducer (Stylet and Working Cannula)


StabiliT Introducer (with bevel tip stylet)

StabiliT Introducer (with diamond tip stylet)

bevel tip stylet

diamond tip stylet

Osteotomes

VertecoR StraightLine Cement Staging Osteotome

PowerCURVE™ Navigating Osteotome

SpineSTAR® Ablation Instrument

SpineSTAR Ablation Instrument

Cables

Activation Element (AE) Cable

Hand Switch Cable

StabiliT Introducer (Stylet and Working Cannula) (Bevel or Diamond Stylet Inserted the Working Cannula)

- The StabiliT Introducer with bevel or diamond tip stylet is used for percutaneous bone access.
 - The device is packaged with a bevel tip stylet and cannula.
 - The Introducer cannula is 3.6 mm in outer diameter.
 - See the Table 1 for the device length.


Table 1. StabiliT Introducer with bevel or diamond tip stylet

Device	Stylet		Size	Working Length in Patient	Gauge	
	Stylet Tip	Extends Beyond Working Cannula				
StabiliT Introducer	Diamond	0.7 cm	Short	10.0 cm	10G	
StabiliT Introducer	Bevel		Long	12.0 cm		
StabiliT Introducer	Diamond					
StabiliT Introducer	Bevel					

VertecoR® StraightLine Cement Staging Osteotome without Stylet (Inserted in the Working Cannula)

- The VertecoR StraightLine Cement Staging Osteotome (SLO) is for scraping or coring of bone in the spine.
 - The device is packaged with a blunt tip Stylet.
 - It is to be used with the StabiliT Introducer.
 - The StraightLine cannula is 3.0mm in outer diameter.
 - See the Table 2 for working length.


Table 2. VertecoR StraightLine Cement Staging Osteotome (SLO)

Device	Size	Extends beyond Working cannula	Working Length in Patient	Gauge
VertecoR StraightLine Cement Staging Osteotome	Short	1.5 cm	11.0 cm	11G
VertecoR StraightLine Cement Staging Osteotome	Long		13.0 cm	

PowerCURVE™ Navigating Osteotome (Inserted in the Working Cannula)

- The PowerCURVE Navigating Osteotome is for scraping or coring of bone in the spine.
 - It is to be used with the StabiliT Introducer.
 - The Shaft is 3.0 mm in outer diameter.
 - See the Table 3 for working length.


Table 3. PowerCURVE Navigating Osteotome

Device	Size	Extends beyond Working cannula	Working Length in Patient	Gauge
PowerCURVE™ Navigating Osteotome	Short	3.0 cm	12.5 cm	11G
PowerCURVE™ Navigating Osteotome	Long		14.5 cm	


Figure 1: The Deployment Handle (A) is the rotating portion that articulates the distal tip of the device. The Tip Direction Indicator (B) indicates the direction in which the tip will articulate.

SpineSTAR Ablation Instrument (Inserted in the Working Cannula)

- The SpineSTAR Ablation Instrument is for ablation of tumor within the vertebral body.
 - It is to be used with the StabiliT Introducer.
 - The SpineSTAR cannula is 3.0 mm in outer diameter.
 - See the Table 4 for working length.

13.0cm working length in patient (short)
15.0cm working length in patient (long)

Table 4. SpineSTAR Ablation Instrument

Device	Thermocouple Configuration (Distal / Proximal)	Size	Extends beyond Working cannula	Working Length in Patient	Gauge
SpineSTAR Ablation Instrument	10 / 15	Short	3.5 cm	13.0 cm	11G
SpineSTAR Ablation Instrument	5 / 10				
SpineSTAR Ablation Instrument	10 / 15	Long		15.0 cm	
SpineSTAR Ablation Instrument	5 / 10				


Figure 2

Figure 2: SpineSTAR (SS) Ablation Instrument (A) has a handle that contains an Articulator (grey color) that curls the Articulating Segment in the direction of the Indicator (white arrow) and white Electrode Knob, which extends the Electrode from the distal end of the Articulating Segment (C). The total maximum distance the Electrode extends beyond the end of the Articulating Segment is 10mm. There are two thermocouples (TC) on the articulating segment of the shaft (B): For 10/15 TC configuration, the Distal TC is 10mm from the insulation tip while Proximal TC is 15mm from the insulation tip; for the 5/10 TC configuration, the Distal TC is 5mm from the insulation tip while Proximal TC is 10mm from the insulation tip. Three black marks on the shaft represent the location of bevel tip and TCs in relation to the distal end of the Working Cannula. The Insulation Tip (C) is radiolucent and represents the center of the ablation zone.

Cables

- The AE Cable connects to the MetaSTAR Generator and SpineSTAR Ablation Instrument.
- The Hand Switch Cable connects to the MetaSTAR Generator. It activates and stops the RF energy during treatment.

HOW SUPPLIED

- The STAR™ Tumor Ablation System and all components are supplied sterile in a peel-open package. The kit is sterilized by gamma irradiation. This kit is intended for single use only. Do not resterilize. Do not use if the package is open or damaged. In the event of damage to the sterile packaging, notify the manufacturer.

WARNINGS:

- This device should only be used by qualified physicians with training in the clinical procedure in which it is being used.
- For safe use of the device, the physician should have specific training, experience, and thorough familiarity with the use and application of this product.
- Do NOT ablate in painful osteoporotic vertebra without tumor.
- Do NOT use the device in fractures due to osteoblastic metastases to the spine (such as purely blastic prostate metastatic lesions).
- Do NOT use this product in dense (sclerotic) bone including traumatic high energy fractures or purely blastic lesions, device damage resulting in patient injury may occur.
- Do NOT use this device in patients without metastatic malignant lesions in a vertebral body.
- Do NOT use this device in patients with multiple myeloma, solitary plasmacytoma, or primary malignant lesions in the index vertebra.
- Do NOT use a device in more than one vertebral body.
- Do NOT touch the Electrode at the tip of the device while power is being applied.
- Do NOT withdraw the device while power is still being applied.
- Do NOT touch the Electrode to metal objects while power is being applied, damage to the device may occur.
- Use standard electrosurgical precautions when using the device in the vicinity of nerves and nerve roots.
- It is essential to maintain strict sterile technique during all phases of handling and use of this product.
- Dispose used product per Local, State and Federal blood borne pathogen controls including Biohazard sharps container and disposal procedures.
- The product is sterilized by gamma radiation. Do NOT use if package is opened or damaged.
- Breakage of the device may require intervention or retrieval.
- All device manipulations (bending and straightening of the Articulating Segment or advancing and retracting the device) should be done slowly and under imaging guidance while carefully monitoring the position of the device and only with the Electrode fully retracted.
- Extension and retraction of the Electrode should be done slowly and under imaging guidance to monitor the position of the tip of the device.
- Ablation procedure must be performed under image guidance. Do NOT perform ablation without first imaging area as it may result in severe injury to patient.

- The device must be positioned such that when ablation is performed per the thermal distribution graphs (see MetaSTAR RF Generator Operator's Manual), nerves and nerve roots are beyond the ablation zone. Failure to utilize the thermal distribution graphs in pre-operative planning may result in severe injury to the patient.
- Do NOT use flammable anesthetics.
- Observe fire precautions at all times. Sparking and heating associated with electro surgery may be an ignition source.
- Do NOT rotate the blue Handle when the SpineSTAR Ablation Instrument is extended beyond the Working Cannula, breakage of the device may occur.
- It is essential to maintain strict sterile technique during all phases of handling and use of this product.
- The Introducer Stylet or the Locking Delivery Cannula should be inserted and engaged with the Working Cannula during Introducer removal or manipulation.
- Do NOT use the PowerCURVE™ Navigating Osteotome or the VertecoR StraightLine Cement Staging Osteotome to scrape or core in more than one vertebra.

PRECAUTIONS:

- It is important to read the Instructions For Use and these precautions prior to device operation.
- Use the device prior to Use By Date noted on the package.
- Do NOT use damaged product. Before use, inspect the device and packaging to verify that no damage has occurred.
- Physicians using the SpineSTAR Ablation Instrument should be familiar with the physiology and pathology of the selected anatomy, and be trained in the performance of the chosen surgical technique.
- The device should be manipulated only while under fluoroscopic (and/or CT) observation with radiographic equipment that provides high quality images.
- Do NOT allow patient contact with grounded metal objects.
- As with other electrosurgical units, electrodes and cables can provide paths for high frequency current. Position cables and electrode to avoid contact with the patient or other leads.
- The location of the Working Cannula in the vertebra should be confirmed using image guidance before advancement as well as during removal of the SpineSTAR Ablation Instrument.
- Do NOT re-sterilize and/or reuse. The device is for single use only. Reconditioning refurbishing, repair, modification, or resterilization of the device to enable further use is expressly prohibited as it will result in device malfunction and/or infection.
- Do NOT clean device with Isopropyl Alcohol or other solvents.
- Do NOT use device in more than one vertebral body.
- Remove the Stylet prior to inserting the StraightLine Osteotome into the Working Cannula. Not doing so will result in inadequate coring of bone.
- The location of the StabiliT Introducer's Working Cannula in the vertebra should be monitored under imaging before and during advancement of the StraightLine and PowerCURVE Navigating Osteotome.

ADVERSE EVENTS:

- As a consequence of electrosurgery, damage to surrounding tissue through iatrogenic injury could occur
- Nerve injury including thermal injury, puncture of the cord or nerve roots potentially resulting in radiculopathy, paresis or paralysis
- Pulmonary embolism
- Hemorrhage
- Hemothorax or pneumothorax
- Hematoma
- Infection including deep or superficial wound infection
- Pain
- Unintended puncture wounds including vascular puncture and dural tear

INTRODUCERS AND OSTEOTOME DIRECTIONS FOR USE:

- Check packaging for damage prior to placing contents in sterile field.
- Remove product from package using standard sterile technique.
- Access the vertebral body using the StabiliT Introducer.
- Under image guidance direct the StabiliT Introducer down the pedicle of the vertebral body. If using fluoroscopy check A/P and lateral images to confirm proper device placement.
- Once the StabiliT Introducer is positioned in the vertebral body, remove the Stylet with a counterclockwise turn, leaving the Working Cannula in place.
- Create an access channel using the SLO and/or the PowerCURVE.
Warning: Use image guidance and follow the IFU to avoid patient injury.
- Remove the Stylet from the SLO and set aside prior to inserting the SLO in the Working Cannula.
- Remove the Stylet from the Working Cannula.
- Using imaging guidance, advance the SLO through the Working Cannula until the tip of the SLO is at the desired location. The luer of the Working Cannula limits the SLO shaft to extend 15 mm beyond the distal end of the Working Cannula.
Caution: The location of the Working Cannula in the vertebra should be monitored before and during advancement of the SLO.
- Note that retraction of the Working Cannula may be required to completely seat the SLO against the Working Cannula.
- When scraping or coring is complete, using imaging guidance, remove the SLO from the Working Cannula.
- Use the PowerCURVE to create access channel(s) in the bone.
- The Articulating Tip is the distal portion of the PowerCURVE. The Tip Direction Indicator points in the direction which the articulating tip bends.
- The Deployment Handle is the rotating portion of the PowerCURVE Handle.
- Turning the Deployment Handle one (1) full 360° turn clockwise will cause the Articulating Tip to fully deploy (bend).
- Turning the Deployment Handle counter-clockwise will cause the Articulating Tip to straighten.
- Do not deploy the Deployment Handle greater than one half (1/2) turn counter clockwise from the neutral position.
- Ensure the Articulating Tip is fully extended in the straight position prior to insertion into the Working Cannula.
Caution: The location of the Working Cannula in the vertebra should be monitored before, and during advancement of the PowerCURVE through the Working Cannula. Insert the PowerCURVE into the Working Cannula until the first laser mark on the shaft is even with the proximal end of the luer on the Working Cannula. Confirm by image guidance that the distal end of the PowerCURVE is at the distal end of the Working Cannula before proceeding. When fully inserted into Working Cannula, the shaft extends approximately 3.0 cm beyond the distal end of the Working Cannula.

19. As the Articulating Tip exits the Working Cannula turn the Deployment Handle to bend the Articulating Tip in the direction of the Tip Direction Indicator on the proximal shaft of the PowerCURVE
20. The PowerCURVE should be advanced to the desired position using image guidance.
21. Care should be taken at all times to NEVER strike the arms of the Deployment Handle, especially when rotated from its starting position.
22. Using image guidance (while stabilizing the Working Cannula), the PowerCURVE™ can be carefully withdrawn and advanced multiple times to scrape or core bone until the desired cavity (size and location) is created.
23. When the Articulating Tip is substantially deployed, the PowerCURVE should not be rotated.
24. The rotation limiting mechanism will slip if the Handle is rotated while the Articulating Tip is substantially articulated.
25. When cavity creation is complete, using imaging guidance straighten the Articulating Tip and return the Deployment Handle to its starting position.
- Caution:** Straightening the Articulating Tip should be done slowly and under imaging guidance while carefully monitoring the position of the tip of the device. Straighten the device by turning the Deployment Handle counter clockwise.
26. Under imaging guidance, remove the PowerCURVE from the Working Cannula.
27. The location of the Working Cannula in the vertebra should be monitored and adjusted if necessary after removal of the PowerCURVE.

SPINESTAR ABLATION INSTRUMENT (SS) DIRECTIONS FOR USE

Warning: All device manipulations should be done/confirmed under imaging guidance (CT or fluoroscopy) that provides high quality images.

Warning: Do not insert, withdraw or articulate or straighten the SpineSTAR (SS) while the Electrode is extended.

1. Connect the Hand Switch Cable to the MetaSTAR RF Generator.
2. Connect the AE Cable to the MetaSTAR RF Generator and the SpineSTAR.
3. Confirm the MetaSTAR RF Generator recognizes the SS – if it does not, the graphical display will illuminate a flashing warning “connect device”.
4. To bend the Articulating Segment - turn the Articulator (grey) clockwise. It will only bend in one direction, the direction of the large white arrow on the blue handle
- Caution:** DO NOT turn the Articulator greater than one half (1/2) turn clockwise.
5. To straighten the Articulating Segment - turn the Articulator counterclockwise.
6. Insert the Introducer to desired location and remove the stylet from the WC.
7. Ensure the Electrode is fully retracted and the Articulating Segment is in straight position prior to inserting the SS into the Working Cannula (WC).
8. In combination with image guidance, use laser marks on SS shaft to verify position of SS and the TCs in relation to end of the WC:
- a. Distal mark on SS shaft is level with WC luer = SS bevel tip is at distal end of WC.
 - b. Middle mark on SS shaft is level with WC luer = SS Distal TC is 2mm beyond distal end of WC.
 - c. Proximal mark on SS shaft is level with WC luer = SS Proximal TC is 2mm beyond distal end of WC.

Note: Distance SS extends beyond end of WC when fully inserted is:

- 30mm when Electrode is fully retracted.
- 35mm when Electrode is fully extended.

9. After inserting the SS bevel tip to the distal end of the WC CONFIRM the MetaSTAR RF Generator graphic display registers $\geq 35^{\circ}\text{C}$ for both proximal and distal TCs before proceeding.
10. As the Articulating Segment of the SS is advanced beyond the WC into the vertebra, rotate the Articulator to bend the Articulating Segment, as necessary, to navigate to the desired location. The degree to which the Articulating Segment will actually bend will depend on bone density.

Caution: When the Articulating Segment is bent, the device handle should not be rotated.

11. Once positioned as desired, extend the Electrode to the desired length by rotating the white Electrode Knob clockwise.

Caution: Do NOT rotate the SS Handle when the Electrode is extended.

Caution: Monitor the location of WC before advancing the SS through it.

12. Prior to initiating ablation confirm (with CT or fluoroscopic imaging) location of the distal end of WC and location of Electrode.

13. Confirm AE Cable is inserted into Cable Connector of SS.
14. Set power level on the MetaSTAR RF Generator to the desired level to create thermal distribution for the major and minor diameter (see in Figure 3). Refer to MetaSTAR RF Generator IFU for Thermal Distribution Diameter Graph showing major and minor diameter at various power levels.

15. Press the Blue Button on the Hand Switch to start ablation. A light on the MetaSTAR RF Generator will illuminate and an audible tone will emit from the MetaSTAR RF Generator.
16. During ablation the MetaSTAR RF Generator graphic display must be observed to monitor temperatures and impedance.

- a) TCs may be used to monitor temperature propagation along the length of the SS.
 - 50°C on the SS at the 5 mm TC represents an ablation zone approximately 10 mm long
 - 50°C on the SS at the 10 mm TC represents an ablation zone approximately 20 mm long
 - 50°C on the SS at the 15 mm TC represents an ablation zone approximately 30 mm long
- b) MetaSTAR RF Generator software is designed to stop RF delivery when any one of following three conditions occur:
 - Proximal TC reaches 50°C .
 - Impedance is too high.
 - Cycle time counts down to zero.

Warning: If proximal TC temperature suddenly drops below 30°C ; STOP THE PROCEDURE and replace the SS.

17. If larger ablation zone is desired, retract the Electrode and reposition the Articulating Segment in desired location. Extend the Electrode and reinitiate ablation.
18. If impedance limit is reached, the MetaSTAR RF Generator will stop RF delivery. Retracting and extending the Electrode or slight manipulation of the Articulating Segment (with the Electrode retracted) may be required to continue ablation. In some cases removal of the SS from the WC and cleaning the Electrode of char may be required.
19. After the ablation is completed, remove the SS per the following steps:
- Fully retract the Electrode.
 - a. If Articulating Segment has been deployed in the anterior third of the vertebra, pull the SS back to the middle third of the vertebra before straightening.
 - b. Confirm the Electrode is still fully retracted before withdrawing the distal end of the SS into WC.

Caution: Stabilize the Working Cannula while removing the SS.

Warning: Do not withdraw the SS while power is being applied.


Figure 3: Ablation Zone and its dimensions

Warning: If the patient complains of SUDDEN ONSET OF PAIN during ablation, STOP THE PROCEDURE and:

- Closely examine the position of the device using CT or fluoroscopic imaging.
- Confirm the path of and placement of device.
- Confirm distal tip of Electrode is in the desired location within vertebra.
- DO NOT continue procedure until proper placement of device has been verified.

Warning: Excessive ablation may create heat within the spine and adjacent neural tissue.

Warning: See MetaSTAR RF Generator Operator’s Manual for selection of power level and ablation time.

SYMBOLS GLOSSARY:

	Caution		Use By Date
	Consult Instructions For Use		
	Sterilized using irradiation		Keep away from sunlight
	Lot Number		Keep dry
	Catalog Number		Manufacturer
	Single Use Device, DO NOT REUSE		Not made with natural rubber latex
	10/11 gauge cannula outer diameter		Do not use if package is damaged
	Longer length instrument		Thermocouple location 5 and 10mm from the Insulation Tip
	Shorter length instrument		Thermocouple location 10 and 15mm from Insulation Tip
	Quantity		Do Not Sterilize

STORAGE & HANDLING

- Handle with care. Store in original packaging in a clean, cool and dry location. Avoid exposure to temperature and humidity extremes.

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Merit Medical Systems, Inc.
U.S.A.
U.S.A. Customer Service