| | |
|--|---|
| en Instructions for use | bg Инструкции за употреба |
| da Brugsanvisning | hr Upute za upotrebu |
| nl Instructie voor gebruik | cs Návod k použití |
| fi Käyttöohjeet | hu Használati útmutató |
| fr Mode d'emploi | no Bruksanvisning |
| de Gebrauchsanleitung | pl Instrukcja użytkowania |
| el Οδηγίες χρήσης | ro Instrucțiuni pentru utilizare |
| it Istruzioni per l'uso | sk Návod na použitie |
| pt Instruções de Uso | tr Kullanma Talimatı |
| ru Инструкции по применению | et Kasutamishuend |
| es Instrucciones de uso | lt Naudojimo instrukcija |
| sv Bruksanvisning till användning | zh 使用說明 |

Exacta™ Percutaneous Sheath Introducer

The following are general instructions for use of the Exacta™ Percutaneous Sheath Introducer kit. Clinicians are strongly recommended to follow this procedure, or an alternative with which they are familiar, or which conforms to their hospital standard practice.

Description

The Merit Exacta Percutaneous Sheath Introducer is a semi-rigid introducer through which a balloon-tipped catheter and central venous catheter can be introduced percutaneously. The introducer is supplied in a kit that includes a scalpel, guidewire, over-needle introducer, and needle introducer, syringe, contamination shield, and Dilator/IObuturator.

Intended Use

Merit Exacta Percutaneous Sheath Introducers are surgically invasive devices that are intended for short term use (no longer than 30 days) to allow introduction of balloon-tipped catheters into the human central circulatory system. This venous system access allows direct measurement of pressures within the chambers of heart, blood sampling, thermolubrication cardiac measurements, infusion of drugs & fluids and transvenous pacing.

All cautions should be read and understood before attempting to use this kit.

CAUTION Use only as directed by a physician. A full list of components is on the package lid.

General Warnings

- Complications associated with percutaneous sheath introduction include vessel wall perforation, pleural and mediastinal injury, air embolism, catheter embolism, thoracic duct laceration, bacteraemia, septicaemia, thrombosis, and inadvertent arterial puncture.
- The patient should be placed in a slight Trendelenburg position for subclavian or jugular introduction.
- Luer connections: as standard practice, the security of luer connections must be checked routinely.
- When removing dressings at or close to catheter sites, care must be taken to avoid severing the catheter.
- Acetone must not come in contact with the catheter as the material may weaken and this may result in leakage or aspiration.
- Exposure of product componentry to topicals containing alcohol is not recommended.
- Do not attempt to re-insert a partially or completely withdrawn introducer cannula.
- Percutaneous puncture of a central vein may be contraindicated in patients with pulmonary hypertension.
- Use of a syringe smaller than 5 mL to irrigate or de-clot an occluded catheter may cause catheter rupture.
- Syringes are supplied for blood aspiration only. **Warnings - Seldinger Technique**

- Do not reinsert the needle into catheter.
- Maintain a firm grip on the guidewire at all times. Do not force the guidewire. If resistance is met, carefully withdraw the guidewire and re-attempt insertion.
- On withdrawing the guidewire, when using the 18 gauge thin-wall needle, caution should be taken to prevent damage or possible severing of the wire on the needle bevel.
- Potential for guidewire breakage. Although the incidence of guidewire breakage is extremely uncommon, physicians must be aware of the potential of guidewire breakage if undue

force is applied to the wire. If resistance is met when attempting to remove the guidewire after Sheath placement, the wire may be kinked within the area of the sheath tip and the vessel. Undue force may cause the wire to break. If resistance is encountered, withdraw the sheath relative to the guidewire (2-3 cm) and attempt to remove the wire. If resistance is still apparent remove the wire and sheath simultaneously.

5. Physicians should be aware that the guidewire can pick up material from the vein. This may prevent the guidewire from being withdrawn through the needle.

Warnings - Sheath Insertion

- Tuohy Borst adaptor: when using the Tuohy Borst adaptor with the Exacta™ Percutaneous Sheath Introducer kit, caution should be taken not to insert introducer needles into the haemostasis valve.
- The Contamination Shield is not intended for use as a sterile barrier.
- Do not insert introducer needles or hypodermic needles through the contamination shield.
- Care should be taken when handling the assembly and when disconnecting the sideport from an infusion tube, to prevent air passing into the patient.
- If use of a Merit Percutaneous Sheath Introducer is delayed or interrupted, it is recommended that the relevant sized Merit obturator is used with the sideport assembly and sheath.

5 Fr Obturator for 6 Fr & 7 Fr Sheath Introducers.

7 Fr Obturator for 8 Fr, 8.5 Fr & 9 Fr Sheath Introducers.

This will maintain sheath patency, ensure against leakage and prevent contamination.


- Ensure that the valve mechanism and contamination shield are not damaged by sharp objects.
- Hospital Protocols should adhere to the following instructions for use of stopcock. Excessive force when connecting, particularly in the presence of lubricious infusates such as fat emulsions, can cause over stressing to stopcock ports. Contact with some organic solvents and infusion solutions can also relieve internal stresses causing the material to crack. To prevent potential cracks causing air embolism or other risks to the patient follow four main principles:

- Avoid over tightening connections.
- Check connections regularly.
- Change stopcocks every 72 hours.
- When lubricious infusates are used change stopcocks every 24 hours.

Venaguide™

- Release guidewire by removing the guidewire cap.
- Straighten guidewire 'J' by retracting into introducer system with thumb.
- Insert into hub of introducer needle and advance guidewire into vein. Advance to required depth.

WARNING - to prevent over insertion stop when last depth mark reaches the venepuncture site.

Patient Preparation

- Disinfect, and drape as required, the proposed puncture site.
- Perform a skin wheal with local anaesthetic at the proposed site.
- Insert the tip of a CritiCath™ or other pulmonary artery catheter into the Tuohy-Borst adaptor end of the contamination shield. Advance the catheter through tubing and hub at distal end. Slide contamination shield past the 50 cm mark on the catheter. Extend full length of the sheath and tighten Tuohy-Borst adaptor at proximal end of sheath.

Sheath Insertion

- Using the introducer needle (or over needle cannula) puncture the vein, aspirate and advance the needle (or cannula) into the vein. Failure to aspirate venous blood through the needle may indicate incorrect venous puncture.
 - Insert the desired end of the guidewire through the Introducer. (The 'J' end can be prepared by sliding the plastic insertion sleeve over the tip to straighten it.)
 - While holding the guidewire in place, remove the Introducer.
 - If required, enlarge the insertion site with the scalpel.
- CARE SHOULD BE TAKEN NOT TO CUT THE GUIDEWIRE.**
- Thread the tip of the dilator over the guidewire until the guidewire exits the hub of sheath/dilator assembly. While firmly holding the guidewire in place, grasp dilator near the skin and advance the sheath/dilator assembly into vessel with a slight twisting motion.
 - Holding the sheath in place, withdraw the guidewire and dilator sufficiently to allow venous blood flow to be aspirated into the side port. Connect to an appropriate intravenous infusion line and flush.
 - Always remove the dilator over the guidewire after dilation of the skin.
 - Holding the sheath in place, remove the guidewire and dilator.
 - Test the catheter balloon according to manufacturer's recommended procedure and insert tip of catheter into the Tuohy-Borst adaptor on sheath/sideport assembly.

13. Advance the catheter to the desired position. Secure catheter in place by tightening Tuohy-Borst adaptor on the sideport assembly. Use of the Tuohy-Borst adaptor may help prevent catheter migration. **CARE SHOULD BE TAKEN TO AVOID OVER-TIGHTENING OF THE TUOHY-BORST ADAPTOR**, as this may occlude the catheter lumen.

14. Suture the sheath in place utilizing suture loop and recessed ring if desired. If the sheath assembly is not to be used immediately, an obturator may be used to help preserve patency.

Physicians should be aware that central venous catheters are intended for use for up to 30 days.

Storage Conditions

Store in cool dry place away from direct sunlight.

STERILE and non-pyrogenic in unopened, undamaged package. For single use only. Check integrity of the individual package before use. Dispose of product after use. Do not resterilize.

Re-use may lead to infection or other illness/injury.

This product contains a phthalate which may cause a risk to {children (male), pregnant or nursing women} in conjunction with prolonged duration of use.

For reordering information or assistance please contact local representative.

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This symbol is applicable for US market only

Rx ONLY

CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a physician.

Non-Pyrogenic
 ikke-pyrogen
 NIET-PYROGEEN
 EI-PYROGEENINEN
 NON PYROGÈNE
 Nicht pyrogen
 ΜΗ ΠΥΡΕΤΟΓΟΝΟ
 APIROGENO
 apirogénico
 АПИРОГЕННЫЙ
 NO PIROGÉNICO
 ICKE-PYROGENT
 неапириногенно
 NEZAPALJIVO
 APYROGENNÍ
 Nem pirogén
 ikke-pyrogen
 NIE PYROGENNY
 APIROGEN
 PIROJENİK OLMAYAN
 Nepyrogenický
 Mitte-pürogeenne
 nepirogeninis or nepirogeninė
 無致熱原

Do Not Resterilize
 Må ikke resteriliseres
 Niet opnieuw steriliseren
 Tuotetta ei saa steriloida uudelleen
 Ne pas restériliser
 Nicht erneut sterilisieren
 Μη επανααστεριρώνετε
 Non risterilizzare
 Não reesterilize
 не стерилизовать повторно
 No reesterilizar
 Får inte steriliseras på nytt
 Да не се стерилизира повторно
 Nemojte ponovno sterilizirati
 Neresterilizujte
 Ne sterilizálja újra
 Skal ikke steriliseres mer enn én gang
 Nie sterylizować ponownie
 A nu se reesteriliza
 Zákaz opakovanej sterilizácie
 Yeniden sterilize etmeyin
 請勿重新消毒
 재소독하지 마십시오.
 再滅菌禁止
 Ärge steriliseerige uuesti
 Pakartotina nesterilizuokite
 Ne sterilizirajte ga
 Nav atkārtoti sterilizējams

Statement is only applicable for US market only
 This product can expose you to chemicals including Di(2-ethylhexyl) phthalate (DEHP), which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov

Merit Medical Singapore
 Pte. Ltd., 198 Yishun Avenue 7,
 Singapore 768926

Merit Medical Ireland Ltd,
 Parkmore Business Park West,
 Galway, Ireland
 European Customer Service
 +31 43 358 82 22

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