## ENGLISH

**BONE MARROW ASPIRATION NEEDLE**

(MALLARME)

**INSTRUCTIONS FOR USE**

**INDICATIONS FOR USE**  
Disposable biopsy needle used to aspirate sternal or iliac bone marrow.

**REUSE PRECAUTION STATEMENT**

For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

**STORAGE**

Keep the devices in their packaging, flat in the original box. Store in a cool dry place (below 26 °C) away from humidity and direct heat. Do not use after expiry date.

**PRECAUTIONS**

The operator should not use the device before being trained and before reading the instructions for use. It is a disposable device. After use, each component of the device may be a potential biohazard. Handle and dispose such devices in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

**INSTRUCTIONS**

- Check expiration date and integrity of the packaging.
- Position the patient according to the technique used (sternal aspiration or iliac aspiration technique).
- Using a sterile technique, prepare the puncture site.
- Infiltrate the skin, especially the periosteum, with local anesthetic.
- A scalpel blade may be used to puncture the skin prior to inserting the needle.
- Introduce the needle through the incision and, using firm pressure, advance the needle. Decreased resistance generally indicates entrance into the marrow cavity.
- Pull the stylet out of the cannula.
- Attach a syringe to the cannula hub.
- Apply negative pressure by quickly drawing the syringe plunger. Disconnect the syringe and remove the aspirated specimen.

**Note:** It is advisable to obtain a bone marrow aspirate in conjunction with a core biopsy specimen, since aspiration alone does not always permit an adequate evaluation of the functional state of the marrow.

|  |  |  |                                  |
|--|--|--|----------------------------------|
|  | Single Use   |  | Do not use if package is damaged |
|  | Temperature limitation   |  | Do not resterilize               |
|  | Caution: consult accompanying documents  |  | Sterilized using Ethylene Oxide  |
|  | For electronic copy scan QR Code, or go to <a href="http://www.merit.com/ifu">www.merit.com/ifu</a> and enter IFU ID Number. For printed copy, call U.S.A or E.U. Customer Service |  |                                  |

|   |   |  |
|---|---|--|
|   |   | <a href="http://www.merit.com">www.merit.com</a> |
| Manufacturer:<br>Merit Medical Systems, Inc.<br>1600 West Merit Parkway, South<br>Jordan, Utah 84095 U.S.A.<br>1-801-253-1600<br>U.S.A Customer Service<br>1-800-356-3748 | EC REP<br>Authorized Representative:<br>Merit Medical Ireland Ltd., Parkmore<br>Business Park West, Galway, Ireland<br>EC Customer Service<br>+31 43 358822 |  |

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