**INSTRUCTIONS FOR USE**

IMPORTANT INFORMATION - PLEASE READ BEFORE USE


CAUTION Federal (USA) Law restricts this device to sale by or on the order of a physician.

INDICATION The StabiliT® Introducer is indicated for percutaneous access to bone.

DESCRIPTION The StabiliT® Introducer is indicated for percutaneous access to bone and consists of a cannula with either a trocar or bevel tip style. The introducer cannula is 3.6mm in outer diameter. See the chart below for the device length.

Device Working Length

StabiliT Introducer (Short)	10 cm
StabiliT Introducer (Long)	12 cm

**WARNINGS**

- The device should only be used by qualified physicians with training in the clinical procedure in which it is being used.
- Avoid use of fluoroscopic guidance to avoid patient injury.
- For safe use of the StabiliT® Introducer, the physician should have specific training, experience, and thorough familiarity with the use and application of this product.
- It is essential to maintain strict sterile technique during all phases of handling and use of this product. The product is sterilized by Ethylene Oxide. Do not use if package is opened or damaged.
- Dispose of used product per Local, State and Federal Blood borne pathogen controls including Biohazard sharps disposal guidelines.
- Do NOT use the product in dead bone; device damage resulting in patient injury may occur. Breakage of the device may require intervention or retrieval.
- The Introducer Stylet should be inserted and engaged with the Working Cannula during Introducer removal or manipulation.
- For single patient use only. Do not reuse, reprocess or sterilize. Reuse, reprocessing or sterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or sterilization also creates a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

PRECAUTIONS

- It is important to read the Instructions for Use and these precautions prior to device operation.
- Use the StabiliT® Introducer prior to Use Date noted on the package to verify no damage has occurred.
- Physicians using the StabiliT® Introducer should be familiar with the physiology and pathology of the selected anatomy, and the performance of the chosen surgical technique.

The StabiliT® Introducer should be manipulated only while under fluoroscopic observation with radiographic equipment that provides high quality images.

**ADVERSE EVENTS**

Adverse events potentially associated with use of the device include:

- Nerve injury including puncture of the cord or nerve roots potentially resulting in radiculopathy, paresis or paralysis
- Emboli
- Hemorrhage
- Hemothorax or pneumothorax
- Hematoma
- Infection including deep or superficial wound infection
- Pain
- Unintended puncture wounds including vascular puncture and dural tear
- Fracture

DIRECTIVES

- Make a skin incision over the selected area for access. Using fluoroscopic imaging guidance, manually advance the StabiliT® Introducer to the selected bone surface checking AP/lateral images to confirm proper placement. Do not advance the cannula without the stylet fully inserted into the cannula.

2. While holding the cannula in place, turn the stylet counter-clockwise to release and remove it from the cannula. The cannula is now ready to accept other instrumentation.

NOTE: Markings on the cannula may be used as reference markers only, they are not intended to replace the use of fluoroscopic observation.

**SYMBOL GLOSSARY**

SYMBOL	DESIGNATION	SYMBOL	DESIGNATION
	Caution: Consult accompanying documents. Read instructions prior to use.		Catalog Number
	Use by: YYYY-MM-DD		Single Use Device, DO NOT REUSE
	Do not use if package is damaged		
	Sterilized using Ethylene Oxide		Bevel Tip
	Keep away from sunlight		Trocar Tip
	Keep dry		Caution: Federal (USA) law restricts this device to sale by or on the order of a physician
	Lot Number		For electronic copy scan QR code or go to www.merit.com/fiu and enter IFU ID Number. For printed copy, call U.S.A. or EU Customer Service.

HOW SUPPLIED

The StabiliT® Introducer is supplied sterile in a peel-open package. The device is sterilized using Ethylene Oxide and is intended for single use only. Do not resterilize. Do not use if the package is open or damaged. Notify manufacturer if damaged.

**STORAGE & HANDLING**

Handle with care. Store in original packaging in a clean, cool, and dry location. Avoid extreme humidity and temperature.

www.merit.com

Manufacturer: Merit Medical Systems, Inc.

EC REP Authorized Representative: Merit Medical Ireland Ltd, Parkmore Business Park West, Galway, Ireland

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