INSTRUCTIONS FOR USE

Intended for use in obtaining core biopsy samples from soft tissue such as breast, kidney, liver, prostate and various soft tissue masses.

Not intended for use in bone.

Do not use if package is damaged.

WARNING: Test firing the needle may damage the cutting cannula and/or stylet tip. Keep fingers or other objects from in front of the needle as it will advance 2.5 cm very rapidly when activated. Never test fire the instrument in an unsafe direction.

STEP 1: PREPARATION OF ACHIEVE FOR USE

To charge the Achieve biopsy device, simply pull back on the charging arm and release twice, first to charge the cutting cannula, the second to charge the inner stylet.

STEP 2: POSITION THE NEEDLE

Position the needle at or near the edge of the region to be sampled.

STEP 3: FIRE THE INSTRUMENT

Two firing sequence methods are available once the Achieve biopsy device is charged.

1. Delay Mode—Depression of the D button releases the inner stylet into the sample area. The operator uses this delay firing method to verify the position of the sample notch in the target area prior to depressing the A button which activates the cutting cannula, capturing the tissue sample.
2. Automatic Mode—The operator can choose the fully sequential mode, omitting the delay feature, simply by depressing the A button first. This releases the stylet and cutting cannula in sequence, capturing the tissue sample.

A = Automatic mode

D = Delay mode (Stylet Only)

STEP 4: BIOPSY SPECIMEN

To retrieve the tissue sample, hold the instrument with the A and D firing buttons facing up, and pull back and release the charging arm once. This action will retract the outer cutting cannula exposing the stylet notch with the tissue sample.

STEP 5: OBTAINING MULTIPLE SAMPLES

After the tissue sample is removed, additional samples from the same patient may be obtained by pulling the charging arm a second time, which will retract the stylet within the cutting cannula, placing the Achieve biopsy device once again in a fully charged position. Repeat Steps 3-5.

In the automatic mode, the Achieve automatic biopsy device is designed to fire the cannula and stylet in rapid succession. Please exercise caution before proceeding as the device will advance 2.5 cm when activated.

CAUTION:Forsingleuseonly.Re-usemayresultinanon-functionalproductorcontributetocrosscontamination.

Only clinicians who have been properly trained in this procedure, including patient preparation, sample preservation and the use of image guidance for percutaneous biopsy, should use this device.

	Use By: YYYY-MM-DD
	Lot Number
	Catalog Number
	Sterilized Using Ethylene Oxide
	Do Not Use if Package is Damaged
	Single Use
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician
	For electronic copy scan QR code or go to www.merit.com/ifu and enter IFU ID Number. For printed copy, call U.S.A or E.U Customer Service.
	Medical Device
Sterile Package	Sterile Package
	Manufacturer
	Date of Manufacture: YYYY-MM-DD

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Merit Medical Ireland Ltd, Parkmore Business Park West, Galway, Ireland EU Customer Service +31 43 358 82 22