HeartSpan®

Transseptal Needles

Transseptal Needle and Stylet Set

Directions for Use

Read instructions for use in its entirety prior to use.

Single use only disposable medical device.

Rx Only: Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

Device Description

The transseptal needle kit consists of an outer needle cannula and an inner stylet. The needle is comprised of flexible thin walled tubing with an ergonomic hub and stopcock attached to the proximal end. The stylet consists of a solid wire that when inserted in the needle protrudes beyond the distal tip of the cannula.

Indications for Use

The transseptal needle is used in conjunction with a transseptal catheter and/or introducer to create the puncture in the atrial septum to allow left heart catheterization procedure to occur through the right atrium.

Contraindications

- Left atrial thrombus or tumor
- Dilated aortic root
- Continual anticoagulation
- Inability to lie flat
- Substantial deformity of the spine or chest
- Marked atrial enlargement
- Distorted anatomy due to congenital heart disease
- Previous intra-septum patch

Warnings

- For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which in turn may result in patient injury, illness or death. This single-use product is not designed or validated to be reused. Reuse may cause a risk of cross-contamination, affect the measurement accuracy, system performance, or cause a malfunction as a result of the product being physically damaged due to cleaning, disinfection, re-sterilization, or reuse.
- Only those physicians who specialize in the practice of invasive cardiology techniques should use this device. The device should be restricted for use by those specialists trained to perform transseptal procedures.
- Maintain continuous pressure monitoring and repeated biplane fluoroscopy during positioning.
- Caution should be used in patients with small left atrium, to avoid left atrial wall puncture.
- The transseptal needle should never be advanced until the catheter is positioned correctly on the atrial septum.

- Always ensure that the transseptal needle has clearly entered the left atrial cavity by confirming distinct left atrial pressure and fluoroscopy of the needle tip before advancing the dilator, sheath or catheter.
- Do not remove a dilator, sheath or catheter that has been inadvertently advanced into the pericardial space until the patient is in surgery.
- Do not reuse this device. The device must be discarded after one use using acceptable medical practices and applicable local, state and federal laws and regulations.
- Follow instructions accompanying the transseptal introducers.

Precautions

- Store in a cool, dark and dry place.
- Inspect all components prior to use.
- If resistance is met while advancing or withdrawing the introducer or guidewire, fluoroscopy should be used to determine the cause.
- Prior to use, ensure the appropriate catheter is being used with the transseptal needle.

Adverse Events

In addition to the complications associated with any cardiac catheterization, the following could occur during a transseptal catheterization:

- Inferior vena cava puncture
- Aorta puncture
- Atrial free wall puncture
- Coronary sinus puncture
- Arterial embolism from thrombus at the puncture site
- Tamponade
- Residual atrial septal defects
- Atrial arrhythmia

Directions for Use

1. To determine the size and location of the left atrial septum, it may be helpful to perform right side angiography.
2. Thoroughly flush the transseptal needle.
3. Inspect all components prior to use for integrity and appropriateness for the particular procedure.
4. The transseptal needle/stylet set is advanced through the transseptal catheter and/or introducer until the tip of the needle is just within the catheter. Note: Ensure the transseptal needle is free to twist and/or rotate without resistance, as it is advanced to this position.
5. Retract the stylet prior to puncturing the interatrial septum.
6. The tip of the transseptal needle is positioned against the interatrial septum.
7. The transseptal needle is advanced through the septum into the left atrium once the correct location is confirmed using fluoroscopy (preferably biplane).
8. The successful puncture into the left atrium is confirmed by observing a left atrial pressure tracing. If an incorrect pressure tracing is observed, a small amount of contrast media is injected to identify the positioning.
9. The transseptal catheter is advanced over the transseptal needle into the left atrium.

10. The transseptal needle is slowly removed.

Shelf Life: 3 years.

	Contents non-pyrogenic and sterilized using ethylene oxide
	Do not resterilize
	Single use only / Do not reuse
	Manufacturer
	Manufacturing date
	Lot number
	Expiration date
	Number of
	Do not use if packaging is damaged!
	Catalog number
	Caution, read instructions for use prior to use
	Product number
	Proximal gage
	Distal gage
	Store in a cool, dark and dry place

Manufacturer:
Merit Medical System, Inc.
1600 West Merit Parkway South Jordan, Utah 84095 U.S.A.
1-801-253-1600 U.S.A
Customer Service 1-800-356-3748

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