BLUNT OBTURATOR

EN	INSTRUCTIONS FOR USE
FR	MODE D'EMPLOI
IT	ISTRUZIONI PER L'USO
DE	GEBRAUCHSANWEISUNG
ES	INSTRUCCIONES DE USO
PT	INSTRUÇÕES DE UTILIZAÇÃO
PTBR	INSTRUÇÕES DE USO
NL	GEBRUIKSINSTRUCTIES
SV	BRUKSANVISNING
DA	BRUGSVEJLEDNING
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
TR	KULLANMA TALİMATLARI
BG	ИНСТРУКЦИИ ЗА УПОТРЕБА
PL	INSTRUKCJE UŻYTKOWANIA
CS	NÁVOD K POUŽITÍ
HU	HASZNÁLATI UTASÍTÁS
RU	РУКОВОДСТВО ПО ЭКСПЛУАТАЦИИ

BLUNT OBTURATOR

INSTRUCTIONS FOR USE

DESCRIPTION

The blunt obturator is an access device with a handle and tapered tip for percutaneous insertion into the desired region to be treated. The blunt obturator is provided sterile and is a single use device.

INDICATIONS FOR USE

The blunt obturator is intended for use with the SAVI® applicator. The blunt obturator is intended to provide a pathway through which the SAVI applicator is placed into the treatment area.

CONTRAINDICATIONS

Use of the blunt obturator in regions that are in generally poor condition (e.g., ulcerated) is not recommended.

WARNINGS

The blunt obturator should only be used by physicians trained in brachytherapy techniques.

The physician is responsible for its proper clinical use.

The blunt obturator is shipped sterile and MUST NOT BE RESTERILIZED.

The blunt obturator is for SINGLE USE only.

Do not use if the package is open or damaged.

PLACEMENT OF BLUNT OBTURATOR

1. Through the skin nick and using ultrasound guidance, advance the obturator and sheath provided in the tray into the cavity. Always align the long axis of the obturator to the long axis of the cavity.

2. Retract the obturator, leaving the sheath in place.

3. Insert the SAVI applicator through the sheath until the distal tip of the device contacts the far wall of the cavity.

4. Carefully retract the sheath, leaving the SAVI applicator in place.

5. Follow Instructions for Use included with the SAVI applicator

DISPOSAL PROCEDURES

It is the user's responsibility to dispose of all parts according to local regulations. Upon disposal of the product or related parts, it is suggested to:

- Disinfect the parts prior to disposal
- Dispose of sharp objects in a safe manner

STORAGE

Store at room temperature. Avoid storing the blunt obturator at conditions of excessive heat or humidity.

COMPONENTS

Description	Quantity
Blunt obturator w/ Sheath	1

Merit Medical Systems, Inc.
6 Journey, Suite 125
Aliso Viejo, CA 92656 USA

SYMBOL	DESIGNATION
[REF]	Catalog Number
[STERILE EO]	Sterilized Using Ethylene Oxide
(X)	Do Not Use if Package is Damaged
(2)	Single Use
R ONLY	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician
(2) STERILE	Do Not Re-sterilize
⚠	Caution: Consult accompanying documents
▲	Manufacturer
MD	Medical Device.
Sterile Package	Sterile Package.

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