# **Q50® PLUS STENT GRAFT BALLOON CATHETER**

## **Instructions for Use**

**STERILE | EO**

Contents sterile unless enclosed package  
has been opened or damaged.  
Sterilized by Ethylene Oxide (EO).

**CAUTION: Federal (USA) law restricts this device to sale by or on the order of  
a physician.**

<b>English (English).....</b>	<b>2</b>
<b>Bulgarian (Български) .....</b>	<b>6</b>
<b>Czech (Česky) .....</b>	<b>11</b>
<b>Danish (Dansk).....</b>	<b>15</b>
<b>Dutch (Nederlands).....</b>	<b>19</b>
<b>French (Français).....</b>	<b>23</b>
<b>German (Deutsch).....</b>	<b>27</b>
<b>Greek (Ελληνικά) .....</b>	<b>31</b>
<b>Hungarian (Magyar).....</b>	<b>36</b>
<b>Italian (Italiano) .....</b>	<b>40</b>
<b>Polish (Polski) .....</b>	<b>44</b>
<b>Portuguese (Português) .....</b>	<b>49</b>
<b>Portuguese-Brazil (Português do Brasil) .....</b>	<b>53</b>
<b>Russian (Русский).....</b>	<b>57</b>
<b>Spanish (Español).....</b>	<b>62</b>
<b>Swedish (Svenska).....</b>	<b>66</b>
<b>Turkish (Türkçe).....</b>	<b>70</b>

### 1.DEVICE DESCRIPTION


The Stent Graft Balloon Catheter is designed for temporary occlusion of large vessels and to assist in the expansion of stent grafts used for the treatment of aortic aneurysms. Ballooning of a stent graft may improve modeling of the graft material and fixation of the stent graft to the vessel wall. Sub-optimal expansion of stent grafts may also be improved by inflating the balloon at the site of the stent. The Stent Graft Balloon Catheter has been tested to temporarily occlude vessels up to 41 mm in diameter.

The Stent Graft Balloon Catheter is an over-the-wire (OTW) tri-lumen catheter with a compliant polyurethane balloon having a maximum diameter of 50 mm at 60 cc. Two lumens inflate and deflate the balloon while one lumen is reserved for guidewire passage. The catheter is offered in two (2) usable lengths, 100 cm and 65 cm.

This device is designed to accommodate a 0.038" diameter (or smaller) guidewire. Two (2) radiopaque marker bands are placed within the balloon at each end [40 mm apart] to facilitate balloon placement prior to inflation. The catheter can accommodate a 0.038" diameter (or smaller) guidewire and is compatible with 12 Fr (or larger) introducer sheaths. The catheter has an extension tube with a stopcock in order to facilitate handling and fluid control. The device models are listed below:

Catalogue Number	Catheter Effective Length
Q50-65P	65 cm
Q50-100P	100 cm

There are risks involved with any medical procedure. Both physician and patient should fully understand those risks associated with surgery, and additional new risks associated specifically with the use of this endoluminal device.

**CAUTION:** Only physicians trained in vascular surgery, interventional radiology or cardiology, and who have completed training or have experience with stent grafts, balloon catheters and associated devices should consider using this device.

**CAUTION:** Read the entire Instructions for Use manual prior to using the device.

### 2. INDICATIONS FOR USE

The Stent Graft Balloon Catheter is intended for temporary occlusion of large vessels or to expand vascular prostheses.

### 3. CONTRAINDICATIONS

The Stent Graft Balloon Catheter is contraindicated in patients who:

- Are contraindicated to contrast media or anticoagulants
- Have an arterial entry site that cannot accommodate a 12 Fr (min) introducer sheath
- Are minors <18 years old
- Are pregnant

### 4. WARNINGS

- The Stent Graft Balloon Catheter is supplied STERILE and for single use only. Do not reprocess or re-sterilize. Reprocessing and re-sterilizing could increase the risk of patient infection and of compromised device performance.
- The catheter should only be manipulated and inflated/deflated while observing under fluoroscopy.
- If resistance is encountered at any time during the insertion procedure, do not force passage or torque the catheter. Resistance may cause damage to device, vessel or stent graft. Carefully withdraw the catheter.
- Do not torque or twist the catheter during insertion or withdrawal.
- The catheter should only be advanced or withdrawn over a guidewire.
- Adhere to balloon inflation parameters outlined in the **Balloon Compliance Chart (Table 1)**. Do not exceed a balloon diameter of 50 mm and do not exceed 60 cc inflation volume at 50 mm balloon diameter. Rupture of balloon may occur. Over-inflation may result in damage to vessel wall and/or vessel rupture, or damage to the stent graft.
- Balloon rupture may occur under certain anatomical, procedural, and/or clinical circumstances. It is therefore recommended to have back-up Stent Graft Balloon Catheters on hand.

- Ensure that the balloon is fully deflated before moving the Stent Graft Balloon Catheter.
- When expanding a vascular prosthesis, there is an increased risk of vessel injury and/or rupture, and possible patient death, if the balloon inflation is not completely within the covered (graft fabric) portion of the prostheses.
- Studies indicate that the danger of micro-embolization increases with increased manipulation and/or duration of the procedure.
- Over inflation of the balloon can cause graft tears and/or vessel rupture. Care should be taken when inflating the balloon in vessels, particularly when inflating in the distalmost area of the stent graft, or in calcified, stenotic, and/or otherwise diseased vessels.
- Do not use in the heart or coronary arteries.

- Do not use in the treatment of dissections.
- Hand injections using a 60 cc syringe are recommended. Do not use a pressure inflation device for balloon inflation.
- Do not use a power injector for injection of contrast medium through distal catheter lumen. Rupture may occur.
- QXMédical is not responsible for mis-sizing, misuse or misplacement of the device.

## 5. PRECAUTIONS

- Preparations should be made and a trained vascular surgical team available in the event conversion to open surgery is required.
- Carefully inspect the package and catheter prior to use to verify no damage occurred during shipment. Do not use if the package or catheter is damaged since the sterility or integrity of the device may be compromised and thus increasing the risk of patient infection and device malfunction.
- Use the catheter prior to the **Use By** date specified on the package.
- Do not attempt to pass the catheter through an introducer sheath smaller than 12 Fr. Damage to the device may occur.
- Only use guidewires with a diameter  $\leq 0.038"$ .
- To avoid kinking and damaging the catheter, advance slowly in small increments until the proximal end of the guidewire emerges from the catheter.
- To avoid damage to the catheter or vessel, do not advance or withdraw the device without a guidewire in place.
- Catheter should not be advanced into a vessel or stent graft having a diameter smaller than the catheter outer diameter. Damage to the device or vessel may occur.
- Carefully monitor the patient's blood pressure throughout the procedure.
- If an obstruction in the vessel (e.g., a tortuous bend, stenosis, calcification, etc.) prevents advancement of the catheter, use standard techniques to try to dilate and/or straighten the vessel before continuing to advance the catheter.
- It is recommended that back-up Stent Graft Balloon Catheters are made available.
- Balloon is highly compliant. Inflate slowly. Do not over-inflate balloon when modeling graft in vessels. Operator should visualize the stent graft at all times during balloon inflation to detect any movement of the stent graft. Use special care in areas of diseased vessels to avoid rupture or vessel trauma.
- The Stent Graft Balloon Catheter is not intended for use as an angioplasty or dilatation balloon.
- The device is not for use in endovascular prostheses contraindicated for use with compliant balloon catheters.
- When aligning the position of the catheter so that the balloon is in proper position for expansion within vessel or stent graft, pay careful attention to the fluoroscope location in order to avoid parallax or other sources of visualization error.
- Do not use the Stent Graft Balloon Catheter for more than 20 inflation/deflation cycles.
- Do not exceed 50 mm balloon diameter and do not exceed 60 cc inflation volume at 50 mm balloon diameter. Adhere to the balloon inflation guidelines outlined in Table 1.
- The Stent Graft Balloon Catheter has been tested to temporarily occlude vessels up to 41 mm in diameter.

## 6. POTENTIAL COMPLICATIONS/ADVERSE EFFECTS

Complications may occur with the use of any stent graft balloon catheter or during any catheterization procedure. Therefore, only physicians trained in vascular surgery, interventional radiology or cardiology, and who have completed training

or have experience with access/support catheters and associated devices should consider using this device. Possible complications associated with this type of procedure include but are not limited to the following:

### CLINICALLY RELATED:

- Vessel perforation, dissection, rupture or injury
- Occlusion at some locations may cause arrhythmia
- Paresthesia
- Drug reactions
- Entry site infection
- Entry site hematoma
- Cardiac events
- Respiratory failure
- General malaise
- Arterial thrombosis and/or embolism
- Hemorrhage
- Stroke
- Aneurysm rupture
- Renal complications
- Death

### DEVICE RELATED:

- Balloon catching inside of graft during catheter movement
- Stent graft migration caused by device
- Balloon rupture
- Inability to inflate/deflate balloon
- Inability to insert guidewire
- Inability to withdraw catheter from introducer

## 7. PACKAGING, STERILIZATION AND STORAGE

The device has been sterilized using ethylene oxide (EO) and is supplied sterile and non-pyrogenic. The package label indicates the **Use By** date.

Do not use the device after the **Use By** date. The device sterility and integrity may be compromised and possibly result in patient infection and device malfunction.

The device should be stored in cool dry place.

## 8. RECOMMENDED ITEMS

Each Stent Graft Balloon Catheter package includes the following:

- Single-use sterile disposable catheter with extension tube and stopcock mounted on a paperboard card.

Materials required but not provided are:

- Introducer Sheaths (12 Fr Minimum size)
- Guidewires (0.038" Maximum diameter)
- 60 cc syringes
- Diluted contrast solution (75% sodium chloride / 25% renograpin) is recommended for balloon inflation/deflation
- Additional Stent Graft Balloon Catheters

It is also recommended that a freely-angled C-Arm with high resolution fluoroscopy, high quality angiography and Digital Subtraction Angiography (DSA) be used during the procedure.

## 9. DEVICE PREPARATION

- 9.1 Carefully inspect the package and catheter prior to use to verify no damage occurred during shipment. Do not use catheter if either the catheter or packaging is damaged or compromised.

- 9.2 If the packaging is free of damage, carefully open the outer pouch and introduce the inner pouch

- (with catheter) to the sterile field using sterile techniques.
- 9.3 Open the inner pouch and carefully disengage the catheter, manifold and extension tube assembly from the paperboard card. Carefully inspect the catheter for any signs of damage. If damaged, please discard and use another catheter.
- NOTE:** Patient and entry site preparation to be performed prior to device preparation. Prepare patient in accordance with standard techniques; including proper administration of anticoagulation and antiplatelet medication. Using standard techniques prepare entry site, including placement of the introducer sheath.
- 9.4 Fill 60 cc syringe with 10 cc to 15 cc of heparinized saline solution.
- 9.5 Attach the syringe to the guidewire lumen port and flush the guidewire lumen.
- 9.6 Remove the syringe from the guidewire lumen port, fill with 15 cc of heparinized saline solution and connect it to the stopcock. Ensure that the stopcock is OPEN and draw vacuum.
- 9.7 While maintaining the balloon under vacuum, gently remove the protective sleeve from the balloon by twisting the sleeve in one direction and pulling it off the balloon.
- 9.8 Partially inflate the balloon and draw a vacuum to remove air from the catheter.
- 9.9 Air removed from the catheter and accumulated in the syringe may be expelled by turning the stopcock OFF, removing the syringe and ejecting the air.
- 9.10 Attempt to remove all the air in the balloon by repeating steps 9.8 and 9.9 as required. Some changes in the catheter orientation may be necessary to vent all the air. Ensure no air is introduced to the catheter when reconnecting the syringe to the stopcock.
- 9.11 Draw vacuum, turn the stopcock OFF and remove the syringe. Set the catheter aside until needed.
- NOTE:** If air bubbles continue to appear in the syringe during catheter preparation or if liquid leakage is seen at any point on the catheter (including the balloon bonds, tip or manifold), please discard the catheter and use a new one.
- ## 10. DEVICE USAGE
- 10.1 Using standard techniques, introduce the guidewire through the introducer sheath, ensuring it is advanced beyond the balloon target location.
- 10.2 Prior to introducing the catheter, determine the amount of diluted contrast (75% sodium chloride / 25% renograpin) needed to inflate the balloon to the desired diameter. Refer to the balloon inflation parameters outlined in the Balloon Compliance Chart (Table 1).
- WARNING:** Over-inflation may result in damage to vessel wall and/or vessel rupture, or damage to the stent graft.
- 10.3 Place the prepared catheter over the previously inserted guidewire by threading the end of the guidewire through the tip of the catheter.
- 10.4 Verify that the balloon is completely deflated and slowly advance through the introducer sheath.
- 10.5 Using fluoroscopic imaging, carefully advance the catheter to the desired location in the vasculature or stent graft using the markerbands (distal & proximal end of the balloon) for visual guidance. When using the Stent Graft Balloon Catheter to expand a vascular prosthesis, use the radiopaque markers to ensure that the entire balloon is positioned within the prosthesis.
- CAUTION:** When aligning the position of the catheter so that the balloon is in proper position for expansion within vessel or stent graft, pay careful attention to the fluoroscope location in order to avoid parallax or other sources of visualization error. When the Stent Graft Balloon Catheter is positioned properly, it may be appropriate to decrease the patient's blood pressure to approximately 80 mm Hg to facilitate ballooning.
- CAUTION:** Failure to lower the blood pressure during thoracic cases may cause balloon movement and may result in graft misplacement.
- 10.6 Fill the 60 cc syringe with the desired amount of diluted contrast solution and attach it to the stopcock ensuring that no air is introduced to the balloon catheter. Open the stopcock and inflate the balloon at the target location. Continuously monitor fluoroscope screen, watching for stent movement. Proper modeling should show very slight outward expansion of stent with balloon inflation. Be careful not to over inflate – stop right after stent expansion is evident. Table 1 below is a guideline for determining the volume of diluted contrast solution required to obtain a given balloon expansion diameter.
- | Recommended Inflation Volume | Balloon Diameter |
|------------------------------|------------------|
| 3 cc                         | 10 mm            |
| 6 cc                         | 20 mm            |
| 16 cc                        | 30 mm            |
| 32 cc                        | 40 mm            |
| 44 cc                        | 45 mm            |
| 60 cc                        | 50 mm*           |
- \* Maximum Inflation Diameter -- DO NOT EXCEED 50 mm Balloon Diameter and DO NOT EXCEED 60 cc inflation at 50 mm balloon diameter.
- CAUTION:** The above chart is only a guide. Balloon expansion should be carefully monitored with the use of fluoroscopy.
- NOTE:** The balloon has been tested to temporarily occlude vessels up to 41 mm in diameter. Fully deflate balloon and verify using fluoroscopy that the balloon is deflated. Move the balloon; inflate and deflate the balloon to smooth wrinkles in the graft material and to enhance contact between the stent graft and the vessel wall.
- 10.8 If necessary, repeat until all target areas of the stent graft have been balloon-modeled. Inflate the balloon with sufficient pressure to firmly embed the stent graft against the vessel.
- CAUTION:** Balloon is highly compliant. Inflate slowly. Do not over-inflate balloon when modeling

- graft in vessels. Operator should visualize the stent graft at all times during balloon inflation to detect any movement of the stent graft. Use special care in areas of diseased vessels to avoid rupture or vessel trauma.
- 10.10 If balloon pressure is lost and/or balloon rupture occurs, deflate the balloon and remove balloon and sheath as a unit.
- 10.11 When balloon inflations are complete, draw a vacuum in the balloon and verify that the balloon is fully deflated.
- 10.12 Maintain vacuum on the balloon and withdraw the Stent Graft Balloon Catheter back through the introducer sheath. Use fluoroscopic imaging to ensure that the stent graft does not move and to track the movement of the Stent Graft Balloon Catheter throughout the withdrawal.
- 10.13 If resistance is met during withdrawal, apply negative pressure with a larger syringe before proceeding. If resistance continues, remove balloon and sheath as a unit.
- WARNING:** Do not torque or twist catheter during insertion or withdrawal.
- 10.14 Close entry site with standard surgical closure technique.

## 11. DISPOSAL

After use, this product is a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable hospital, local, state and federal laws and regulations.

## WARRANTY DISCLAIMER AND LIMITATION OF LIABILITY

ALTHOUGH THE STENT GRAFT BALLOON CATHETER HAS BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, QXMÉDICAL, LLC HAS NO CONTROL OVER THE CONDITIONS UNDER WHICH THE STENT GRAFT BALLOON CATHETER IS USED. QXMÉDICAL, LLC THEREFORE DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE STENT GRAFT BALLOON CATHETER INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE AND ANY IMPLIED WARRANTY ARISING FROM COURSE OF PERFORMANCE, DEALING, USAGE OR TRADE. QXMÉDICAL, LLC SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE STENT GRAFT BALLOON CATHETER, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND QXMÉDICAL, LLC TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE STENT GRAFT BALLOON CATHETER.

Descriptions and specification appearing in QXMédical, LLC's printed materials, including this publication, are informational only and are not made or given as a warranty of the Stent Graft Balloon Catheter in any way. The exclusions and limitations set forth above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this Warranty Disclaimer is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this Warranty Disclaimer shall not be affected, and all rights and obligations shall be construed and enforced as if this Warranty Disclaimer did not contain the particular part or term held to be invalid.

Manufactured by QXMédical, LLC.

<p><b>Manufacturer:</b></p> <p><b>QXMédical, LLC</b> 2820 Patton Road St. Paul, MN 55113 U.S.A. <a href="http://www.qxmedical.com">www.qxmedical.com</a></p>
<p><b>Authorized Representative:</b></p> <p><b>Emergo Europe</b> Prinsessegracht 20 2514 AP The Hague The Netherlands Tel: +31 (0)70 345 8570</p>
<p><b>Worldwide Distributor:</b></p> <p><b>Merit Medical Systems, Inc.</b> 1600 West Merit Parkway, South Jordan, Utah 84095 U.S.A. 1-801-253-1600</p>

Q50 is a registered trademark of QXMédical, LLC.

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QXMédical, LLC tarafından üretilmiştir.
Q50, QXMédical, LLC'nin tescilli ticari markasıdır.