# BIG60 ALPHA™

EN	INSTRUCTIONS FOR USE
FR	MODE D'EMPLOI
IT	ISTRUZIONI PER L'USO
DE	GEBRAUCHSANWEISUNG
ES	INSTRUCCIONES DE USO
PT	INSTRUÇÕES DE UTILIZAÇÃO
PTBR	INSTRUÇÕES DE USO
NL	GEBRUIKSAANWIJZING
SV	BRUKSANVISNING
NO	BRUKSANVISNING
DA	BRUGSANVISNING
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
TR	KULLANMA TALİMATLARI
PL	INSTRUKCJA UŻYTKOWANIA
CS	NÁVOD K POUŽITÍ
BG	ИНСТРУКЦИИ ЗА УПОТРЕБА
HU	HASZNÁLATI ÚTMUTATÓ
RU	ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ
LV	LIETOŠANAS INSTRUKCIJA
LT	NAUDOJIMO INSTRUKCIJOS
RO	INSTRUCȚIUNI DE UTILIZARE
SK	NÁVOD NA POUŽITIE
HR	UPUTE ZA UPORABU
ET	KASUTUSJUHEND
SL	NAVODILA ZA UPORABO

# BIG60 ALPHA™

## INSTRUCTIONS FOR USE

### DESCRIPTION

The BIG60 Alpha™ Inflation Device is a 60ml disposable inflation device capable of producing a maximum pressure of 12 ATM/BAR, fitted with a threaded plunger assembly with a lock/release bar, a flexible high-pressure tube, and a three-way medium pressure stopcock. The BIG60 Alpha™ Inflation Device is designed to generate and monitor inflation pressures over a range of -.68 to 12ATM/BAR (-10 to 176 psi with an accuracy of 3% ATM/BAR).

### INDICATIONS FOR USE/CLINICAL BENEFIT

This inflation device is used to inflate and deflate an angioplasty balloon or other interventional device, and to measure the pressure within the balloon.

**NOTE:** This device has not been cleared for dispensing fluid into the body.

### PRECAUTIONS

- Inspect the BIG60 Alpha™ Inflation Device and packaging for damage prior to use. Do not use product if opened or damaged. Confirm the product is consistent with the package label. Contact Customer Service to report and replace damaged product.
- For single patient use only. Do not reuse, reprocess, or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness, or death. Reuse, reprocessing, or resterilization may also create a risk of contamination of the product and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the product may lead to injury, illness, or death of the patient.
- Graduations on device barrel are for reference only and are not intended for exact measurement.
- In the EU, any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the applicable Member State.

### POTENTIAL COMPLICATIONS


The potential exists for complications which may include, but are not limited to the following:

- Obstruction non vascular
- Pulmonary event
- Soft tissue injury

### INSTRUCTIONS FOR USE

#### DEVICE PREPARATION

1. Prior to use, free the plunger tip by twisting the device plunger/handle 360 degrees clockwise.
2. To prep device, simply aspirate up to 60ml of contrast solution into the syringe by pulling back on the device handle.  
**CAUTION:** Inspect the device tubing and stopcock (if used) to ensure there is no air in the system.
3. Remove any excess air by orienting the device upwards, pushing the handle against a table or other solid surface until the device plunger tip is oriented with the black arrow printed on the device barrel (approximately 35ml).

### **ATTACHING THE INFLATION DEVICE TO THE BALLOON**

**NOTE:** Refer to the manufacturer's directions accompanying the balloon dilatation catheter or other interventional device for specific information on use, maximum inflation pressure, precautions, and warnings for that device.

1. Prepare and test the balloon catheter according to the catheter manufacturer's directions for use.
2. Create a fluid-fluid connection between the balloon and the inflation syringe extension tube, connect the luer connectors securely.
3. Squeeze the trigger and pull back on the plunger handle to apply a vacuum to the balloon.

### **BALLOON INFLATION AND DEFLATION**

1. To inflate the balloon, squeeze the trigger to allow the plunger to return to a resting position (0 ATM/BAR or PSI). Release grip on the trigger, which will lock the plunger into position.

2. To increase pressure, rotate handle clockwise until the desired pressure is achieved.

**CAUTION:** Do not exceed 12ATM.

**NOTE:** Loss of pressure may indicate a leak in the system.

3. To deflate the balloon (rotate the handle counterclockwise to relieve pressures above 6 ATM) squeeze the trigger and pull the handle back to generate the desired negative pressure. Release the grip on the trigger to lock the plunger in the negative pressure position.

### **SAFE DISPOSAL**

- Follow local guidelines and practices for the disposal of infected waste products.

**STORAGE:** Store in a cool, dry place.

### **WARRANTY**

The manufacturer warrants that reasonable care has been used in the design and manufacture of this device. This warranty is exclusive and manufacturer makes no other representations or warranties of any kind to customers, its end users, or to any third parties with respect to the device and hereby expressly disclaims any and all other warranties, express or implied, statutory or otherwise, including, but not limited to, infringement and the implied warranties of merchantability and fitness for a particular purpose, even if manufacturer is aware of such purpose. Handling and storage of this device, as well as other factors relating to the patient, diagnosis, treatment, implant procedures, and other matters beyond the control of the manufacturer, directly affect the device and the results obtained from its use. The manufacturer's obligation under this warranty is limited to the replacement of the device. Under no circumstances shall manufacturer be liable to customer or any other person or entity for any punitive, special, incidental or consequential damages directly or indirectly arising from the use of this device. The manufacturer neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device. This warranty shall not apply, and manufacturer assumes no liability with respect to, devices that have been (i) modified, changed, altered, misused, mishandled, repaired, reused, reprocessed, refurbished or resterilized; (ii) subjected to improper maintenance, testing or storage, accident, tampering, or inadequate protection against shock, vibration, excessively high or low temperatures, overpressure, or physical, environmental or electrical stress; (iii) been used outside the approved "Indications for Use" as cleared by the relevant competent authority, used contrary to the use outlined in the device specifications, or in an application or environment for which such device was not designed or contemplated; or (iv) distributed or used contrary to applicable federal, state, local or regulatory standards.

**Rx ONLY CAUTION:** Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

	Catalog number
	Batch code
	Date of Manufacture: YYYY-MM-DD
	Use by date: YYYY-MM-DD
	Single use.
	Do not resterilize
	Non-pyrogenic
	Do Not Use If Package is Damaged and Consult Instruction for Use
<b>Rx ONLY</b>	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Sterilized using ethylene oxide
	Temperature limitations
	Caution
	Consult Instructions for Use For electronic copy scan QR code, or go to <a href="http://www.merit.com/ifu">www.merit.com/ifu</a> and enter IFU ID. For printed copy, call U.S.A. or EU Customer Service
	Medical Device
	Single Sterile Barrier System
	Unique Device Identifier
	Authorized Representative in European Community
	Manufacturer

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Merit Medical Systems, Inc.
South Jordan, Utah 84095 U.S.A.
U.S.A. Customer Service 1-800-356-3748
Merit Medical Ireland Ltd