Cultura™ TRANSPORT SYSTEM

Sterile Flocked Collection Device with Transport Vial and Biohazard Specimen Bag

EN	INSTRUCTIONS FOR USE
FR	MODE D'EMPLOI
IT	ISTRUZIONI PER L'USO
DE	GEBRAUCHSANWEISUNG
ES	INSTRUCCIONES DE USO
PT	INSTRUÇÕES DE UTILIZAÇÃO
PTBR	INSTRUÇÕES DE USO
NL	GEbruIKSAANWIJZING
SV	ANVÄNDNINGSSINSTRUKTIONER
NO	BRUKSANVISNING
DA	BRUGSANVISNING
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
TR	KULLANMA TALİMATLARI
PL	INSTRUKCJA UŻYTKOWANIA
CS	NÁVOD K POUŽITÍ
BG	ИНСТРУКЦИИ ЗА УПОТРЕБА
HU	HASZNÁLATI UTASÍTÁS
RU	ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ
LV	LIETOŠANAS PAMĀCĪBA
LT	NAUDOJIMO INSTRUKCIJOS
RO	INSTRUCȚIUNI DE UTILIZARE
SK	NÁVOD NA POUŽÍVANIE
HR	UPUTE ZA UPORABU
ET	KASUTUSJUHEND
SL	NAVODILA ZA UPORABO

Cultura™ TRANSPORT SYSTEM

English

The Viral Transport Medium (VTM) consists of Hank's Balanced Salt Solution (HBSS) enriched with proteins and sugars with a neutral pH and pH indicator. The VTM contains antibiotics to inhibit overgrowth of bacteria, yeast and fungi, maintain cellular integrity and encourage preservation of viruses when properly stored and has been evaluated for transport of SARS-CoV-2. Prior to use, vials should be stored at 2-25°C. After collection, the transport tube containing the specimen can be stored for up to 120 hours at 2-25°C, which provides convenience for transportation to the laboratory and storage. Samples can be stored at -70°C for archival purposes.

INDICATIONS FOR USE

The Merit Cultura™ Transport System is intended for transport to the analysis laboratory of clinical specimens for subsequent diagnostic/identification techniques.

CONTRAINDICATIONS

None known

LIMITATIONS

1. The Cultura VTM is intended to be used with specimen collection swabs. The use of medium, tubes, or swabs from any other source may compromise performance.
2. PCR studies were performed in triplicate using the Cultura VTM and swabs with positive and negative SARS-CoV-2 controls to confirm nucleic acid detection.
3. Testing has been performed in triplicate to show the suitability of the VTM for the preservation of nucleic acids (RNA and DNA) for downstream nucleic acid extraction and molecular assay testing and experiments when VTM and collection samples are stored per instructions provided.
4. Although the VTM is formulated to encourage preservation of viruses, virus viability testing has not been performed.

Results obtained largely depend on proper and adequate specimen collection as well as the promptness with which the specimens are transported to the laboratory and analyzed.

WARNINGS

1. This product is for single use only; reuse may cause a risk of infection and/or inaccurate results.
2. Do not re-pack.
3. Not suitable for any other application than intended use.
4. The use of this product in association with a rapid diagnostic kit or with diagnostic instrumentation should be previously validated by the user.
5. Instructions for use must be followed carefully. The manufacturer cannot be held responsible for any unauthorized or unqualified use of the product.
6. To be handled by trained personnel only.
7. It must be assumed that all specimens contain infectious micro-organisms; therefore all specimens must be handled with appropriate precautions. After use, tubes and swabs must be disposed of according to laboratory regulations for infectious waste.

CONTENTS

The Cultura™ Transport System is ready for use and requires no further preparation. The packaging includes:

- 1 Vial with 3mL of Viral Transport Medium

PRECAUTIONS

1. Carefully follow the instructions for use.
2. Product can not be reused.
3. Sterile gloves and protective clothing and eyewear should be worn when collecting and handling microbiology specimens.

STORAGE

Store product in a cool, dry place.

Storage temperature for the vial with reagent is 2-25°C. Do not overheat or freeze prior to use.

ADVERSE REACTIONS

None known

RX ONLY

Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician trained and/or experienced in the use of this device, specific therapeutic or diagnostic procedure.

INSTRUCTIONS FOR USE

1. Collect the sample. During sampling, the swab end (tip) shall only come in contact with the suspected infection, so as to reduce contamination risks.
2. After patient sampling, immediately place the tip of the swab into the tube containing the VTM.
3. Break the swab off into vial by aligning breakpoint with upper edge of vial and bend until swab breaks off at mark, leaving the swab tip in the transport tube. Discard the stick and recap the transport tube containing the swab sample.
4. After collection, standard operating protocols for clinical specimen handling and preservation should be followed. The transport tube containing the swab with nucleic acids from the specimen can be stored for up to 120 hours at 2-25°C, which provides convenience for transportation to the laboratory and storage. It can be stored at -70°C for archival purposes.
5. After the stored sample's target RNA or DNA is extracted by commercial extraction kits, various gene detection and molecular assay tests can be performed.

Specimens collected for clinical investigations should be collected and handled following published manuals and guidelines. Specific requirements for the shipment and handling of specimens should be in full compliance with state and federal regulations. Shipping of specimens within medical institutions should comply with internal guidelines of the institution.

CAUTIONS

1. In the laboratory, wear protective gloves and other protection commensurate with universal precautions when handling clinical specimens. Observe biosafety recommendations when handling or analyzing patient samples.
2. Condition, timing, and volume of specimen collected for clinical investigation are significant variables in obtaining reliable results. Follow recommended guidelines for specimen collection.
3. Check the package before using it. If damaged, do not use.
4. Vial container may become brittle at cold temperatures and may crack if dropped. Handle with care.

WASTE DISPOSAL

VTM must be disposed of in accordance with applicable regulations. Dispose of any used or contaminated disposable materials following procedures for infectious or potentially infectious products. It is the responsibility of each laboratory to handle waste and effluents produced according to their nature and degree of hazard and to treat and dispose of them (or have them treated and disposed of) in accordance with applicable regulations.

REAGENTS

The VTM formulation includes proteins for virus stabilization, antibiotics and antimycotics to prevent overgrowth of bacterial and fungal flora and a buffer solution to maintain a neutral pH.

Component	Concentration
HBSS Solution	10 g/l
pH Indicator	< 1 g/l
FBS	20 g/l
Antibiotics	< 1 g/l

 Single Use		 Date of Manufacture
 Do not use if package is damaged and consult instructions for use	 For In Vitro Diagnostic Use	
 Storage Temperature, (2-25°C)	 Use by Date	

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Merit Medical Ireland Ltd