INSTRUCTIONS FOR USE
MODE D'EMPLOI
ISTRUZIONI PER L'USO
GEBRAUCHSANWEISUNG
INSTRUCCIONES DE USO
INSTRUÇÕES DE UTILIZAÇÃO
INSTRUÇÕES DE USO
GEBRUIKSINSTRUCTIES
ANVÄNDNINGSSINSTRUKTIONER
BRUKSANVISNING
BRUGSANVISNING
ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
KULLANMA TALIMATI
INSTRUKCJA UŻYTKOWANIA
NÁVOD K POUŽITÍ
ИНСТРУКЦИИ ЗА УПОТРЕБА
HASZNÁLATI UTASÍTÁS
ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ
LIETOŠANAS PAMĀCĪBA
NAUDOJIMO NURODY
INSTRUCȚIUNI DE UTILIZARE
NÁVOD NA POUŽITIE
UPUTE ZA UPORABU
KASUTUSJUHEND
NAVODILA ZA UPORABO
KÄYTTÖOHJEET


INSTRUCTIONS FOR USE

DESCRIPTION

The Prelude Roadster™ Guide Sheath Introducer is designed to perform as a guiding sheath and/or standard sheath introducer. The sheath utilizes a flat coiled wire to provide kink resistance as well as flexibility. The sheath comes in 4-8 Fr and in lengths of 30, 45, 65, and 90 cm. The sheath will be offered in multiple shapes: straight, multipurpose, hockey stick and renal. The distal 35cm of the sheath is hydrophilic coated. The sheath is radiopaque and has a marker band located approximately 5mm from the distal tip of the sheath. NOTE: The tip of the sheath is also radiopaque so position of the marker band to the tip can be visually verified under fluoroscopy imaging.


The product consists of the following components:

- One (1) Sheath Introducer
- One (1) Vessel Dilator
- One (1) Hemostasis Valve Adapter (HVA)

*Image above shows the Prelude Roadster Guide Sheath Straight

The Prelude Roadster Guide Sheath has a marker band located approximately 5mm from the distal tip.

INDICATIONS FOR USE

The Prelude Roadster Guide Sheath Introducer is indicated to be used for the introduction of interventional and diagnostic devices into the peripheral (and coronary) vasculature.

INTENDED PURPOSE


The Prelude Roadster Guide Sheath Introducer is intended to provide access and facilitate the percutaneous introduction of various devices into veins and/or arteries while maintaining hemostasis for a variety of diagnostic and therapeutic procedures.

CONTRAINDICATIONS

Radial access is contraindicated if there is an abnormal Allen's Test, radial pulse, or insufficient dual arterial supply.

WARNINGS

- Extreme care should be taken when manipulating a sheath and wire combination to prevent possible intravascular tissue damage. If resistance is felt during advancement, manipulation, or removal, stop immediately and confirm wire position under fluoroscopy.
- Do not use device with a power injector.
- Prior to beginning radial artery access, an assessment such as the Allen's test should be performed to assess the presence/adequacy of dual arterial circulation to the hand.
- Do not leave the introducer in place for extended periods of time without a catheter or obturator to support the cannula wall.
- Appropriate flushing protocols should be utilized to prevent thrombus formation during procedural use.
- The sheath is radiopaque and has a marker band located approximately 5mm from the distal tip of the sheath. When inflating a balloon at, or close to the sheath tip, do not inflate it inside the distal end of the sheath

CAUTIONS

- Read instructions prior to use.
- Store in a cool, dry place.
- This device is intended for single use only. Do not reuse or resterilize.
- This device is sterile if package is unopened or undamaged.
- This device is non-pyrogenic.
- This device should be used by clinicians with adequate training in the use of the device.
- Follow appropriate anticoagulant protocols if required during the procedure.
- Prior to use, ensure the sheath and dilator are sized appropriately for the vessel and devices being used

CLINICAL BENEFIT STATEMENT

The Prelude Roadster Guide Sheath provides an indirect clinical benefit to the patient by facilitating access and the percutaneous introduction of various devices into a patient's vasculature, while maintaining hemostasis for a variety of diagnostic and therapeutic procedures.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

For a copy of this device's current European Summary of Safety and Clinical Performance (SSCP), please go the European database on medical devices (Eudamed), where it is linked to the basic UDI-DI. <https://ec.europa.eu/tools/eudamed>.

Pending availability of the Eudamed site, the SSCP may also be accessed at the following link: <http://www.merit.com/sscp/>

PATIENTS

The Prelude Roadster Guide Sheath product line is intended for patients that require percutaneous access for interventional therapies and diagnostic procedures per physician assessment.

USERS

These products are typically used by interventional radiologists, interventional cardiologists, diagnostic only cardiologists and vascular surgeons.

POTENTIAL COMPLICATIONS

Potential complications include, but are not limited to: air embolism, infection, hematoma, bleeding, perforation or laceration of the vessel wall, thrombus formation, pseudo aneurysm formation, guide wire embolization, vessel spasm, and risks normally associated with percutaneous diagnostic and/or interventional procedures. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established

INSPECTION PRIOR TO USE

Product is sterile if package is unopened and undamaged. Prior to use, carefully examine sheath to verify that the sterile package or product has not been damaged in shipment. Prior to and during use, inspect the sheath carefully for bends or kinks which may have occurred. Do not use if sheath or packaging are damaged.

REUSE PRECAUTION STATEMENT

For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

INSTRUCTIONS FOR USE

The following instructions provide technical direction but do not replace the necessity of formal training in the use of the device. The techniques and procedures described do not represent all medically acceptable protocols, nor are they intended as a substitute for the clinician's experience and judgment in treating any specific patient.

1. Identify and prepare the insertion site using proper aseptic technique and local anesthesia as required.
2. Remove the Prelude Roadster Guide Sheath components from package using proper aseptic technique.
3. Ensure the valve and the sheath are connected before use. Do not over-tighten the luer fitting of the hemostasis valve adapter.
4. Flush all components with heparinized saline or suitable isotonic solution. Be sure to wet the outer surface of the sheath introducer to activate the

hydrophilic coating with wet gauze or place in saline basin prior to use. The sheath should not be used in a dry state.

NOTE: Be sure to use in a hydrated state prior to use.

WARNING: After flushing side port, turn stopcock to off position to maintain flush in the side port and prevent bleed back upon insertion into the vessel.

WARNING: Do not wipe outer surface of the sheath introducer with dry gauze.

5. Insert vessel dilator into Prelude Roadster Guide Sheath through the center of the hemostasis valve (HVA) and snap into place.

WARNING: Dilator must be securely snapped into place to avoid damage to the vessel.

6. Insert appropriate access needle into vessel. While holding the access needle, place the flexible end or J end of the guide wire through access needle into vessel.

NOTE: Refer to product labeling for appropriate guide wire compatibility with the system components.

WARNING: Do not advance the guide wire if resistance is met. Determine the cause of resistance before proceeding.

7. Hold guide wire in place while removing access needle. Apply manual pressure above puncture site during needle removal and until the introducer/ dilator assembly is placed.

WARNING: Do not withdraw the guide wire through the needle because it may damage the guide wire.

8. Insert the introducer/dilator assembly over the guide wire into the vessel using a rotating motion, advance the introducer/dilator assembly through the tissue into the vessel.

WARNING: Ensure that the surface of the sheath is hydrated prior to insertion; the sheath should not be used in a dry state.

NOTE: Placing sheath in saline basin prior to use will ensure adequate hydration before placing into the vessel.

WARNING: During insertion, hold assembly near distal tip while passing over the guide wire and into the vessel to avoid buckling.

9. After introducer/dilator assembly has been placed into vessel and target destination has been reached, detach the dilator from the HVA by bending the dilator hub down slightly (this will unsnap the dilator hub from the introducer HVA cap). While holding the sheath, carefully remove the dilator leaving the sheath introducer in the vessel.

10. Aspirate from the side port extension to remove any potential air or debris. After aspiration, flush the side port with a suitable solution.

WARNING: Stopcock handle must be turned to the off position (toward the sheath hub) to prevent inadvertent blood loss.

11. Hold the sheath in place when inserting, positioning, or removing the devices to prevent un-intentional sheath movement. Always exchange or remove devices slowly through the sheath.


NOTE: Use caution when inserting and removing selected device(s) (wires, catheters, etc.) into Prelude Roadster Guide Sheath.

12. REMOVAL: The sheath should be removed at the end of the procedure or when clinically indicated. **Prior to removal of the sheath, reinsert the dilator and secure dilator to the sheath. Remove the sheath/dilator as a whole unit.** Compression on the vessel, above the puncture site, should be started as the sheath is slowly removed. Non-occlusive compression or closure device should be used to achieve hemostasis once the sheath is removed.

NOTE: Ensure all devices are removed at the end of the procedure

NOTE: Collected fibrin at the tip of the sheath may be aspirated via the side arm tubing prior to removal of the sheath.

13. DISPOSAL: After use, dispose of product and packaging in accordance with hospital protocol.

SYMBOL	DESIGNATION
	Use by date: YYYY-MM-DD
	Batch code
	Catalog number
	Sterilized using ethylene oxide
	Do Not Use If Package is Damaged and Consult Instruction for Use
	Single use
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Do not resterilize
	Caution
	Date of Manufacture: YYYY-MM-DD
	Medical Device
	Single Sterile Barrier System
	Authorized Representative in European Community
	Manufacturer
	Consult Instructions for Use. For electronic copy scan QR Code, or go to www.merit.com/ifu and enter IFU ID Number. For printed copy, available within seven days, call U.S.A or E.U. Customer Service
	Unique Device Identifier
	Max Guide Wire
	2D barcode

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Prelude Roadster is a registered and unregistered trademark of Merit Medical Systems, Inc
Merit Medical Systems, Inc.
South Jordan, Utah 84095 U.S.A.
U.S.A Customer Service 1-800-356-3748
Merit Medical Ireland Ltd,