BioFlo DuraMax[®]

with ENDEXO[®] Technology

Chronic Hemodialysis Catheter

EN	A Patient's Guide
FR	Guide du patient
IT	Guida per il paziente
DE	Patientenleitfaden
ES	Guía para el paciente
PT	Guia do Doente
PTBR	Um guia para o paciente
NL	Een patiëntengids
SV	En patientguide
NO	En pasientveiledning
DA	En vejledning til patienter
EL	Οδηγός για τον ασθενή
TR	Hasta Kılavuzu
RU	Руководство для пациентов
PL	Przewodnik dla pacjenta
CS	Příručka pro pacienta
BG	Ръководство за пациента
HU	Útmutató betegeknék
LV	Pacienta rokasgrāmata
LT	Paciento vadovas
RO	Ghid pentru pacient
SK	Príručka pre pacienta
HR	Vodič za pacijenta
ET	Patsiendi juhend
SL	Vodič za bolnike
FI	Potilasopas
KO	환자 가이드
MK	Водич за пациентот
AR	دليل المريض

BioFlo DuraMax[®]

with ENDEXO[®] Technology

Chronic Hemodialysis Catheter

As a patient or companion to someone who has received a BioFlo DuraMax[®] with ENDEXO[®] Technology Chronic Hemodialysis Catheter, here are some valuable things to know:

- The BioFlo DuraMax with ENDEXO Technology Chronic Hemodialysis Catheter provides access to the bloodstream in order to perform hemodialysis.
- Be aware of your dialysis schedule and know when, where, and how often you typically go for treatment.
- If traveling, be sure to set up appointments at facilities local to your destination so that you can maintain your current dialysis schedule, unless otherwise specified by your healthcare provider.
- Information about your BioFlo DuraMax with ENDEXO Technology Chronic Hemodialysis Catheter is available by calling the Merit Medical Customer Service at +1-800-356-3748.

CAUTION: Federal Law (USA) restricts this device for sale by or on the order of a physician.

Important instructions about Implant Card

- You will receive an Implant Card when receiving your BioFlo DuraMax with ENDEXO Technology Chronic Hemodialysis Catheter.
- Always carry your BioFlo DuraMax with ENDEXO Technology Chronic Hemodialysis Catheter implant card with you. The implant card is conveniently sized to fit in a wallet.
- The implant card has important information about your catheter that will be beneficial for your healthcare provider.

The image shows a red rectangular card with rounded corners. At the top left, the text "BioFlo DuraMax®" is written in white, with "with ENDEXO Technology" in smaller white text below it. To the right of this is the Merit Medical logo, which consists of a stylized "M" icon followed by the words "MERIT MEDICAL" in white. In the center of the card, the text "Patient Implant Card" is written in white. Below this, there is a large white rectangular area with a thin black border. Inside this area, the text "PLACE THE REMOVABLE STICKER FROM THE PRODUCT LABEL HERE" is written in black. At the bottom left of the white area, there is a QR code. To the right of the QR code, the text "Scan or go to URL <https://cloud.merit.com/catalog/IFUs/408039002.pdf> or a hard copy may be requested from Merit Customer Service: USA 1-800-356-3748 or EC +31 43 3588222" is written in black. At the bottom right of the white area, there is a small black icon of a factory or building. To the right of this icon, the text "Merit Medical Systems, Inc. 1600 West Merit Parkway South Jordan, Utah 84095 U.S.A." is written in black. Below this, there is a small black box with the text "EC REP" in white. To the right of this box, the text "Merit Medical Ireland Ltd Parkmore Business Park West Galway, Ireland" is written in black.

Catèter para hemodialisis crónica / Cathéter d'hémodialyse chronique / Dauerkatheter für die Hämodialyse / Catetere per emodialisi di mantenimento / Katheter voor chronische hemodialyse / Kronisk hæmodialysekateter / Kronisk hemodialyskateter / Καθετήρας χρόνιας αιμοκάθαρσης / Cateter para hemodiálise crónica / Állandó hemodializis-katéter / Katéřr pro chronickou hemodialýzu / Cewnik do hemodializy długoterminowej / Kateter for kronisk hemodialyse / 만성 혈액투석 카테터 / Kronik Hemodializ Kateteri / Ilgalaikis hemodializės kateteris / Regulaarse hemodialüüsi kateeter / Cateter de hemodializă cronică / Katéter na chronickú hemodialýzu / Kateter za kroničnu hemodijalizu / Катетер за хронична хемодијализа / Kateter za hroničnu hemodijalizu / Катетер для хронического гемодиализа / Hemodialyysikatetri kroonilised sairauden hoitoon / زمزملي لڳالا لیسڤ راسدو / Hroniskās hemodializes katetrs / Катетър за хронична хемодиализа / Kateter za kronično hemodializo

408039001 001 ID 2024-03-21

Information about your BioFlo
DuraMax with ENDEXO Technology
Chronic Hemodialysis Catheter is
available by calling the
Merit Medical U.S.A .Customer
Service 1-800-356-3748

BioFlo DuraMax[®] with ENDEXO[®] Technology Chronic Hemodialysis Catheter

Intended Use:

The BioFlo DuraMax with ENDEXO Technology Chronic Hemodialysis Catheter is intended to provide Long-Term vascular access for Hemodialysis in adults.

Catheters greater than 40 cm are intended for femoral vein insertion.

This patient guide is intended to be educational and is not a substitute for the directions for use provided with the device.

* Merit and the Merit logo are trademarks and/or registered trademarks of Merit Medical Systems, Inc., an affiliate or subsidiary. All other trademarks are property of their respective owners.

BioFlo DuraMax

with ENDEXO Technology

Chronic Hemodialysis Catheter

Table of Contents

Introduction.....	6
What is a BioFlo DuraMax with ENDEXO Technology Chronic Hemodialysis Catheter?	6
What do BioFlo DuraMax with ENDEXO Technology Chronic Hemodialysis Catheters consist of?	7
How is the BioFlo DuraMax with ENDEXO Technology Chronic Hemodialysis Catheter inserted?	8
How is your hemodialysis catheter used and accessed?.....	10
Care and Maintenance of your BioFlo DuraMax with ENDEXO Technology Chronic Hemodialysis Catheter	11
Identification.....	11
Commonly Asked Questions	11
Important information your healthcare provider should know	13
Patient Check List	14
Commonly used terms.....	15
Additional Instructions.....	15
Notes:.....	16

Introduction

Your healthcare provider has recommended the placement of your catheter in order that you can begin treatment to best manage your healthcare needs. While other options for long term vascular access for hemodialysis may be available to you, the placement of the BioFlo Duramax with ENDEXO Technology Chronic Hemodialysis Catheter makes it possible for you to begin your treatment plan as soon as your healthcare provider recommends. Other options for long term vascular access may require additional time to implant or create, and time for healing may be required before use. In some cases, your hemodialysis catheter may be used to manage your healthcare needs until a different type of vascular access is placed and ready for use. It is important to periodically discuss all options for long-term vascular access with your healthcare provider as these needs may change over time.

NOTE: A Summary of Safety and Clinical Performance on this device, is available on Eudamed at ec.europa.eu/tools/eudamed, where it is linked to the basic UDI-DI. The basic UDI-DI for the device is 0884450BUDI651Q6. Pending availability of the Eudamed site, the SSCP may also be accessed at the following link: <http://www.merit.com/sscp/>

What is Hemodialysis and why is it needed?

Normal functioning kidneys filter excess water and waste from your blood which keeps the composition of the blood stable. If your kidneys are damaged, dysfunctional, or missing, these harmful wastes build up. A hemodialysis machine, which functions in place of a kidney, can be used to cleanse these wastes from your blood. The hemodialysis catheter is inserted into a vein and allows for the exchange of blood to and from the hemodialysis machine.

What is a BioFlo DuraMax with ENDEXO Technology Chronic Hemodialysis Catheter?

The BioFlo DuraMax with ENDEXO Technology Chronic Hemodialysis Catheter is a tunneled hemodialysis catheter. A hemodialysis catheter is a hollow, soft tube that has two openings—one to send your blood to the hemodialysis machine and the other to return the cleansed blood back to your body. Each opening has a clamp and a hub with a cap. The clamp is kept closed and the endcap is kept on whenever the catheter is not being used. The term “tunneled” means a portion of the catheter is implanted under the skin before it is inserted into a vein. Tunneling the catheter helps prevent it from moving, which is important for long term hemodialysis treatments.

What do BioFlo DuraMax with ENDEXO Technology Chronic Hemodialysis Catheters consist of?

The BioFlo Duramax with ENDEXO Technology Chronic Hemodialysis Catheter is a long tube constructed of polyurethane material which contains an additional additive that has been shown to be effective in reducing thrombus formation (blood clotting) on the surface of the catheter. Special tests were conducted in a laboratory setting to determine the number of platelets that accumulate on the device when it is exposed to blood. Platelets are one part of “blood clotting” and this clotting can lead to difficulty in removing blood from the catheter and/or returning blood to your body through the catheter. Catheter shaft is made of Carbothane 3585A (6.43g), with 20% barium sulfate for radiopacity (1.61g), 2% Endexo plastic polymer (0.164g), and 0.2% teal colorant (0.0164g). It has permanent exposure to blood.


The tube is divided into separate sections so that the blood that is being removed from the body does not mix with the blood being returned to the body. This is important for your treatment to be successful.

The end of the catheter that will be outside of your body has a “y” shape. This “y” shape consists of two separate tubes called “extension legs”. The extension legs each have hubs and clamps described further here:

- The hub is the very end of the catheter. The hub is the part of the catheter that is utilized to connect tubing and administer treatment to the patient. The hub should always have a cap on it when the catheter is not in use to prevent anything from getting into the catheter such as debris or bacteria.
- The clamps can be set to an open or closed position. If the clamps are open, they allow treatment to flow through the catheter. If it is closed, it prevents this. When not in use, the clamps should remain in the closed position. Each clamp is labeled with its associated priming volume. This is the amount of fluid required to fill the length of the catheter.

The catheter also features the following:

- A rotating suture wing that allows the clinician to suture the catheter into place. This helps prevent movement of the catheter, keeping the tip that is inside the body in the appropriate position for treatment.
- A “cuff” that encircles the catheter and is intended to remain in the “tunnel” between the site where the catheter exits the skin and the site where the catheter enters the vessel. The cuff provides a barrier along the tunnel that may help to prevent bacteria from entering the tunnel and also may promote tissue adherence to the catheter. This adherence may help to keep the catheter from moving in or out of the insertion site.


How is the BioFlo DuraMax with ENDEXO Technology Chronic Hemodialysis Catheter inserted?

Insertion of the hemodialysis catheter is a brief procedure done under local or general anesthesia, or sedation. The skin on the chest is first cleaned to remove germs. There are a variety of techniques to successfully place a hemodialysis catheter. A common technique involves puncture of the vessel where the catheter will be placed with a small needle. Then, a special device is placed in the vessel that will be used to feed the catheter into the vessel. Before the catheter is placed, a small incision about 8-10cm below the vessel puncture site is made. The catheter will be tunneled under the skin from the small incision to the vessel puncture site, where the catheter is then inserted into the vein. The tip of the catheter is then positioned near the heart as this is the location in the body with the

highest blood flow. The proper catheter placement is typically confirmed using fluoroscopy or chest radiograph. The catheter may be sutured in place or a special adhesive securement device may be used to hold the catheter in place. The site will then be covered with a dressing. The healing process from this procedure takes approximately one week but may vary from patient to patient. Follow your healthcare provider's instructions on the care of your incision as well as instructions on when your catheter can be utilized.

Dialysis Catheter Insertion Site

*This part of the catheter is
tunneled under your skin.*

How is your hemodialysis catheter used and accessed?

Your hemodialysis catheter will be used for hemodialysis according to your healthcare team's instructions. Typically, hemodialysis will occur three times per week. At each dialysis session, your catheter will be accessed by attaching the dialysis tubing to your catheter. Careful attention should be paid to make sure that the hubs of your catheter are cleansed well before the catheter is used. Each hospital or facility has procedures that will be followed to prevent infection of your catheter.

Once the dialysis catheter is connected to the hemodialysis machine, the blood will be pulled from the catheter into the machine, followed by the filtration of impurities in the blood. The filtered blood will then be returned from the machine into the opposite lumen of the catheter. Dialysis sessions typically take several hours. Once dialysis is complete, the catheter will be flushed with saline and a locking solution that helps to keep the catheter flowing during dialysis sessions.

Typically, the exit site dressing will be changed as part of your dialysis visit. The catheter should remain secured for entire time it is implanted.

- Do NOT swim or shower.
- Do NOT soak your dressing while bathing.
- If your dressing becomes wet or loose, it should be changed. Discuss the plan for unscheduled dressing changes with your healthcare provider to determine if you or your caregiver will need to receive instructions. Alternately, your healthcare provider may ask that you return to your healthcare setting for any unscheduled dressing changes.
- Do NOT use acetone on any part of the catheter tubing. Exposure to this agent may cause catheter damage.
- Do not use sharp instruments near the extension tubing or catheter lumen.
- Do not use scissors to remove dressing.
- Catheter will be damaged if clamps other than what is provided with this kit are used.
- In the event a clamp breaks, replace the catheter at the earliest opportunity.
- Examine catheter lumen and extensions before and after each treatment for damage.
- Clamping of the tubing repeatedly in the same location may weaken tubing. Avoid clamping near the luers and hub of the catheter.
- Use only Luer Lock (threaded) connectors with this catheter.
- Repeated over tightening of bloodlines, syringes, and caps will reduce connector life and could lead to potential connector failure.

Care and Maintenance of your BioFlo DuraMax with ENDEXO Technology Chronic Hemodialysis Catheter

Your health care provider will give you specific instructions on how to care for your catheter, but you should follow these general guidelines:

Identification

How to determine that your dialysis catheter is a BioFlo DuraMax with ENDEXO Technology Chronic Hemodialysis Catheter

There are several ways to determine that you have a BioFlo DuraMax with ENDEXO Technology Chronic Hemodialysis Catheter. Upon receiving your catheter, your healthcare provider provided you with an implant card. This lets your healthcare provider know that you have a BioFlo DuraMax with ENDEXO Technology Chronic Hemodialysis Catheter. The name of the catheter is also printed on the catheter itself, above the suture wings but below the clamps. You should always keep your implant card with you and present it to your healthcare provider prior to treatment.

Commonly Asked Questions

Q How long will I have my BioFlo DuraMax with ENDEXO Technology Chronic Hemodialysis Catheter?

A The appropriate time to remove the catheter varies from patient to patient. It is important to discuss your personal needs and preferences with your care team to determine when the catheter will be removed.

Q How do I take care of my BioFlo DuraMax with ENDEXO Technology Chronic Hemodialysis Catheter?

A General instructions for care are provided in this brochure; however, your care team may provide additional instructions according to the protocols of the institution.

Q Will my catheter affect my daily activities?

A You should become familiar with the instructions in this brochure as well as any additional instructions regarding your daily activities from your care team. Keep in mind that the catheter site should remain clean and dry.

Q Will I have to change my dressing or flush my catheter?

A Generally, dressing changes and routine catheter flushing will be performed by your care team. If your dressing becomes soiled, wet, or is no longer adhering to your skin, the dressing should be changed. Your care team will instruct you if your catheter needs to be flushed between dialysis sessions for any reason.

Q Is my catheter safe in the MRI setting?

A The BioFlo DuraMax with ENDEXO Technology Chronic Hemodialysis Catheter does not contain any metal and therefore, there is no contraindication to having an MRI study with the catheter in place.

Q What do I do if a cap is loose or comes off of my catheter?

A First, make sure that the clamp on the extension leg is closed. Then, notify your care team immediately for further instructions on replacing the missing or loose cap.

Q What do I do if my catheter has a leak or a hole in it?

A While this is very unlikely, you should clamp the catheter such that clamp separates the hole or leak from the rest of the catheter if possible. The clamp position will be closest to the catheter exit site so that the catheter portion closest to your body is closed off from the hole or leak.

It is important to notify your healthcare provider(s) immediately about the following, as well as any new symptoms that you are concerned about:

- Aching discomfort, swelling or pain in the shoulder, neck or arm.
- Pain, redness, soreness or inflammation over or around the catheter site after the site has healed.
- Any unusual drainage from the catheter site or the needle puncture site (usually in the neck).
- Burning sensation or pain during any infusion or use of the catheter.
- Chills, shortness of breath, dizziness or a fever.
- Catheter leaking or damage to the catheter (see additional information in Commonly Asked Questions)

Important information your healthcare provider should know

Talk to your Healthcare Providers

As a patient who has a newly implanted hemodialysis catheter it is important to discuss potential complications that may occur as a result of the device placement, or the use of the catheter for long term access to the bloodstream. Information about potential complications associated with the device is available in the product packaging that is provided with the device. This information is also available by contacting your healthcare provider, or Merit Medical at +1-800-356-3748. Risks and benefits of the implant procedure and device placement should be discussed with your care team.

IMPORTANT:

Any serious incident which has occurred with the use of this device should be reported to your healthcare provider immediately. It is also important for product-related concerns or complaints to be shared with Merit Medical at www.merit.com

Product-related concerns or complaints can be submitted by patients and/or healthcare providers.

Notice for patients in Australia:

Please report product concerns/complaints to the Therapeutic Goods Administration (TGA) at <https://www.tga.gov.au/reporting-adverse-events>.

Notification of adverse event or complaint to Merit Medical or the regulatory body in your country is not a replacement for medical opinion.

Patient Check List

- ☐ You should keep your implant card with you and present it to your healthcare provider prior to treatment. If the healthcare provider needs additional information, they may call the Merit Medical Customer Service number at +1-800-356-3748.
- ☐ Speak up and share your concerns with your healthcare provider.
- ☐ Additional patient information and educational tools may be accessed by visiting www.merit.com/ifu.

Commonly used terms

Catheter

A flexible tube inserted through the skin into a blood vessel or cavity to draw out body fluid or infuse fluid.

Dialysis

A form of renal replacement therapy where a semipermeable membrane – in hemodialysis the membrane of the dialyzer – is used to selectively filter solute from the patient’s blood into the dialysate.

Hemodialysis (HD)

A treatment method where the patient’s blood flows outside the body through disposable bloodlines into a special filter, the dialyzer. The dialysis solution carries away waste products and excess water, and the cleaned blood is returned to the patient. The process is controlled by a hemodialysis machine that pumps blood, adds anticoagulants, regulates the purification process, and controls the mixing of the dialysis solution and its flow rate through the system. A patient typically receives three treatments per week, lasting from three to four hours.

Recommended Flushing Protocols

(to be completed by your healthcare provider)

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The image shows a red patient implant card for BioFlo DuraMax. At the top left, the text 'BioFlo DuraMax®' is written in white, with 'with ENDEXO Technology' below it. To the right is the Merit Medical logo. In the center, 'Patient Implant Card' is printed. Below this, a large white rectangular area is marked with red corner brackets and contains the instruction 'PLACE THE REMOVABLE STICKER FROM THE PRODUCT LABEL HERE'. At the bottom left, there is a QR code and a URL: 'Scan or go to URL https://cloud.merit.com/catalog/IFUs/408039002.pdf or a hard copy may be requested from Merit Customer Service: USA 1-800-356-3748 or EC +31 43 3588222'. At the bottom right, there is a silhouette of a mountain range, followed by the address 'Merit Medical Systems, Inc. 1600 West Merit Parkway South Jordan, Utah 84095 U.S.A.' and a box containing 'EC REP' next to the text 'Merit Medical Ireland Ltd Parkmore Business Park West Galway, Ireland'.

Catèter para hemodialisis crónica / Cathéter d'hémodialyse chronique / Dauerkatheter für die Hämodialyse / Catetere per emodialisi di mantenimento / Katheter voor chronische hemodialyse / Kronisk hæmodialysekateter / Kronisk hemodialyskateter / Καθετήρας χρόνιας αιμοκάθαρσης / Cateter para hemodiálise crónica / Állandó hemodializis-katéter / Katéřr pro chronickou hemodialýzu / Cewnik do hemodializy długoterminowej / Kateter for kronisk hemodialyse / 만성 혈액투석 카테터 / Kronik Hemodializ Kateteri / Ilgalaikis hemodializės kateteris / Regulaarse hemodialüüsi kateeter / Cateter de hemodializă cronică / Katéter na chronickú hemodialýzu / Kateter za kroničnu hemodijalizu / Катетер за хронична хемодијализа / Kateter za hroničnu hemodijalizu / Катетер для хронического гемодиализа / Hemodialyysikatetri kroonilised sairauden hoitoot / زمزملي لڳالا لیسڤ راسق / Hroniskās hemodializes katetrs / Катетър за хронична хемодиализа / Kateter za kronično hemodializo

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The image shows a red rectangular card with rounded corners. At the top left, the text "BioFlo DuraMax®" is written in white, with "with ENDEXO Technology" in smaller white text below it. To the right of this text is the Merit Medical logo, which consists of a stylized "M" icon followed by the words "MERIT MEDICAL" in a bold, sans-serif font. Below the company name, the words "Patient Implant Card" are printed in a large, white, sans-serif font. In the center of the card, there is a white rectangular area with a thin black border. Inside this area, the text "PLACE THE REMOVABLE STICKER FROM THE PRODUCT LABEL HERE" is printed in a small, black, sans-serif font. At the bottom left of the white area is a square QR code. To the right of the QR code, there is a block of text in a black, sans-serif font: "Scan or go to URL [https://cloud.merit.com/catalog/IFUs/408039002.pdf](\"https://cloud.merit.com/catalog/IFUs/408039002.pdf\") or a hard copy may be requested from Merit Customer Service: USA 1-800-356-3748 or EC +31 43 3588222". To the right of this text is a small black icon of a factory or building. Below the icon is a small white box with the text "EC REP" in black. To the right of this box is another block of text in a black, sans-serif font: "Merit Medical Systems, Inc. 1600 West Merit Parkway South Jordan, Utah 84095 U.S.A.". Below this text is another block of text in a black, sans-serif font: "Merit Medical Ireland Ltd Parkmore Business Park West Galway, Ireland".

Catèter para hemodialisis crónica / Cathéter d'hémodialyse chronique / Dauerkatheter für die Hämodialyse / Catetere per emodialisi di mantenimento / Katheter voor chronische hemodialyse / Kronisk hæmodialysekateter / Kronisk hemodialyskateter / Καθετήρας χρόνιας αιμοκάθαρσης / Cateter para hemodiálise crónica / Állandó hemodializis-katéter / Katéřr pro chronickou hemodialýzu / Cewnik do hemodializy długoterminowej / Kateter for kronisk hemodialyse / 만성 혈액투석 카테터 / Kronik Hemodializ Kateteri / Ilgalaikis hemodializės kateteris / Regulaarse hemodialüüsi kateeter / Cateter de hemodializă cronică / Katéter na chronickú hemodialýzu / Kateter za kroničnu hemodijalizu / Катетер за хронична хемодијализа / Kateter za hroničnu hemodijalizu / Катетер для хронического гемодиализа / Hemodialyysikatetri kroonilised sairauden hoitoon / زمزملي لڳالا لیسڊ راسڊو / Hroniskās hemodializes katetrs / Катетър за хронична хемодиализа / Kateter za kronično hemodializo

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The image shows a red patient implant card for BioFlo DuraMax. At the top left, the text 'BioFlo DuraMax®' is written in white, with 'with ENDEXO Technology' below it. To the right is the Merit Medical logo. In the center, the text 'Patient Implant Card' is displayed. Below this, a white rectangular area contains the instruction 'PLACE THE REMOVABLE STICKER FROM THE PRODUCT LABEL HERE'. At the bottom left of this white area is a QR code. To the right of the QR code is the text: 'Scan or go to URL <https://cloud.merit.com/catalog/IFUs/408039002.pdf> or a hard copy may be requested from Merit Customer Service: USA 1-800-356-3748 or EC +31 43 3588222'. At the bottom right of the white area is the Merit Medical Systems, Inc. logo and address: 'Merit Medical Systems, Inc. 1600 West Merit Parkway South Jordan, Utah 84095 U.S.A.'. Below this is the text: 'Merit Medical Ireland Ltd Parkmore Business Park West Galway, Ireland'.

Catèter para hemodialísis crónica / Cathéter d'hémodialyse chronique / Dauerkatheter für die Hämodialyse / Catetere per emodialisi di mantenimento / Katheter voor chronische hemodialyse / Kronisk hæmodialysekateter / Kronisk hemodialyskateter / Καθετήρας χρόνιας αιμοκάθαρσης / Cateter para hemodiálise crónica / Állandó hemodialízis-katéter / Katétr pro chronickou hemodialýzu / Cewnik do hemodializy długoterminowej / Kateter for kronisk hemodialyse / 만성 혈액투석 카테터 / Kronik Hemodializ kateteri / Ilgalaikis hemodializės kateteris / Regulaarse hemodialüüsi kateeter / Cateter de hemodializă cronică / Katéter na chronickú hemodialýzu / Kateter za kroničnu hemodijalizu / Катетер за хронична хемодијализа / Kateter za hroničnu hemodijalizu / Катетер для хронического гемодиализа / Hemodialyysikatetri kroonilised sairauden hoitoon / نمرزملی ایگلا لیسڈر راطسہ / Hroniskās hemodializes katetrs / Катѣть за хронична хемодиализа / Kateter za kronično hemodializo

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[illegible]

Manufacturer:
Merit Medical Systems, Inc.
1600 West Merit Parkway,
South Jordan, Utah 84095 U.S.A.
1-801-253-1600
U.S.A Customer Service 1-800-356-3748

Authorized Representative:
Merit Medical Ireland Ltd,
Parkmore Business Park West, Galway, Ireland
EC Customer Service +31 43 3588222

408039002_002 ID 2025-02-27

© 2023 Merit Medical Systems, Inc., or its affiliates.
All rights reserved.

Persons depicted in this brochure are models and
included for illustrative purposes only.