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Catalog number	Product Description	Components		Specification
MX-1100L-01	StabiliT MX Fracture Kit with PowerCURVE (Base, Long)	4457	StabiliT Bone Cement and Saturate Mixing System	-
		DF-5000A	Master Syringe Assembly (Master Syringe, Coupler and Elbow)	-
		1426	Locking Cement Delivery Cannula (11G) with Diamond-tipped Stylet	工作長度 (無探針): 14.5cm
		1472	StabiliT Introducer (10G) with Diamond-tipped Stylet	總長度: 12.5cm
		PWR-2011L	PowerCURVE Navigating Osteotome	工作長度: 17.1cm
		IN8KVCF	StabiliT Touch Delivery Syringe	操作壓力範圍: 0 - 35 ATM
MX-1100S-01	StabiliT MX Fracture Kit with PowerCURVE (Base, Short)	4457	StabiliT Bone Cement and Saturate Mixing System	-
		DF-5000A	Master Syringe Assembly (Master Syringe, Coupler and Elbow)	-
		0975	Locking Cement Delivery Cannula (11G) with Diamond-tipped Stylet	工作長度 (無探針): 12.5cm
		1493	StabiliT Introducer (10G) with Diamond-tipped Stylet	總長度: 10.5cm
		PWR-2011S	PowerCURVE Navigating Osteotome	工作長度: 15.1cm
		IN8KVCF	StabiliT Touch Delivery Syringe	操作壓力範圍: 0 - 35 ATM

Catalog number	Product Description	Components		Specification
MX-2100L-01	StabiliT MX Fracture Kit with PowerCURVE (Full, Long)	4457	StabiliT Bone Cement and Saturate Mixing System	-
		DF-5000A	Master Syringe Assembly (Master Syringe, Coupler and Elbow)	-
		1426	Locking Cement Delivery Cannula (11G) with Diamond- tipped Stylet	工作長度 (無探針): 14.5cm
		1472	StabiliT Introducer (10G) with Diamond-tipped Stylet	總長度: 12.5cm
		1467	StabiliT Introducer (10G) with Bevel-tipped Stylet	總長度: 12.5cm
		1545	VertecoR StraightLine Osteotome	套管總長度: 16.6cm
		PWR-2011L	PowerCURVE Navigating Osteotome	工作長度: 17.1cm
		IN8KVCF	StabiliT Touch Delivery Syringe	操作壓力範圍: 0 - 35 ATM
MX-2100S-01	StabiliT MX Fracture Kit with PowerCURVE (Full, Short)	4457	StabiliT Bone Cement and Saturate Mixing System	-
		DF-5000A	Master Syringe Assembly (Master Syringe, Coupler and Elbow)	-
		0975	Locking Cement Delivery Cannula (11G) with Diamond- tipped Stylet	工作長度 (無探針): 12.5cm
		1493	StabiliT Introducer (10G) with Diamond-tipped Stylet	總長度: 10.5cm
		1488	StabiliT Introducer (10G) with Bevel-tipped Stylet	總長度: 10.5cm
		1011	VertecoR StraightLine Osteotome	套管總長度: 14.6cm
		PWR-2011S	PowerCURVE Navigating Osteotome	工作長度: 15.1cm
		IN8KVCF	StabiliT Touch Delivery Syringe	操作壓力範圍: 0 - 35 ATM

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


Long/ short StabiliT Introducer 示意圖

Long/ short PowerCURVE Navigating Osteotome 示意圖


Figure 2: PowerCURVE Navigating Osteotome

Device	Exit Marker	Working Length
PowerCURVE Navigating Osteotome, Short	12.0 cm	15.1 cm
PowerCURVE Navigating Osteotome, Long	14.0 cm	17.1 cm