Manufacturer:

Merit Medical Systems, Inc.

1600 West Merit Parkway, South Jordan,

Utah 84095 U.S.A. 1-801-253-1600

U.S.A. Customer Service 1-800-356-3748

www.merit.com

Tenor™

STEERABLE GUIDEWIRE

Instructions for use.....	English	2
Instructions d'utilisation	Français (French)	4
Instrucciones de utilización	Español (Spanish).....	6
Instruções de uso.....	Português-Brazil (Portuguese)	8

DESCRIPTION

This Tenor Steerable Guidewire (0.014"/0.36 mm or 0.018"/0.46 mm) is a stainless steel guidewire with a polymer distal tip and a 3 cm radiopaque distal tip that is shapeable for 1 cm. The distal 50 cm is coated with a hydrophilic coating. The guidewire length is depicted on the product label.

INDICATIONS FOR USE

The Tenor Steerable Guidewire is intended to facilitate the placement of catheters within the peripheral vasculature for various diagnostic and interventional procedures.

WARNING

- Contents supplied sterile.
- Do not use if pouch is opened or damaged. Verify that the package integrity has been maintained to ensure the sterility of the device.
- Do not use in case of any surface irregularities, bends, or kinks. Any damage of the guidewire may change its characteristics likely to affect its performance
- Use the device prior to the "Use Before" date noted on the package.
- For single patient use only. Do not reuse, reprocess, or resterilize. Reusing, reprocessing or resterilizing may compromise the structural integrity of the device and or lead to device failure, which in turn may result in patient injury, illness, or death. Reusing, reprocessing or resterilizing may also create a risk of contamination of the device and or cause patient infection or cross infection including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
- Discard the product after use.
- This device should be used only by physicians thoroughly trained in percutaneous, intravascular techniques and procedures in relevant areas of the anatomy.
- Do not attempt to move the guidewire without observing the guidewire tip. Always maintain visualization of the guidewire under fluoroscopy.
- Do not push, pull, or rotate the wire against resistance. If resistance is met, discontinue movement of the guidewire, determine the reason for resistance, and take appropriate action before continuing. Movement of the catheter or guidewire against resistance may result in separation of the catheter or guidewire tip, damage to the catheter, or vessel perforation.

- The hydrophilic coating has an extremely lubricious surface only when properly hydrated.

PRECAUTIONS

- Do not expose guidewires to extreme temperatures or solvents.
- Prior to use, carefully examine the guidewire and packaging to verify proper function and integrity.
- Extreme care should be taken when shaping the guidewire distal tip. Over-manipulation of the guidewire distal tip may cause damage. Damaged guidewires should not be used.

COMPLICATIONS

Procedures requiring percutaneous catheter introduction should not be attempted by physicians unfamiliar with the possible complications listed below. Possible complications may include, but are not limited to, the following:

- Hemorrhage
- Infection
- Ischemia
- Perforation of vessel or arterial wall
- Puncture site hematoma
- Thrombus formation
- Vasospasm
- Vessel wall dissection

PREPARATION FOR USE

- Open carefully the pouch and remove the hoop from the pouch
- Flush the hoop with heparinized saline prior to guidewire removal.
- Remove the guidewire from the hoop.
- If the guidewire is to remain unused at any time during the procedure, be sure to rehydrate with heparinized saline prior to reinsertion.

INSTRUCTIONS FOR USE

- The distal tip of the guidewire can be carefully shaped to the desired tip configuration, using standard practices.
- Insert the steerable guidewire, flexible end first, into the catheter lumen using a guidewire insertion tool. Test the guidewire for free movement within the catheter. Then, advance the catheter into the guide catheter. During this step, exercise caution not to allow the guidewire to protrude beyond the catheter tip, this may damage it.
- To aid in rotating or steering the guidewire, secure the supplied torque device to the proximal end of the guidewire.

- In order to aid in the selective placement of the catheter into a particular vessel, gently rotate the proximal end of the guidewire as it is advanced forward.
- Use accepted angiographic techniques to steer the guidewire to the intended location.

Warning: Always maintain visualization of the guidewire under fluoroscopy, ensuring that the tip is moving freely when torque is applied.

- When the desired guidewire position is achieved, secure the guidewire in place while tracking the catheter over it and to the treatment location.
- Once the micro-catheter is in position, gently remove the guidewire prior to any intervention.

CAUTION:

⚠ Only: Federal (USA) law restricts this device to use by or on the order of a licensed physician.

CONTRAINDICATIONS:

There are no known contraindications for the use of steerable guidewires.

STORAGE:

Store the Tenor Steerable Guidewire in a cool, dark, dry area.

COMPATIBILITY:

The Tenor Steerable Guidewire is compatible with catheters which use 0.014” or 0.018” guidewires in intravascular procedures.

Merit Medical assumes no liability with respect to instruments reused, reprocessed or re-sterilized, and makes no warranties, expressed or implied, including, but not limited to, merchantability or fitness for intended use with respect to such instrument.

Information on packaging:

Symbol	Designation
	Manufacturer: Name & Address
	Use by date: year-month
	Batch code
	Catalogue number
	Do not resterilize
	Do not use if package is damaged
	Keep away from sunlight
	Keep dry
	Do not re-use
	Caution
	Consult Instructions For Use
	Sterilized using ethylene oxide

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